

Protocol No.: 20256031708-BC1

**Clinical Trial on Performance and Safety of
Disposable Powered Articulating Endoscopic Linear
Cutter Stapler in Total-Thoracoscopic Anatomic
Lobectomy (segmentectomy)**

(Prospective, multi-center, stratified block randomization, incomplete blinding, parallel positive control, non-inferiority testing)

Study Device: Disposable Powered Articulating Endoscopic Linear Cutter Stapler

Product Specifications: See the IFU for details.

Regulatory Category of Study Medical Device: Class III medical device

Requiring Clinical Trial Review and Approval: Yes ☐ No ☒

Similar Product(s) Existing in China: Yes ☒ No ☐

Protocol Version No.: 1.1

Protocol Version Date: 2019.03.25

Leading Site: Northern Jiangsu People's Hospital

Coordinating Investigator: Yu Shusheng

Sponsor: Fengh Medical Co., Ltd.

Agent :/

Northern Jiangsu People's
Hospital
★
National Drug Clinical Trial
Institute

Fengh Medical Co., Ltd.
★
3202811906439

Fengh Medical Co.,

Ltd.
★
3202811906439

Confidentiality Statement

All information contained in this report is proprietary to Fengh Medical Co., Ltd. Therefore, it is only provided to the investigators, co-investigators, ethics committees, regulatory authorities and other relevant medical institutions for review. Without the written approval of the sponsor, no information should be given to any third party not involved in the study, except for necessary explanation when signing the informed consent form with the subjects who may participate in the study.

Notes:

1. To carry out clinical trials of medical devices, clinical trial protocols shall be developed.
2. The clinical trial protocol is jointly designed and developed by the investigators of the clinical study institution and the sponsor. The sponsor shall sign the clinical trial protocol approved by both parties, and clinical study agreement or contract with the clinical study institutions and investigators.
3. The clinical trial protocol shall be approved by the ethics committee of this institution.
4. The clinical study of investigational medical devices listed in the *Catalogue of Class III Medical Devices Requiring Clinical Trial Review and Approval* shall not be implemented until approved by NMPA.
5. Clinical trial institutions and investigators shall jointly formulate the case number and duration of clinical study for each disease according to clinical statistical methods, characteristics of investigational medical devices, existing data and evidence with the sponsor, thereby ensuring to achieve the intended purpose of the study and reduce the waste of resources.
6. For a multi-center clinical trial, only the leading site should be filled in for the "clinical study institution" on the cover, and other institutions should be listed in the protocol.
7. For a multi-center clinical trial, the coordinating investigator shall be filled in for the "investigator" on the cover.

Content

Confidentiality Statement	1
Abstract of Trial Protocol	1
List of Abbreviations	6
Trial Flow Chart	7
Text of Trial Protocol	8
1. Sponsor Information	8
1.1 Sponsor Name	8
1.2 Sponsor Address	8
1.3 Sponsor Contact	8
1.4 Sponsor-Related Qualification Documents	8
2. List of All Clinical Trial Institutions and Investigators in Multi-Center Clinical Trial	8
3. Purpose and Content of Clinical Trial	9
3.1 Purpose	9
3.2 Content	9
4. Background Data of Clinical Trial	9
4.1 Pulmonary Lobectomy	9
4.2 Study Device	10
5. Product Model, Structural Composition, Operating Principle, Mechanism of Action and Trial Scope	11
5.1 Product Model/Specification	11
5.2 Product Structural Composition	15
5.3 Product Performance	15
5.4 Trial Scope	16
6. Product Applicable Scope, Contraindications and Precautions	16
6.1 Applicable Scope	16
6.2 Contraindications	16
6.3 Cautions	16
6.4 Precautions	16
6.5 Application Method	16
7. Reference Medical Device-Related Information	16
7.1 Stapler Information	16
7.2 Reload Information	17
8. Overall Design	18

8.1 Test Design	18
8.2 Test Process	35
8.3 Monitoring Plan	40
9. Statistical Considerations	41
9.1 Statistical Design, Methods and Analysis Procedure	41
9.2 Sample Size Calculation	41
9.3 Significance Level and Power of Clinical Trial	42
9.4 Expected Dropout Rate	42
9.5 Qualification/Disqualification Criteria for Clinical Trial Results	43
9.6 Criteria and Reasons for Terminating Trial Based on Statistical Reasons	43
9.7 Statistical Methods of Data, Together with Processing Methods of Missing, Unused or Wrong Data and Unreasonable Data	43
9.8 Procedures for Reporting Deviations from Original Statistical Plan	44
9.9 Selection Criteria and Reasons for Subjects Included in Analysis	44
9.10 Exclusion of Special Information During Hypothesis Validation and Reasons	45
10. Data Management	45
10.1 Database Design	45
10.2 Data Receipt and Data Entry	45
10.3 Management of Data Verification and Query	45
10.4 Database Lock and Data Archiving	46
11. Feasibility Analysis	46
11.1 Probability Analysis of Success	46
11.2 Probability Analysis of Failure	47
11.3 Conclusion	47
12. Quality Control of Clinical Trial	47
12.1 Quality Assurance	47
12.2 Quality Control	48
13. Ethical Issues and Informed Consent of Clinical Trial	49
13.1 Ethical Considerations	49
13.2 Approval of Study Protocol	49
13.3 Informed Consent Process and Informed Consent Form Text	50
13.4 Personal Data and Data Protection	50

14. Provisions for Reporting Adverse Events and Device Defects	51
14.1 Adverse Event	51
14.2 Serious Adverse Event	61
14.3 Reporting and Handling of Device Defects	62
15. Device Management	62
15.1 Study Device Management	62
15.2 Device Packaging	62
15.3 Recovery of Remaining Device	62
16. Deviations from Clinical Trial Protocol and Provision for Amendments to Clinical Trial Protocol	62
16.1 Deviations from Clinical Trial Protocol	62
16.2 Amendments to Clinical Trial Protocol	66
17. Direct Accessed Source Data/Source Documents	66
18. Finance and Insurance	66
19. Contents Covered in Clinical Trial Report	67
20. Confidentiality	67
21. Agreement on Publication of Study Results	68
22. Responsibilities Assumed by Each Party	68
22.1 Responsibilities of Sponsor	68
22.2 Responsibilities of Clinical Trial Institution and Investigators	69
23. References	71
24. Supplementary Provisions	71
Investigator Statement	74

Abstract of Trial Protocol

Trial Name	Clinical Trial on Performance and Safety of Disposable Powered Articulating Endoscopic Linear Cutter Stapler in Total-Thoroscopic Anatomic Lobectomy (segmentectomy)
Sponsor	Fengh Medical Co., Ltd.
Device Name	Study device: Disposable Powered Articulating Endoscopic Linear Cutter Stapler Control device: ECHELON Flex Powered Articulating Endoscopic Linear Cutters
Regulatory Category	Class III medical device
Trial Scope	The trial is intended for patients accepting anatomic lobectomy (segmentectomy).
Clinical Trial Indication	Anatomic lobectomy (segmentectomy)
Trial Design	Prospective, multi-center, stratified block randomization, incomplete blinding, parallel positive control, non-inferiority testing
Trial Purpose	By comparing the effectiveness and safety of Disposable Powered Articulating Endoscopic Linear Cutter Stapler (product under application) manufactured by Fengh Company with the similar product (ECHELON Flex Powered Articulating Endoscopic Linear Cutters manufactured by Johnson & Johnson Company) in thoroscopic anatomical lobectomy (segmentectomy), it is demonstrated that the product under application can be used for pulmonary tissue resection and anastomosis, and the clinical trial meets the requirements of <i>Good Clinical Practice for Medical Devices</i> and <i>Guidelines for Technical Review of Endoscopic Stapler Registration</i> , which can be used for product registration application.
Sample Size	164 cases 82 cases in both the test group and control group, respectively
Center Number	5
Follow-up period:	30 days
Test Time	(1) Enrollment date of the first subject: August 2018 (2) Duration of subject recruitment: 6 months (3) Completion date of the last subject: March 2019
Primary Effectiveness Evaluation Indicator	(1) Success rate of device cutting and anastomosis

Secondary Effectiveness Evaluation Indicator	(1) Operation time
	(2) Anastomosis time
	(3) Incidence of chest tube drainage greater than 7 days
	(4) Conversion rate to thoracotomy
	(5) Number of staple line repair needles
Safety Evaluation Indicator	(1) Incidence of staple line dehiscence
	(2) Incidence of staple line bleeding
	(3) Incidence of air leakage
	(4) Incidence of staple line infection
	(5) Incidence of SAEs
	(6) Incidence of AEs
	(7) Incidence of device defects
Operational Performance Evaluation	(1) Stapler operational performance evaluation
Inclusion Criteria	(1) Subjects aged 18-70 years (inclusive), male or female;
	(2) Subjects who are intended to undergo thoracoscopic anatomical lobectomy (segmentectomy);
	(3) Subjects or their guardians are able to understand the purpose of the study, show adequate compliance with the trial protocol, and sign an informed consent.
Exclusion Criteria	(1) Subjects with a contraindication to video-assisted thoracoscopic surgery;
	(2) Subjects have preoperative imaging findings indicating severe pleural adhesions and significantly calcified hilar lymph nodes;
	(3) Subjects have a platelet (PLT) $< 60 \times 10^9/L$ or international normalized ratio (INR) > 1.5 ;
	(4) Subjects with forced expiratory volume in one second (FEV1)/predicted value $\leq 50\%$, or forced expiratory volume in one second (FEV1)/forced vital capacity (FVC) $\leq 60\%$;
	(5) Subjects with ejection fraction $\leq 50\%$;
	(6) Subjects with major organ failure or other serious diseases (including clinically relevant cardiovascular disease or myocardial infarction within 12 months prior to enrollment; severe neurological or psychiatric history;

	<p>preoperative severe infection that must be controlled by medication; active disseminated intravascular coagulation; high risk of thrombosis);</p> <p>(7) The subject is a pregnant or lactating woman;</p> <p>(8) Subjects participating in other drug or device clinical trials within 1 month before the trial;</p> <p>(9) Other situations not suitable for inclusion judged by the investigators.</p>
Device Usage Specification	See the IFU for details.
Concomitant medication/concomitant treatment	If a subject receives any medication or treatment for an adverse event during the study for any reason, the investigators are required to record such occurrence on the case report form (CRF) for easy access.
Inspection Items	<p>Medical history and vital signs, etc.:</p> <p>(1) Demographic data: gender, date of birth;</p> <p>(2) History of present illness: admission time, etiology;</p> <p>(3) Past medical history: history of allergy, dysfunction of vital organs (heart, lung, liver and kidney), surgery and previous medication history (within 1 month);</p> <p>(4) Vital signs: height, weight, blood pressure, pulse, body temperature, respiration, BMI;</p> <p>Laboratory examinations:</p> <p>(1) Blood routine: hemoglobin (Hgb), red blood cell count (RBC), white blood cell count (WBC), neutrophil count (NEUT), neutrophil ratio (NEUT%), platelet count (PLT);</p> <p>(2) Coagulation blood routine: prothrombin time (PT), activated partial thromboplastin time (APTT), international normalized ratio (INR);</p> <p>(3) Blood biochemistry: aspartate aminotransferase (AST), alanine aminotransferase (ALT), serum total bilirubin (STB), serum albumin (ALB), serum creatinine (Scr), blood urea nitrogen (BUN);</p> <p>(4) Cardiac examination: ECG, cardiac score;</p> <p>(5) Pulmonary function tests: forced expiratory volume in one second (FEV1)/predicted, forced expiratory volume in one second (FEV1)/forced vital capacity (FVC).</p> <p>(6) Arterial blood gas analysis: pH, oxygen partial pressure, carbon dioxide</p>

partial pressure, oxygen saturation (SpO₂)

- (7) Blood pregnancy test: applicable to women of childbearing age (18 ~ 50 years old);

Preoperative imaging:

- (1) CT/MRI examination: whether there is pleural adhesion and degree, whether there are obvious calcified lymph nodes in the hilum;

Post-operative imaging:

- (1) Angiography: whether there is anastomotic bleeding;
(2) Thoracoscopy: whether there is air leak, staple line bleeding;

Other tests:

- (1) Intraoperative observation items: start time of operation, end time of operation, start time of anastomosis, end time of anastomosis, resected lobe or segment of lung, intraoperative blood loss, whether the transection and anastomosis are successful, whether the cutting staple line is complete, whether there is leakage of cutting staple line, whether there is bleeding of cutting staple line, whether there is air leakage of staple line, whether there is conversion to thoracotomy, number of staple line repair needles, whether there is dehiscence of staple line, whether there are other complications, evaluation of stapler operating performance;
(2) Device use: model and specification of electric stapler, model of staple cartridge, device defect;
(3) Postoperative follow-up observation: whether there is anastomotic line bleeding, whether there are other complications, and the duration of chest tube drainage (days);
(4) Routine observation items: adverse events/serious adverse events, concomitant and concomitant medication, protocol deviations.

Statistical analysis is performed using SAS 9.1 software.

Statistical Analysis

Perform statistical description on all data, including demographic data, baseline data, all efficacy indicators and all safety data. Measurement data include the mean, standard deviation, minimum, maximum, median, 25th percentile and 75th percentile; enumeration data include frequency and corresponding percentage. The primary effectiveness evaluation indicators are evaluated using the confidence interval method.

For secondary evaluation indicators, the difference test is adopted. Two independent samples t test is adopted for measurement data.

Difference test is used for safety evaluation indicators. Two independent samples t test is used for measurement data. Pearson χ^2 test or Fisher's exact test is used for enumeration data. The adverse events and device defects occurred in this trial are described in list.

Difference test is used to evaluate the operating performance of the stapler.

Pearson χ^2 test or Fisher exact test is used to analyze enumeration data.

Center effect is considered during the analysis process. General linear model is used for measurement data. CMH is used for enumeration data.

List of Abbreviations

Abbreviation	English	Chinese
AE	Adverse Events	不良事件
CRF	Case Report Form	病例报告表
EC	Ethics Committee	伦理委员会
FAS	Full Analysis Set	全分析集
GCP	Good Clinical Practice	医疗器械临床试验质量管理规范
ICF	Informed Consent Form	知情同意书
PPS	Per-Protocol Set	符合方案集
SAE	Serious Adverse Events	严重不良事件
SOP	Standard Operating Procedure	标准操作规程
SS	Safety Set	安全数据集
VATS	Video-Assisted Thoracic Surgery	电视辅助胸腔镜外科

Trial Flow Chart

Visit Time Visit Item	Visit 1	Visit 2	Visit 3	Visit 4
	Before surgery (Day -7 ~ 0)	During surgery (Day 0)	After surgery (Day 7 ± 2 days)	After surgery (Day 30 ± 5 days)
Informed consent	▲			
Subject screening	▲			
Demographic data	▲			
History of present illness	▲			
Past medical history	▲			
Vital signs	▲	▲	▲	▲
Blood routine	▲		▲	
Coagulation routine	▲		▲	
Blood biochemistry	▲		▲	
Cardiac examination	▲			
Pulmonary function test	▲			
Arterial blood gas analysis	▲			
Preoperative imaging	▲			
Blood pregnancy test ¹¹	▲			
Randomization	▲			
Intraoperative observation record		▲		
Use of device		▲		
Postoperative imaging ¹²			▲	▲
Postoperative follow-up observation			▲	▲
Concomitant medication	▲	▲	▲	▲
Adverse events	▲	▲	▲	▲
Device defect		▲		

Text of Trial Protocol

1. Sponsor Information

1.1 Sponsor Name

Feng Medical Co., Ltd.

1.2 Sponsor Address

Address/production address: D3 No. 6 Dongsheng West Road, Jiangyin National High-tech Zone, 214437
Jiangsu, P.R. China

1.3 Sponsor Contact

Contact: Liu Jiao

Tel.: (86) 510-81695550

Mobile: 13815132106

E-mail: Liu.jiao_fengmedical@163.com

1.4 Sponsor-Related Qualification Documents

Social credit code: 91320281583765063L

Registration inspection report: 2017ZC2908, 2017ZC2908-EMC, 2017QW0247

2. List of All Clinical Trial Institutions and Investigators in Multi-Center Clinical Trial

Code number	Name of clinical trial institution	Investigator	Title	Contact information
01	Northern Jiangsu People's Hospital	Shu Yusheng	Chief Physician	Tel: 18051061999
				E-mail: shuyusheng65@163.com
02	Taizhou Hospital of Zhejiang Province	Ma Dehua	Associate Chief Physician	Tel: 13305760862
				E-mail: madh@enzemed.com
03	Second Hospital of Ningbo	Zhao Guofang	Chief Physician	Tel: 13957420885
				E-mail: guofzhao@hotmail.com
04	Second Affiliated Hospital of Shantou University Medical College	Xu Xiaohua	Associate Chief Physician	Tel: 13509888111
				E-mail: xxh8111@sina.com

Code number	Name of clinical trial institution	Investigator	Title	Contact information
05	Shenzhen Second People's Hospital	Qian Youhui	Chief Physician	Tel: 13602585160
				E-mail: 13602585160@163.com

3. Purpose and Content of Clinical Trial

3.1 Purpose

By comparing the effectiveness and safety of Disposable Powered Articulating Endoscopic Linear Cutter Stapler (product under application) produced by Fengh Company with the similar product (ECHELON Flex Powered Articulating Endoscopic Linear Cutters manufactured by Johnson & Johnson Company) in VAT anatomical lobectomy (segmentectomy), it is proved that the product under application can be used for pulmonary tissue resection and anastomosis, and the clinical trial meets the requirements of *Good Clinical Practice for Medical Devices*, *Guidelines for Design of Clinical Trials of Medical Devices* and *Guidelines for Technical Review of Endoscopic Stapler Registration*, which can be used for product registration application.

3.2 Content

This clinical trial is a prospective study. The patients scheduled for thoroscopic anatomical lobectomy (segmentectomy) were selected as the study subjects. The subjects were randomly divided into the test group and control group. Subjects in the test group received the Disposable Powered Articulating Endoscopic Linear Cutter Stapler manufactured by Fengh Company, and the subjects in the control group received the similar product (ECHELON Flex Powered Articulating Endoscopic Linear Cutters) manufactured by Johnson & Johnson Company. The subjects were blinded and were followed up for 4 visits within 1 month. The collected data were analyzed and compared to evaluate the effectiveness and safety of the Disposable Powered Articulating Endoscopic Linear Cutter Stapler manufactured by Fengh Company in VAT anatomical lobectomy (segmentectomy).

4. Background Data of Clinical Trial

4.1 Pulmonary Lobectomy

With the development of video-assisted thoroscopic surgery (VATS) technology, some pulmonary parenchymal lesions that previously required thoracotomy for diagnosis and treatment could be completed by thoroscopic wedge resection, lobectomy, and pneumonectomy.

Lobectomy is indicated for peripheral lung cancer, irreversible lesions confined to the lobes. Under general anesthesia, an incision is made between the ribs to expose the lungs; the diseased lobe is removed after

examination of the thoracic cavity, and a drainage strip is placed to drain the air, effusion, and hemothorax in the pleural cavity; the thoracic cavity is then closed and the skin is sutured.

Lobectomy is indicated for lesions confined to the lobes but more than one pulmonary segment.

- (1) Cavitary lesions without significant absorption or enlargement after the tuberculous cavities are treated with standardized chemotherapy for 12 to 18 months. Patients with positive sputum, especially drug-resistant tuberculosis cases. Secondarily infected patients combined with hemoptysis and recurrent attacks. Patients who can not be excluded from cancerous cavities.
- (2) Tuberculoma: with a diameter greater than 3 cm, without change or enlargement after regular chemotherapy, unable to be excluded from the tumor.
- (3) After regular chemotherapy for 12 months for a large caseous lesion, sputum remained positive with hemoptysis.
- (4) Tuberculous stenosis of the lobar bronchi, resulting in atelectasis and pulmonary consolidation.
- (5) Bilateral lesions, while the main lesion is concentrated in one lobe and can be removed in stages and fractions.
- (6) Localized lesions in the lung caused by atypical acid-fast bacilli should be operated because there are no effective drugs and they are easy to develop and reactivate.
- (7) Early-stage squamous cell lung cancer (T₁N₀M₀).
- (8) Benign lung diseases, including pulmonary cysts, pulmonary inflammatory pseudotumors, etc.

Contraindications for lobectomy.

- (1) Patients whose lesions are unstable as shown by chest X-ray since 3 months.
- (2) Patients with irregular chemotherapy and insufficient course of medication.
- (3) Patients with moderate and severe cardiopulmonary insufficiency and poor general condition who can not tolerate surgery.
- (4) Patients complicated with other diseases (e.g., diabetes, hyperthyroidism), which can not be effectively controlled.

4.2 Study Device

Since VATS is increasingly and widely used in the treatment of lobectomy, the stapler cooperating with thoracoscopic surgery is also widely used in clinical practice. Stapler is an alternative manual suture device used in medicine. Its main operating principle is to use titanium nail for transection or anastomosis of tissues, which is similar to book stapler. According to different applicable scopes, it can be divided into skin stapler, linear stapler, circular stapler, circular hemorrhoidal stapler, circumcision stapler, vascular stapler, hernia stapler, etc.

Compared with the traditional manual suture, the advantages of anastomosis suture include: ① simple and

convenient operation, saving operation time; ② disposable use, avoiding cross-infection; ③ tight and moderate suture using titanium nail or stainless steel nail (skin stapler); ④ few side effects and effectively reducing surgical complications.

In order to provide clinicians with a more convenient thoracoscopic linear stapler, Fengh Company develops one kind of Disposable Powered Articulating Endoscopic Linear Cutter Stapler according to market demands. This linear stapler is driven by electric motor, which can better reduce the surgical intensity of doctors and save surgical time.

The Disposable Powered Articulating Endoscopic Linear Cutter Stapler designed and developed by Fengh Company according to the clinical treatment characteristics have passed the test by Jiangsu Institute for Medical Device Testing of NMPA (original CFDA) [2017ZC2908, 2017ZC2908-EMC, 2017QW0247].

This clinical trial is to test the safety and effectiveness of the device by evaluating the success rate of cutting and anastomosis of the device after use, so as to provide trial basis for the product registration.

5. Product Model, Structural Composition, Operating Principle, Mechanism of Action and Trial Scope

5.1 Product Model/Specification

Table 1 Model and specification of Electric stapler

Model/specification of electric stapler	Arm length (L3)	Length of reload-resisting base (L1)	Overall length (L)
FSMS30	192 ± 20	66.3 ± 2.5	538.6 ± 20
FSMS45	192 ± 20	78.3 ± 2.5	547.6 ± 20
FSMS60	192 ± 20	90.3 ± 2.5	562.6 ± 20
FSMM30	252 ± 20	66.3 ± 2.5	538.6 ± 20
FSMM45	252 ± 20	78.3 ± 2.5	547.6 ± 20
FSMM60	252 ± 20	90.3 ± 2.5	562.6 ± 20
FSML30	352 ± 20	66.3 ± 2.5	538.6 ± 20
FSML45	352 ± 20	78.3 ± 2.5	547.6 ± 20
FSML60	352 ± 20	90.3 ± 2.5	562.6 ± 20
FSAS30	192 ± 20	66.3 ± 2.5	538.6 ± 20
FSAS45	192 ± 20	78.3 ± 2.5	547.6 ± 20
FSAS60	192 ± 20	90.3 ± 2.5	562.6 ± 20
FSAM30	252 ± 20	66.3 ± 2.5	538.6 ± 20
FSAM45	252 ± 20	78.3 ± 2.5	547.6 ± 20
FSAM60	252 ± 20	90.3 ± 2.5	562.6 ± 20
FSAL30	352 ± 20	66.3 ± 2.5	538.6 ± 20
FSAL45	352 ± 20	78.3 ± 2.5	547.6 ± 20
FSAL60	352 ± 20	90.3 ± 2.5	562.6 ± 20
FNMS30	192 ± 20	66.3 ± 2.5	538.6 ± 20
FNMS45	192 ± 20	78.3 ± 2.5	547.6 ± 20
FNMS60	192 ± 20	90.3 ± 2.5	562.6 ± 20
FNMM30	252 ± 20	66.3 ± 2.5	538.6 ± 20
FNMM45	252 ± 20	78.3 ± 2.5	547.6 ± 20
FNMM60	252 ± 20	90.3 ± 2.5	562.6 ± 20

Clinical Trial on Performance and Safety of Disposable Powered Articulating Endoscopic Linear Cutter Stapler in Total-Thoracoscopic Anatomic Lobectomy (segmentectomy)

Model/specification of electric stapler	Arm length (L3)	Length of reload-resisting base (L1)	Overall length (L)
FNML30	352 ± 20	66.3 ± 2.5	538.6 ± 20
FNML45	352 ± 20	78.3 ± 2.5	547.6 ± 20
FNML60	352 ± 20	90.3 ± 2.5	562.6 ± 20
FNAS30	192 ± 20	66.3 ± 2.5	538.6 ± 20
FNAS45	192 ± 20	78.3 ± 2.5	547.6 ± 20
FNAS60	192 ± 20	90.3 ± 2.5	562.6 ± 20
FNAM30	252 ± 20	66.3 ± 2.5	538.6 ± 20
FNAM45	252 ± 20	78.3 ± 2.5	547.6 ± 20
FNAM60	252 ± 20	90.3 ± 2.5	562.6 ± 20
FNAL30	352 ± 20	66.3 ± 2.5	538.6 ± 20
FNAL45	352 ± 20	78.3 ± 2.5	547.6 ± 20
FNAL60	352 ± 20	90.3 ± 2.5	562.6 ± 20
DSMS30	192 ± 20	66.3 ± 2.5	538.6 ± 20
DSMS45	192 ± 20	78.3 ± 2.5	547.6 ± 20
DSMS60	192 ± 20	90.3 ± 2.5	562.6 ± 20
DSMM30	252 ± 20	66.3 ± 2.5	538.6 ± 20
DSMM45	252 ± 20	78.3 ± 2.5	547.6 ± 20
DSMM60	252 ± 20	90.3 ± 2.5	562.6 ± 20
DSML30	352 ± 20	66.3 ± 2.5	538.6 ± 20
DSML45	352 ± 20	78.3 ± 2.5	547.6 ± 20
DSML60	352 ± 20	90.3 ± 2.5	562.6 ± 20
DSAS30	192 ± 20	66.3 ± 2.5	538.6 ± 20
DSAS45	192 ± 20	78.3 ± 2.5	547.6 ± 20
DSAS60	192 ± 20	90.3 ± 2.5	562.6 ± 20
DSAM30	252 ± 20	66.3 ± 2.5	538.6 ± 20
DSAM45	252 ± 20	78.3 ± 2.5	547.6 ± 20
DSAM60	252 ± 20	90.3 ± 2.5	562.6 ± 20
DSAL30	352 ± 20	66.3 ± 2.5	538.6 ± 20
DSAL45	352 ± 20	78.3 ± 2.5	547.6 ± 20
DSAL60	352 ± 20	90.3 ± 2.5	562.6 ± 20
DNMS30	192 ± 20	66.3 ± 2.5	538.6 ± 20
DNMS45	192 ± 20	78.3 ± 2.5	547.6 ± 20
DNMS60	192 ± 20	90.3 ± 2.5	562.6 ± 20
DNMM30	252 ± 20	66.3 ± 2.5	538.6 ± 20
DNMM45	252 ± 20	78.3 ± 2.5	547.6 ± 20
DNMM60	252 ± 20	90.3 ± 2.5	562.6 ± 20
DNML30	352 ± 20	66.3 ± 2.5	538.6 ± 20
DNML45	352 ± 20	78.3 ± 2.5	547.6 ± 20
DNML60	352 ± 20	90.3 ± 2.5	562.6 ± 20
DNAS30	192 ± 20	66.3 ± 2.5	538.6 ± 20
DNAS45	192 ± 20	78.3 ± 2.5	547.6 ± 20
DNAS60	192 ± 20	90.3 ± 2.5	562.6 ± 20
DNAM30	252 ± 20	66.3 ± 2.5	538.6 ± 20
DNAM45	252 ± 20	78.3 ± 2.5	547.6 ± 20
DNAM60	252 ± 20	90.3 ± 2.5	562.6 ± 20
DNAL30	352 ± 20	66.3 ± 2.5	538.6 ± 20
DNAL45	352 ± 20	78.3 ± 2.5	547.6 ± 20
DNAL60	352 ± 20	90.3 ± 2.5	562.6 ± 20

Table 2 Model and specification of cartridge components

Unit: mm

Model/specification of reload	Reload color	Reload number	Anastomosis length (L2)
FMCC30	Gray	36	35.2 ± 2
FMCC45	Gray	70	49.3 ± 2
FMCC60	Gray	88	61.3 ± 2
FMCW30	White	36	35.2 ± 2
FMCW45	White	70	49.3 ± 2
FMCW60	White	88	61.3 ± 2
FMCT30	Brown	36	35.2 ± 2
FMCT45	Brown	70	49.3 ± 2
FMCT60	Brown	88	61.3 ± 2
FMCB30	Blue	36	35.2 ± 2
FMCB45	Blue	70	49.3 ± 2
FMCB60	Blue	88	61.3 ± 2
FMCG30	Green	36	35.2 ± 2
FMCG45	Green	70	49.3 ± 2
FMCG60	Green	88	61.3 ± 2
FMCP30	Purple	36	35.2 ± 2
FMCP45	Purple	70	49.3 ± 2
FMCP60	Purple	88	61.3 ± 2
FMCY30	Gold	36	35.2 ± 2
FMCY45	Gold	70	49.3 ± 2
FMCY60	Gold	88	61.3 ± 2
FMCX30	Black	36	35.2 ± 2
FMCX45	Black	70	49.3 ± 2
FMCX60	Black	88	61.3 ± 2
FACC30	Gray	36	35.2 ± 2
FACC45	Gray	70	49.3 ± 2
FACC60	Gray	88	61.3 ± 2
FACW30	White	36	35.2 ± 2
FACW45	White	70	49.3 ± 2
FACW60	White	88	61.3 ± 2
FACT30	Brown	36	35.2 ± 2
FACT45	Brown	70	49.3 ± 2
FACT60	Brown	88	61.3 ± 2
FACB30	Blue	36	35.2 ± 2

Clinical Trial on Performance and Safety of Disposable Powered Articulating Endoscopic Linear Cutter Stapler in Total-Thoroscopic Anatomic Lobectomy (segmentectomy)

Model/specification of reload	Reload color	Reload number	Anastomosis length (L2)
FACB45	Blue	70	49.3 ± 2
FACB60	Blue	88	61.3 ± 2
FACG30	Green	36	35.2 ± 2
FACG45	Green	70	49.3 ± 2
FACG60	Green	88	61.3 ± 2
FACP30	Purple	36	35.2 ± 2
FACP45	Purple	70	49.3 ± 2
FACP60	Purple	88	61.3 ± 2
FACY30	Gold	36	35.2 ± 2
FACY45	Gold	70	49.3 ± 2
FACY60	Gold	88	61.3 ± 2
FACX30	Black	36	35.2 ± 2
FACX45	Black	70	49.3 ± 2
FACX60	Black	88	61.3 ± 2
DMCC30	Gray	36	35.2 ± 2
DMCC45	Gray	70	49.3 ± 2
DMCC60	Gray	88	61.3 ± 2
DMCW30	White	36	35.2 ± 2
DMCW45	White	70	49.3 ± 2
DMCW60	White	88	61.3 ± 2
DMCT30	Brown	36	35.2 ± 2
DMCT45	Brown	70	49.3 ± 2
DMCT60	Brown	88	61.3 ± 2
DMCB30	Blue	36	35.2 ± 2
DMCB45	Blue	70	49.3 ± 2
DMCB60	Blue	88	61.3 ± 2
DMCG30	Green	36	35.2 ± 2
DMCG45	Green	70	49.3 ± 2
DMCG60	Green	88	61.3 ± 2
DMCP30	Purple	36	35.2 ± 2
DMCP45	Purple	70	49.3 ± 2
DMCP60	Purple	88	61.3 ± 2
DMCY30	Gold	36	35.2 ± 2
DMCY45	Gold	70	49.3 ± 2
DMCY60	Gold	88	61.3 ± 2

Clinical Trial on Performance and Safety of Disposable Powered Articulating Endoscopic Linear Cutter Stapler in Total-Thoroscopic Anatomic Lobectomy (segmentectomy)

Model/specification of reload	Reload color	Reload number	Anastomosis length (L2)
DMCX30	Black	36	35.2 ± 2
DMCX45	Black	70	49.3 ± 2
DMCX60	Black	88	61.3 ± 2
DACC30	Gray	36	35.2 ± 2
DACC45	Gray	70	49.3 ± 2
DACC60	Gray	88	61.3 ± 2
DACW30	White	36	35.2 ± 2
DACW45	White	70	49.3 ± 2
DACW60	White	88	61.3 ± 2
DACT30	Brown	36	35.2 ± 2
DACT45	Brown	70	49.3 ± 2
DACT60	Brown	88	61.3 ± 2
DACB30	Blue	36	35.2 ± 2
DACB45	Blue	70	49.3 ± 2
DACB60	Blue	88	61.3 ± 2
DACG30	Green	36	35.2 ± 2
DACG45	Green	70	49.3 ± 2
DACG60	Green	88	61.3 ± 2
DACP30	Purple	36	35.2 ± 2
DACP45	Purple	70	49.3 ± 2
DACP60	Purple	88	61.3 ± 2
DACY30	Gold	36	35.2 ± 2
DACY45	Gold	70	49.3 ± 2
DACY60	Gold	88	61.3 ± 2
DACX30	Black	36	35.2 ± 2
DACX45	Black	70	49.3 ± 2
DACX60	Black	88	61.3 ± 2

5.2 Product Structural Composition

- Powered stapler consists of stapler body, arm rod, reload-resisting base, reload base, battery pack or power adapter.
- The cartridge component is composed of reload, reload-pushing piece and anastomotic nail, etc.

5.3 Product Performance

The main technical indicators of electric stapler include convenient replacement of stapler components, reliable positioning. In addition, the stapler is equipped with a safety protection device for empty staple

cartridge to maintain its reliability. The product shall be sterile.

5.4 Trial Scope

The product is intended for thoracoscopic anatomical lobectomy (segmentectomy).

6. Product Applicable Scope, Contraindications and Precautions

6.1 Applicable Scope

The Disposable Powered Articulating Endoscopic Linear Cutter Stapler is intended for transection, resection, and/or creation of anastomoses. The instrument has application in multiple open or minimally invasive general, gynecologic, urologic, thoracic, and pediatric surgical procedures. It can be used with staple line or tissue buttressing materials. The instrument may also be used for transection and resection of liver parenchyma (hepatic vasculature and biliary structures), pancreas, kidney and spleen.

6.2 Contraindications

- ◆ Do not use the instruments on the aorta.
- ◆ Do not use the instruments on ischemic or necrotic tissue.
- ◆ Do not use any linear cutter on major vessels without making provision for proximal and distal control.
- ◆ Tissue thickness should be carefully evaluated before firing any stapler. Refer to Reload Product Codes for tissue compression requirement (closed staple height) for each staple size. If tissue cannot comfortably compress to the closed staple height, or easily compresses to less than the closed staple height, the tissue is contraindicated as it may be too thick or too thin for the selected staple size.
- ◆ These instruments are not intended for use when surgical stapling is contraindicated.

6.3 Cautions

See the IFU for details.

6.4 Precautions

See the IFU for details.

6.5 Application Method

See the IFU for details.

7. Reference Medical Device-Related Information

7.1 Stapler Information

Product name	ECHELON Flex Powered Articulating Endoscopic Linear Cutters
Specification/model	PSE45A, PCE45A, PLE45A, PSE60A, PCE60A, PLE60A
Registration No.	GXZJ 20173226240

Structural composition	The product is a sterile device for single patient use, consisting of closing rod, red firing lock, firing handle, anvil release button, lithium battery pack, battery pack-releasing piece, hole cover plate for manual override maintenance, knife-reversing switch, knob, articulation fins, reload alignment slot, anvil jaw and reload jaw. The lithium battery pack must be installed prior to use. The product has been sterilized by irradiation. This product is provided for single use.
Applicable scope	The Disposable Powered Articulating Endoscopic Linear Cutter Stapler is intended for transection, resection, and/or creation of anastomoses. The instrument has application in multiple open or minimally invasive general, gynecologic, urologic, thoracic, and pediatric surgical procedures. It can be used with staple line or tissue buttressing materials. The instrument may also be used for transection and resection of liver parenchyma (hepatic vasculature and biliary structures), pancreas, kidney and spleen.
Manufacturer	Ethicon Endo-Surgery, LLC
Agent	Johnson & Johnson Medical (Shanghai) Ltd.

7.2 Reload Information

Product name 1	echelon 60 ENDOPATH STAPLER Endoscopic Linear Cutter Reloads
Specification/model 1	ECR60W
Registration No. 1	GXZJ 20153080094
Structural composition 1	The linear cutter stapler reload is composed of cartridge, anastomotic nail and reload-protecting board, in which the material of anastomotic nail is Ti3Al2.5V. The reload is sterilized by radiation and provided for single use.
Product name 2	ECHELON 45 ENDOPATH Stapler Endoscopic Linear Cutter and Reloads
Specification/model 2	ECR45M, EC45, SC45, ECLG45, ECR45W, ECR45B, ECR45D, ECR45G
Registration No. 2	GXZJ 20163654883
Structural composition 2	The product is composed of anvil, trigger, firing trigger, handle and cartridge component. The reload is sterilized by radiation and provided for single use.
Product name 3	ECHELON FLEX Articulating Linear Cutters and Reloads Articulating Endoscopic Linear Cutter Reloads
Specification/model 3	EC45A, SC45A, EC45AL, EC60A, SC60A, LONG60A, ECR60W, ECR 60B, ECR60D, ECR60G, ECR60T, ECR60M
Registration No. 3	GXZJ 20173656323

Structural composition 3	This product is composed of stapler and cartridge. Its components include cartridge, closing rod, firing rod, rotary knob, articulating piece, manual knife-reversing switch, trip count indicator, knife direction indicator, anvil release button, closing rod, and anvil. The stapler is made of stainless steel, and the anastomotic needle is made of titanium alloy (Titanium3A12.5V). The product is sterilized by radiation.
Applicable scope	It is used for transection, resection and/or anastomosis. It can be used in a variety of open or minimally invasive general, gynecologic, urologic, thoracic, and pediatric surgical procedures. It can be used with staple line or tissue buttressing materials. These devices can also be used for transection and resection of liver parenchyma tissues (hepatic vasculature and biliary structures), pancreas, kidney and spleen.
Manufacturer	Ethicon Endo-Surgery, LLC
Agent	Johnson & Johnson Medical (Shanghai) Ltd.

8. Overall Design

8.1 Test Design

8.1.1 Primary objective

To evaluate whether the primary effectiveness evaluation indicator "device cutting and anastomosis success rate" of the Disposable Powered Articulating Endoscopic Linear Cutter Stapler manufactured by Fengh Company is non-inferior to the similar product manufactured by Johnson & Johnson Company in total-thoroscopic anatomic lobectomy (segmentectomy).

8.1.2 Secondary objectives

To evaluate whether there are differences in the secondary evaluation indicators (operation time, anastomosis time, incidence of chest tube drainage for more than 7 days, rate of conversion to thoracotomy, number of staple line repair needles) when the Disposable Powered Articulating Endoscopic Linear Cutter Stapler manufactured by Fengh Company is used in VAT anatomical lobectomy (segmentectomy) compared with the similar product manufactured by Johnson & Johnson Company.

To evaluate whether there are differences in the safety evaluation indicators (incidence of staple line dehiscence, incidence of staple line bleeding, incidence of air leakage, incidence of staple line infection, incidence of SAE (serious adverse event), incidence of AE (adverse event), incidence of device defects) between the Disposable Powered Articulating Endoscopic Linear Cutter Stapler manufactured by Fengh Company and the similar product manufactured by Johnson & Johnson Company in total-thoroscopic anatomic lobectomy (segmentectomy).

To evaluate whether there are differences in the operational performance evaluation indicators (operational performance evaluation of the stapler) between the Disposable Powered Articulating Endoscopic Linear Cutter Stapler manufactured by Fengh Company and the similar product manufactured by Johnson & Johnson Company in total-thoracoscopic anatomic lobectomy (segmentectomy).

8.1.3 Selection of test method and reasons

This clinical trial is a "prospective, multi-center, stratified-block randomization, incomplete blinding, parallel positive control, non-inferiority test".

- (1) Prospective: A research method that tracks from the present to the future. Prospective studies can clarify the causal relationship; uniform diagnostic, detection, and evaluation criteria can be available for the obtained data, the data processing therefore is controllable.
- (2) Multi-center: 5 clinical trial institutions were selected to conduct this clinical trial. The application of multi-center trial can allow to obtain more cases than those of single-center trial in the same time, the clinical trial time therefore can be shortened; since the multi-center trial is completed by different regions, different trial centers and many clinical researchers, the conclusion therefore is often widely representative.
- (3) Stratified block randomization: This trial adopts the stratified block randomization method (with the study site as the stratification factor) and the randomization envelope to randomize the subjects who are screened successfully. The randomized-grouping time is that before the surgery. Randomized grouping ensures that, except for treatment factors, other non-treatment factors that may have confounding effect are kept consistent in each group as far as possible, so as to ensure the balance of baseline of each group; randomized-grouping allows each subject to have equal opportunity to be assigned to the test group or control group.
- (4) Incomplete blinding: Since the investigational medical device and the control medical device are different in appearance, while the subjects are under anesthesia during the operation, the grouping arrangement can not be known. Therefore, the method of incomplete blinding is used, that is, the investigators knows the grouping result, while the subject does not know the grouping result. Blinded data review is also required at the same time. Incomplete blinding can minimize the bias caused by "placebo effect", and make the evaluation of adverse reactions by subjects more objective.
- (5) Parallel positive control: The investigational medical device has the mature application technology with many similar products that are marketed at home and abroad. It has been used in clinical practice for many years, with exact effect and predictable occurrence of adverse events. Therefore, parallel positive control is adopted in this clinical trial.

- (6) Non-inferiority test: To evaluate whether the main effectiveness evaluation indicator "success rate of cutting and anastomosis" of the Disposable Powered Articulating Endoscopic Linear Cutter Stapler manufactured by Fengh Company is non-inferior to the similar product manufactured by Johnson & Johnson Company in total-thoroscopic anatomic lobectomy (segmentectomy). The purpose of the non-inferiority test is to validate that if the effectiveness and safety of the investigational medical device is lower than those of the control device. The difference should be less than the pre-defined non-inferiority margin, that is, the difference should be acceptable in clinical practice.

8.1.4 Measures for reducing and avoiding bias

- (1) Before the start of the study, the sponsor shall provide relevant training for the investigators participating in the study to ensure that the investigators fully understand the study process and are familiar with the operation of the study device; during the study, the investigators shall operate in strict accordance with the operating methods and procedures in the study protocol, and the clinical trial monitor shall do a good job in quality control and monitoring to ensure that the investigators operate and implement in strict accordance with the study protocol. The above measures shall be implemented throughout the implementation phase of the study to reduce errors or operational errors.
- (2) Subjects shall be screened in strict accordance with the clinical diagnostic criteria, as well as the inclusion and exclusion criteria of the trial protocol to reduce selection bias.
- (3) VAT anatomical lobectomy (segmentectomy) is difficult to operate. Thoracoscopic anatomical lobectomy (segmentectomy) requires high operating skills and should be carried out in centers with rich experience in thoracoscopic surgery. In view of the learning curve of thoracoscopic anatomical lobectomy (segmentectomy), doctors can be more proficient based on at least 50 cases. The operating surgeon who is recommended to perform thoracoscopic lobectomy (segmentectomy) should be experienced in at least 100 cases of thoracoscopic anatomic lobectomy (segmentectomy).
- (4) When the clinical trial is completed, the data storage and collation shall be done well. When data problems are found, the data analyst shall check and confirm the data through the data clarification form to avoid recording errors.

8.1.5 Investigational medical device and control medical device or control diagnosis and treatment method

See the corresponding IFU for the diagnosis and treatment method with the investigational medical device and the reference medical device. The investigators shall truthfully record the usage for reference.

When subjects experience adverse events related to the stapler, the investigators should truthfully judge the correlation between adverse events and the test device.

8.1.5.1 Contraindications of video-assisted thoracoscopic surgery

- (1) Active stage of tuberculosis.
- (2) The pulmonary pleura is extensively adherent and the pleural space is atretic.
- (3) The general condition is poor and the infection is not effectively controlled.
- (4) Patients with poor pulmonary function who can not tolerate intraoperative one-lung ventilation anesthesia and pneumonectomy.

8.1.5.2 Preparation

The requirements for anesthesia are the same as or even higher than those for common thoracotomy; and thoracotomy is required if thoracoscopic treatment can not be adopted. Preoperative diagnostic tests should be conducted to further understand the condition of the lesion, know the presence or absence of concomitant diseases and general condition, so as to help determine whether to adopt the thoracoscopic surgery, surgical method and surgical date to make a better decision.

8.1.5.3 Medical history and physical examination

Thorough and accurate disease history collection is the first step in preoperative preparation. In addition to inquiring about the history of major diseases, the past history of the subject should also be detailed. In particular, it is necessary to understand whether there is a history of thoracic tuberculosis, pleurisy, chest trauma, surgery and hydropneumothorax. Since these diseases will alter the normal anatomy of the pleural space, the thoracoscopic surgery therefore is difficult or impossible. In addition, attention should be paid to the health of the patient's cardiovascular and respiratory systems to understand whether intraoperative one-lung ventilation and surgical trauma can be tolerated. Systematic preparation is another important part of preoperative preparation, which is helpful for the diagnosis and detection of concomitant diseases and shall be focused on.

8.1.5.4 Laboratory examination

Preoperative routine blood, urine and stool examination shall be conducted to find abnormalities. Reexamination shall be timely carried out and the causes shall be found. If low hemoglobin indicates anemia and malnutrition; a small number of blood transfusions can be given to make hemoglobin close to normal levels in a short period of time, which improves the safety of surgery. Abnormalities in water electrolyte tests should be corrected before surgery. Liver and kidney function tests can help understand the functional status of liver and kidney organs and surgical tolerance. Since intrathoracic hemorrhage itself is not easy to stop spontaneously, the subject's coagulation function should be examined before surgery. If the subject has a history of coronary heart disease in the near future before surgery, myocardial enzyme examination should be performed. Routine blood glucose tests can detect some asymptomatic patients with type II diabetes.

8.1.5.5 Auxiliary examination

Chest X-ray examination is a quick and effective method to understand the intrathoracic condition. Chest X-ray, chest radiography, tomography and CT can be performed according to the needs of the disease. CT examination should generally be done to provide more information about the disease and intrathoracic structure for researchers and provide a more sufficient basis for preoperative surgical planning. Preoperative ECG examination is helpful to understand the heart condition, and echocardiography and myocardial imaging can also be done when necessary. Ultrasonography can examine pleural effusion and abdominal organs. Preoperative pulmonary function test is very important. In addition to the determination of total pulmonary function, it is also necessary to determine the split pulmonary function when conditions permit, so as to have an objective estimate of the subject's surgical tolerance.

8.1.5.6 Mental preparation

The necessity and possible discomfort of the surgery should be patiently introduced to the subjects to relieve the subjects' ideological concerns and enhance the confidence of surgical treatment, which can help subjects to better cooperate with the thoracoscopic surgery.

8.1.5.7 Preoperative routine preparation

The preparation one day before thoracoscopic surgery and on the day of surgery is the same as the preparation before thoracotomy for similar diseases. If the thoracoscopic Operation time is expected to be short, a catheter may not be indwelled. The blood matching volume is also relatively less, and a small amount of blood or no blood matching can be performed according to specific circumstances. Since the thoracoscopic surgery needs to be timely converted into thoracotomy in case of lesions or complications that can not be treated with thoracoscopy during thoracoscopic surgery, thoracotomy instruments should also be prepared for emergency use.

8.1.5.8 Preoperative respiratory function estimation

Although thoracoscopic surgery has little effect on the pulmonary function of subjects, the pulmonary complications after thoracoscopic surgery mainly stem from the low pulmonary function of subjects. The subjects' preoperative pulmonary function should be known.

8.1.5.9 Anesthesia and surgical positions

General anesthesia is performed using a double-lumen endotracheal tube with unilateral lung strengthening ventilation. The patient is in healthy lateral decubitus position, the operating bed is in "jackknife position" 30 °, the upper limb of the operating side is suspended on the anesthesia head frame, and the investigator operates on the side of the patient. After lobectomy, positive pressure ventilation of the affected lung is performed to fully expand the residual lung to avoid the presence of thoracoscopic illegible localized atelectasis.

8.1.5.10 Surgical incision

Thoracoscopic lobectomy surgical incision, including 1 thoracoscopic cannula incision about 1.5 cm in length, 1~3 operating cannula incisions about 1.5 cm in length, or/and 1 chest wall-assisted small incision 5~7 cm in length.

- (1) Thoracoscopic incision: generally selected between the 8th intercostal space and axillary midline. The choice of incision location varies slightly from subject to subject and from lobe to lobe.
- (2) Operating cannula incision: 1 ~ 2 (sometimes 3) operating cannula incisions are generally used; their positions can be determined after thoracoscopic exploration of the thoracic cavity, based on the principle of facilitating surgical operation. The retractor operation hole is generally selected near the posterior axillary line of the 7th and 8th intercostal spaces.
- (3) Small incision on chest wall: The location of small incision is generally selected between the anterior and posterior axillary line of the 5th intercostal space. In addition, it can be selected according to the surgical needs and resection of different lobes. The selection of small incision should generally follow the principles of proximity to the pulmonary hilum, less chest wall injury, and relatively beautiful scar of the incision.

8.1.5.11 Treatment of hilar vessels and bronchi

- (1) Anatomic lobectomy is shown under thoracoscopic visualization. The investigator lifts the adventitia with an atraumatic endoscopic grasper in one hand and endoscopic dissecting hook or endoscopic scissors in the other hand, divides the adventitia, bluntly separates the vessels, overhangs the traction line, lifts the vessels, ligates the pulmonary vessels with an endoscopic knotter or cuts the pulmonary vessels with an endoscopic vessel suture cutter under direct vision through a small chest wall incision, and the larger vessels can be treated by ligating them at the proximal end of the vessel and then cutting them with a cutting stapler. In order to reduce pulmonary congestion and pulmonary section oozing, the pulmonary artery can be firstly treated and the pulmonary veins can be postconditioned. Bronchial handling is usually performed after pulmonary vascular handling. The peribronchial tissues are lifted with endoscopic grasping forceps, the endoscopic dissecting hook is separated to the distal end of the resected bronchus, the peribronchial lymph nodes are removed, electrocoagulation is performed for hemostasis, the bronchial artery is clamped with endoscopic titanium staples, the bronchus is completely freed, the cutting surface is selected, the endoscopic bronchectomy stapler is inserted, the atretic bronchus is compressed. After the reload is examined as exact and complete, the bronchus is cut. The diseased lung is loaded into the specimen bag and removed through a small incision. Electrocoagulation or argon is adopted for hemostasis. The thoracic cavity is flushed, and the two-lung ventilation is changed. The lung is expanded to examine whether there is air leakage at the bronchial stump. The surgery is ended after the drainage tube is placed.

(2) Lobectomy with massive stapling of the hilum without retraction of the ribs. Thoracoscopy showed that the hilum was sharply and bluntly dissected to expose the lobar vessels and bronchi, the lung grasper lifted the lobe to be removed, the hilum was bypassed with a rubber catheter, the catheter was extracted from a chest wall incision, the thick tip of the catheter was placed on an endoscopic stapler to fire the stapler and staple the hilum. The stapler was then removed, the clip was reloaded, advanced through the incision into the thoracic cavity, the hilum was again jammed distal to the first stapling, and the stapler was fired. The lobe was excised from the distal end of the stapler, packed in a specimen bag, and removed from the small incision. Water was injected into the thoracic cavity to inflate the lung to check whether the stapling of bronchial stump was true or not. If there is no active bleeding and there is no air leakage at bronchial stump, place a closed drainage tube to end the operation.

8.1.5.12 Removal of lung specimens

The specimens shall be carefully removed through the incision by placing them in a specimen bag or sterile glove.

8.1.5.13 Lymphadenectomy

The lymph nodes and surrounding adipose tissue of each group shall be carefully dissected.

8.1.5.14 Chest tube placement and incision closure

Two upper and lower thoracic drainage tubes are placed for upper lobectomy, and lower thoracic drainage tubes are placed for lower lobectomy, both of which are placed through the endoscopic hole or operation hole. The upper thoracic drainage tube could be placed on the top of the chest along the lateral side of the lung; the placement of thoracic drainage tube should be performed under direct thoracoscopic vision. Intercostal nerve closure is performed with 1% bupivacaine in the upper and lower ribs of each incision to reduce postoperative pain. Under direct thoracoscopic assistance, the pleura and intercostal muscles are closed, followed by sequential suture of the chest wall tissue and finally suture of the cannula incision.

8.1.5.15 Precautions

- (1) Lobectomy should be based on the results of various examinations to preserve good lung tissue as much as possible.
- (2) The section should ensure that there is no lesion as far as possible. The tracheal stump should preferably be sutured with absorbable suture and embedded with the surrounding tissues to reduce the occurrence of stump fistula.
- (3) Pulmonary vascular injury should be prevented. Accidental vascular injury during surgery is common in cases of chronic pulmonary suppuration, with severe perivascular adhesions or interlobar fissure hypodifferentiation. In this case, it can be operated under direct vision through a small thoracic incision, the hilum is dissected with conventional thoracotomy surgical instruments, and the pulmonary vessels

are freed. Once major bleeding occurs, the bleeding site should be immediately compressed with gauze pieces and converted to thoracotomy in a timely manner.

- (4) The surgery procedures should be flexibly controlled. The anatomical characteristics of each lobe are different, and the operation sequence of thoracoscopic anatomical lobectomy is also different. Left lower lobe or right middle lobe and lower lobe resection is basically the same as conventional thoracotomy. The artery is firstly processed, the vein is then processed, and the bronchus are finally processed; for right upper lobe pneumonectomy, the posterior segmental artery of the upper lobe is firstly dissected. After cutting off, the upper lobe lung is pulled forward to process the upper lobe bronchus, the anterior branches of the upper lobe artery are then processed, and the upper lobe branches of the upper pulmonary vein are finally processed; for left upper pneumonectomy, the lingual artery through the interlobar fissure is generally separated and cut, the anterior and posterior branch arteries upwards are separated and processed, and then the upper pulmonary vein and bronchus are processed, respectively. The basic principle of VATS pneumonectomy is basically the same as that of lobectomy, but the pulmonary vessels and bronchi are thick. Vascular treatment advocates proximal ligation and then cutting stapler.
- (5) Anti-tuberculosis drugs must be continuously taken to for at least half a year after surgery.

Note: The above treatment principles are only for the reference of the investigators, and the investigators can deal with them according to the actual situation.

8.1.6 Selection of Subjects

8.1.6.1 Clinical trial indications

At present, the Disposable Powered Articulating Endoscopic Linear Cutter Stapler can be used for "a variety of open or minimally invasive surgeries in general surgery, obstetrics and gynecology, urology, thoracic surgery and pediatrics. They can be used together with staple line or tissue supporting materials. They also can be used for transection and excision of liver parenchyma (hepatic vascular system and bile duct structure), pancreas, kidney and spleen.

However, according to the *Guidelines for Technical Review of Endoscopic Stapler Registration*, "the surgical site and/or tissue type and thickness selected for clinical trial shall cover the intended clinical applicable scope of the product under application", and the cutting and anastomosis of lung tissues are common uses of the device. Therefore, thoracic surgery is selected for the clinical trial.

At the same time, clinicians will use linear staplers when performing VAT anatomical lobectomy (segmentectomy) according to clinical research. Therefore, the indication of this clinical trial is: "**anatomic lobectomy (segmentectomy).**"

8.1.6.2 Inclusion criteria

Subjects meeting all of the following conditions only can be included:

- (1) Subjects are aged 18-70 years (inclusive), male or female;
- (2) Subjects are scheduled to accept thoracoscopic anatomical lobectomy (segmentectomy);
- (3) Subjects or their guardians are able to understand the purpose of the study, show adequate compliance with the trial protocol, and sign an ICF.

8.1.6.3 Exclusion criteria

Subjects meeting any of the following conditions can not be included:

- (1) Subjects are with video-assisted thoracoscopic surgery-related contraindications;
- (2) Subjects' preoperative imaging findings indicate severe pleural adhesions and significantly calcified hilar lymph nodes;
- (3) Subjects with platelet (PLT) $< 60 \times 10^9/L$ or international normalized ratio (INR) > 1.5 ;
- (4) Subjects with forced expiratory volume in one second (FEV1)/predicted value $\leq 50\%$, or forced expiratory volume in one second (FEV1)/forced vital capacity (FVC) $\leq 60\%$;
- (5) Subjects with ejection fraction $\leq 50\%$;
- (6) Subjects with major organ failure or other serious diseases (including clinically relevant cardiovascular disease or myocardial infarction within 12 months prior to enrollment; severe neurological or psychiatric history; preoperative severe infection that must be controlled by medication; active disseminated intravascular coagulation; high risk of thrombosis);
- (7) Subjects who are pregnant or lactating women;
- (8) Subjects participating in other drug or device clinical trials within 1 month before the trial;
- (9) Other situations not suitable for inclusion judged by the investigators.

8.1.6.4 Study termination criteria

Study termination refers to that the clinical study has not been completed according to the protocol, and all clinical studies are stopped halfway. The purpose of terminating the clinical study is to protect the rights and interests of the subjects and ensure the quality of clinical study. Except for Item 6 below, the sponsor and the investigators shall jointly negotiate whether to terminate the trial:

- (1) If serious safety problems occur during the clinical study, the clinical study should be terminated in a timely manner;
- (2) If the product is found to have no clinical value during the clinical study, the clinical study shall be terminated;
- (3) In the process of clinical study, it is found that the clinical study protocol has major errors and it is difficult to evaluate the product effect; or the protocol has serious deviation in the implementation, in

which it is difficult to evaluate the product effect if the study is continuously to be conducted, , the clinical study shall be terminated;

- (4) The sponsor requests to terminate the study (for reasons, e.g., financial, etc.);
- (5) The clinical trial institutions and investigators fail to comply with the laws and regulations related to clinical trial and clinical trial protocol, and fail to make correction after being pointed out, resulting in serious situations or failure to change.
- (6) NMPA orders to terminate the clinical study for some reason(s).

After the study is terminated, the investigators should continue to provide appropriate treatment for the subjects from the perspective of protecting the rights and interests of the subjects, and inform the subjects of the treatment received during the trial and related treatment in detail. The subject data continues to be acceptable for evaluating the safety of the product.

The investigators should timely feed back the termination of the subject study to the sponsor.

8.1.6.5 Enrollment time

It is estimated that the first subject will be enrolled in August 2018. The total enrollment time is 6 months, and the trial end time for the last subject will be March 2019.

8.1.6.6 Expected duration of clinical trial and the reasons for determination

The expected enrollment time for the first subject is August 2018. The enrollment time for subjects is 6 months. The time for the last subject to complete the trial is March 2019. The time for statistical analysis and preparation of study summary report after the completion of clinical trial is 2 months, and the total clinical study period is 8 months.

8.1.6.7 Expected duration of participation of each subject

The expected screening period from signing the ICF to enrollment for each subject is 7 days, and the postoperative follow-up period is about 30 days. The total subject observation period is 37 days. If a subject experiences an adverse event during the clinical trial, the outcome of the adverse event needs to be observed.

8.1.6.8 Number of subjects required for clinical trial

The total sample size is 164. See Section 9.2 for details.

8.1.6.9 Drop-out criteria and treatment

Definition of drop-out: All the subjects who have signed the ICF and passed screening to enter the clinical trial, no matter when or when they withdraw from the study, will be called drop-out cases as long as they have not completed the observation period specified in the protocol.

Drop-out criteria:

- (1) The subject is unwilling or impossible to continue the clinical trial for any reason, and the clinical trial is terminated by putting forward a request for withdrawal from the clinical trial to the investigator;

- (2) The subject is lost to follow-up due to no treatment or examination although the subject does not explicitly withdraw from the clinical trial;
- (3) If the subject experiences a serious adverse event (SAE), the investigator should withdraw the subject from the clinical trial;
- (4) Poor subject compliance.

Treatment of dropout cases: For all the dropout cases, it is necessary to fill in the summary form of completion of clinical trial in the case report form and analyze the reasons for dropout. When a subject dropped out, the investigator should contact the subject as far as possible, ask the reason and record the last clinical symptom, and complete the evaluation items that could be completed. When withdrawing from the study due to adverse events, the investigator should take appropriate treatment measures according to the actual situation of the subject.

When a subject or his/her guardian actively requests to withdraw from the trial, the subject or his/her guardian should be asked that whether the subject's previous data can be continuously used in the clinical trial, and whether the results can be recorded in the original records. When no enquiry is performed, the data are available by default.

8.1.6.10 Supplementary procedures for subjects

During the clinical trial, if the number of subjects fails to meet the statistical requirements due to dropout for any reason, or the investigator considers it necessary to supplement subjects, the reasons shall be analyzed firstly, and then the investigator, sponsor and statistician shall discuss whether to supplement or not. After determining the number of supplemented cases, the information shall be submitted to relevant departments for review.

8.1.6.11 Study completion criteria

Study completion refers to that the clinical study has included sufficient subjects according to the requirements of clinical protocol and completed the treatment, observation and evaluation for all subjects, and the collected data can be used to evaluate the safety and effectiveness of the study product.

8.1.7 Effectiveness evaluation method

8.1.7.1 Main effectiveness evaluation indicators, method and time for evaluation, recording and analysis

(1) Success rate of device cutting and anastomosis

Justification for selection: Through the mechanical transmission device, the stapler places two or multiple rows of staples prearranged in parallel and staggered rows in the component into the tissues that are aligned and need to be anastomosed. After passing through the tissues, the staples are blocked by the anvil jaw and

bend inward, forming a B-shaped staggered arrangement to anastomose the tissues together. Whether the stapler can successfully cut and staple the lung tissue is very critical to the video-assisted thoracoscopic surgery. If the cutting and stapling of lung tissue fails, the subject may need to be converted into thoracotomy, which will increase the risk of the subject. Therefore, in this clinical trial, the "success rate of device cutting and anastomosis" is selected as the primary endpoint, and the effectiveness of this product is evaluated by calculating the number of subjects required to be converted into thoracotomy due to the stapler.

Evaluation time: During the surgery.

Evaluation site: Tissues receiving linear stapler for pulmonary anastomosis.

Evaluation method: The transection and anastomosis of lung tissues are performed during the operation. The device is withdrawn after successful triggering. The cutting staple line is carefully checked for integrity, air leakage and bleeding. If all staple lines are complete, without air leakage or bleeding, it is judged that the cutting and anastomosis of the device is successful. If there is persistent air leakage and bleeding at the anastomosis site, and **conversion to thoracotomy is required**, it is judged as device cutting and anastomosis failure.

Formula: Device cutting and anastomosis success rate = number of subjects with successful anastomosis in each group/total number of subjects in this group × 100%

Precautions: ① During VAT anatomical lobectomy (segmentectomy), the subject may be converted into thoracotomy due to his complicated condition. At this time, the investigators need to determine whether the subject's conversion to thoracotomy is caused by the stapler. If it is caused by the stapler, it is judged as device cutting and anastomosis failure; ② During intraoperative lung tissue cutting and anastomosis, the investigators need to perform video shooting for the cutting and anastomosis process, so as to perform necessary traceability for the original data. The video should at least contain the whole process of lung tissue cutting and anastomosis.

8.1.7.2 Secondary effectiveness evaluation indicators, method and time for evaluation, recording and analysis

(1) Operation time

Justification for selection: The main purpose of stapler is to reduce the surgeon's operation intensity and save operation time. Therefore, the operation time is selected as one of the indicators to evaluate the product performance. Meanwhile, with reference to the *Guidelines for Technical Review of Endoscopic Stapler Registration* issued by NMPA (original CFDA) on March 13, 2017, the "operation time" is selected as the secondary evaluation indicator for such products.

Evaluation time: During the surgery.

Mode of evaluation: The investigators record the time taken from thoracotomy to chest closure to the

nearest minute.

(2) Anastomosis time

Justification for selection: The main purpose of stapler is to reduce the surgeon's operation intensity and save operation time. Therefore, the stapling time is selected as one of the indicators to evaluate the product performance. Meanwhile, with reference to the *Guidelines for Technical Review of Endoscopic Stapler Registration* issued by NMPA (original CFDA) on March 13, 2017, the "stapling time" is selected as the secondary evaluation indicator for such products.

Evaluation time: During the surgery.

Evaluation method: The investigators record the time from the subject **first** cutting and anastomosis of lung tissues to the time when the device is inserted into the thoracic cavity to the completion of the clamping, cutting and withdrawal from the thoracic cavity. The time is accurate to second (s).

(3) Incidence of chest tube drainage greater than 7 days

Justification for selection: Staple line healing is an important indicator for effectiveness and safety evaluation of such products. Therefore, the healing rate of staple line is selected as one of the indicators to evaluate the product performance. However, since the healing of staple line is not easy to be observed, and postoperative chest tube drainage is a routine operation of VAT lobectomy (segmentectomy), the removal of chest tube drainage indicates that the staple line has gradually healed. Meanwhile, the *Guidelines for Technical Review of Endoscopic Stapler Registration* issued by NMPA (original CFDA) on March 13, 2017 is referred to. Eventually, the "incidence of chest tube drainage greater than 7 days" is used as an alternative indicator for evaluation.

Evaluation time: After surgery.

Mode of evaluation: The investigators not involving in the subject's treatment judge that whether the subject could remove the chest drainage tube, and record the time when the subject finally removes the chest drainage tube. The time is accurate to day (natural day).

Calculation formula: Incidence of chest tube drainage for more than 7 days = number of subjects in each group with chest tube drainage for more than 7 days/total number of subjects in this group × 100%.

Precautions: ① In order to avoid the evaluation bias caused by subjective reasons, whether the subjects can remove the chest drainage tube should be independently judged by the investigators who do not participate in the treatment of subjects; ② For the subjects with chest tube drainage for more than 7 days, the possible reasons for chest tube drainage for more than 7 days should be analyzed to determine whether it is related to the device.

(4) Conversion rate to thoracotomy

Justification for selection: The need for conversion to thoracotomy due to failed thoracoscopic surgery

increases the risk to the subject. Therefore, the "conversion rate to thoracotomy" is selected as one of the effectiveness indicators of the product. At the same time, with reference to the *Guidelines for Technical Review of Endoscopic Stapler Registration* issued by NMPA (original CFDA) on March 13, 2017, the "conversion rate to thoracotomy" is selected as an effectiveness evaluation indicator for such products.

Evaluation time: During the surgery.

Evaluation method: The reasons for conversion to thoracotomy during thoracoscopic surgery are truthfully recorded. The common cases for conversion to thoracotomy are pulmonary artery bleeding, pulmonary vein bleeding, lymph node metastasis or lymph node calcification and adhesion, pulmonary artery anatomical deformity, and linear cutter stapler failure.

Formula: Rate of conversion to thoracotomy = number of eligible subjects in each group/total number of subjects in this group $\times 100\%$

(5) Number of staple line repair needles

Justification for selection: When the stapler is used for pulmonary tissue anastomosis, the situation that the staple line needs to be repaired may occur. When the number of repair needles becomes larger, it may increase the risk to the subject. Therefore, the "number of stitches repaired with staple line" is selected as one of the indicators to evaluate the effectiveness of the product. At the same time, with reference to the *Guidelines for Technical Review of Endoscopic Stapler Registration* issued by NMPA (original CFDA) on March 13, 2017, the "number of stitches repaired by staple line" is selected as an effectiveness evaluation indicator for such products.

Evaluation time: During the surgery.

Evaluation method: The number of staple line repair needles of each subject during operation is recorded.

Precautions: ① Anastomotic repair is required when the anastomotic ring is incomplete, leaking and bleeding. Anastomotic repair may also be necessary when there is a dehiscence of the staple line during surgery. Sometimes the investigators may perform anastomotic repair in order to strengthen the anastomotic site. Therefore, the reasons for repair should be recorded; ② If the investigators use other devices for repair, it is necessary to record the product name and the quantity used, and finally replace the number of repair needles with the quantity used.

8.1.8 Safety evaluation method

(1) Incidence of staple line dehiscence

Justification for selection: When stapler is used for pulmonary tissue anastomosis, it is necessary to select appropriate reload. Dehiscence of the staple line may increase the risk to the subject. Therefore, the "incidence of staple line dehiscence" is selected as one of the safety indicators of the product under evaluation. At the same time, with reference to the *Guidelines for Technical Review of Endoscopic Stapler*

Registration issued by NMPA (original CFDA) on March 13, 2017, the "incidence of staple line dehiscence" is selected as a safety evaluation indicator for such products.

Evaluation time: During the surgery.

Method of evaluation: Staple line dehiscence due to "staple popping" during surgery is recorded.

Formula: Incidence rate of staple line dehiscence = number of eligible subjects in each group/total number of subjects in the group $\times 100\%$

(2) Incidence of staple line bleeding

Justification for Selection: Staple line bleeding is a common complication of stapling procedures in lung tissue and increases the risk to the subject. Therefore, "the incidence of staple line bleeding" is selected as one of the indicators to evaluate the product safety. Meanwhile, with reference to the *Guidelines for Technical Review of Endoscopic Stapler Registration* issued by NMPA (original CFDA) on March 13, 2017, the "incidence rate of staple line bleeding" is selected as a safety evaluation indicator for such products.

Evaluation time: from the end of surgery to the end of follow-up.

Evaluation mode: Respiratory tract bleeding occurs after operation, and it is determined whether the bleeding is the source of anastomotic line by observing thoracic drainage tube, thoracoscopy and angiography. If the source is not clear, it is called respiratory tract bleeding.

Formula: Incidence rate of staple line bleeding = number of eligible subjects in each group/total number of subjects in the group $\times 100\%$

(3) Incidence of air leakage

Justification for selection: Air leakage is a common complication of lung tissue anastomosis, which will increase the risk of subjects. Therefore, the "occurrence rate of air leakage" is selected as one of the safety indicators for product evaluation. At the same time, with reference to the *Guidelines for Technical Review of Endoscopic Stapler Registration* issued by China Food and Drug Administration on March 13, 2017, the "incidence rate of air leakage" is selected as a safety evaluation indicator for such products.

Evaluation time: During the surgery.

Evaluation mode: After the completion of cutting and anastomosis of lung tissues, the investigators drip normal saline to the staple line to see whether there are bubbles, which is used as the basis for judging air leakage.

Formula: Incidence of air leakage = number of subjects with air leakage/total number of subjects in the group $\times 100\%$

(4) Incidence of Staple Line Infection

Justification for Selection: Staple line infection is a common complication of lung tissue stapling procedures and poses an increased risk to the subject. Therefore, "the incidence of staple line infection" is

selected as one of the indicators to evaluate the product safety.

Evaluation time: from the end of surgery to the end of follow-up.

Mode of evaluation: The clinician judges whether the subject has a staple line infection based on the symptoms, signs, endoscopy, angiography, etc.

Formula: Incidence rate of staple line infection = eligible subjects in each group/total number of subjects in the group $\times 100\%$

(5) Incidence of SAEs

Evaluation time: During the whole clinical trial process.

Evaluation method: Serious adverse event refers to events occur in the process of clinical trial that causes death or serious deterioration of health, including fatal disease or injury, permanent damage to the body structure or body function, requiring hospitalization or prolonging the hospitalization, requiring medical or surgical intervention to avoid the permanent damage to the body structure or body function.

Formula: SAE incidence = number of subjects with SAE in this group/number of all subjects in this group $\times 100\%$.

Precautions: The adverse events related to the staple line during clinical trial shall be considered as serious adverse events. Since "staple line repair" may be associated with "staple line dehiscence", the "staple line repair" is not repeatedly considered a serious adverse event.

(6) Incidence of AE (adverse event)

Evaluation time: During the whole clinical trial process.

Formula: Incidence of adverse events = number of subjects with adverse events in this group/number of all subjects in this group $\times 100\%$.

Definition: An adverse event is any untoward medical occurrence in a clinical investigation, whether or not related to the device. It is necessary to analyze the postoperative vital signs and laboratory tests (blood routine, biochemical indicators, etc.).

(7) Incidence of device defects

Definition: An test device defect that occurs during the course of a clinical trial. "Staple pop" due to human factors (e.g., incorrect reload selection) is not a defect of the device.

Measurement time: During the whole process of clinical trial.

Calculation formula: Incidence of device defect = number of device defects/total number of devices $\times 100\%$.

8.1.9 Operational performance evaluation

(1) Stapler operating performance evaluation

Measurement time: During the stapler operation process.

Evaluation method: ① Whether the single-handed operation of stapler is smooth; ② whether the stapler can fire normally; ③ whether the back-off is smooth after the completion of anastomosis, whether there is tissue entrapment, whether there is resistance; ④ whether the stapler can fire continuously for many times is reliable; ⑤ whether the stapler can completely cut off the tissue; ⑥ whether the battery supply is stable; ⑦ whether the electric cutting is smooth.

8.2 Test Process

8.2.1 Test Flow Chart

Visit Time Visit Item	Visit 1	Visit 2	Visit 3	Visit 4
	Before surgery (Day -7 ~ 0)	During surgery (Day 0)	After surgery (Day 7 ± 2 days)	After surgery (Day 30 ± 5 days)
Informed consent	▲			
Subject screening	▲			
Demographic data	▲			
History of present illness	▲			
Past medical history	▲			
Vital signs	▲	▲	▲	▲
Blood routine	▲		▲	
Coagulation routine	▲		▲	
Blood biochemistry	▲		▲	
Cardiac examination	▲			
Pulmonary function test	▲			
Arterial blood gas analysis	▲			
Preoperative imaging	▲			
Blood pregnancy test ¹	▲			
Randomization	▲			
Intraoperative observation record		▲		
Use of device		▲		
Postoperative imaging ²			▲	▲
Postoperative follow-up observation			▲	▲
Concomitant medication	▲	▲	▲	▲
Adverse events	▲	▲	▲	▲
Device defect		▲		

¹ Only for women of childbearing potential.

² In case of anastomotic leakage and bleeding, the investigators may perform relevant examinations for the subjects according to actual needs.

8.2.2 Enrollment

The investigators screened the subjects according to the inclusion and exclusion criteria after obtaining their informed consent; on the day of conducting total-thoracoscopic anatomic lobectomy (segmentectomy), the subjects were randomly divided into the test group and the control group according to the random assignment instructions.

8.2.3 Randomization

The trial was stratified by clinical trial institutes with the randomization method of stratified blocks. The seed number and block length were determined with SAS9.1 statistical software programming. The subjects were divided into the test group and the control group in a ratio of 1:1. The random grouping arrangement of at least 164 subjects was generated, and the corresponding sequentially numbered, opaque and sealed randomization envelopes were made. The randomization envelopes were kept and allocated by relevant personnel who were not involved in the screening and treatment of subjects and authorized by the investigator.

The random number is a 4-digit sequential number consisting of "center number" and "serial number". The first digit is the number of the clinical trial institution, and the last 3 digits are the "serial number" assigned by randomization. When the serial number is less than 3 digits, add "0" to make up 3 digits. If the randomization number "1001" is the first randomization number of the clinical trial institute with the institute number of 01, it corresponds to the first subject randomized at the institute.

When a subject was determined to be enrolled, the investigators or other authorized designee requested randomization assignment from the randomization envelope manager. The randomization envelope keeper took out the randomization envelope in the order of random number, filled in the screening number and initials of the subject to be assigned to the designated place on the back of the envelope, and signed the name and date to complete the allocation. The investigators or other authorized designee then intervened based on the results of the randomization assignment (i.e., test groups).

8.2.4 Blinding method

The investigators could not be blinded due to the difference in appearance between the test device and the control device, this trial therefore adopted an incomplete blinding method, and the subjects were not informed of the randomization results (i.e., the subjects did not know their own treatment assignment).

8.2.5 Inspection items

8.2.5.1 Medical history and vital signs

- (1) Demographic data: gender, date of birth;
- (2) History of present illness: admission time, cause;
- (3) Past medical history: allergy history, dysfunction history of vital organs (heart, lung, liver and kidney), surgery history, and previous medication history ³ (within 1 month);
- (4) Vital signs: height, weight, blood pressure, pulse, temperature, respiration, BMI ⁴;

8.2.5.2 Laboratory examinations ⁵

- (1) Blood routine: hemoglobin (Hgb), red blood cell count (RBC), white blood cell count (WBC), neutrophil count (NEUT), neutrophil ratio (NEUT%), platelet count (PLT);
- (2) Coagulation blood routine: prothrombin time (PT), activated partial thromboplastin time (APTT), international normalized ratio (INR);
- (3) Cardiac examination: ECG, cardiac ejection fraction;
- (4) Blood biochemistry: aspartate aminotransferase (AST), alanine aminotransferase (ALT), serum total bilirubin (STB), serum albumin (ALB), serum creatinine (Scr), blood urea nitrogen (BUN);
- (5) Pulmonary function test: forced expiratory volume in one second (FEV1)/predicted value, forced expiratory volume in one second (FEV1)/forced vital capacity (FVC);
- (6) Arterial blood gas analysis: pH, oxygen partial pressure, carbon dioxide partial pressure, oxygen saturation (SpO₂);
- (7) Blood pregnancy test: applicable to women of childbearing age (18 ~ 50 years old);

8.2.5.3 Preoperative imaging ⁶

- (1) CT/MRI examination: whether there is pleural adhesion and is to what extent, whether there is hilar significantly calcified lymph node;

8.2.5.4 Postoperative imaging examination ⁷

- (1) Angiography: whether there is staple line bleeding;
- (2) Thoracoscopy: whether there is air leakage, whether there is staple line bleeding;

8.2.5.5 Other examinations

- (1) Intraoperative observation items: operation starting time, operation ending time (24 h clock h/min) ⁸, anastomosis starting time, anastomosis ending time (s), resected lobe or segment ⁹, intraoperative blood

³ Mainly record drugs that affect healing (e.g., antineoplastic agents, corticosteroids) or coagulation.

⁴ Body mass index (BMI) = weight (kg) ÷ height² (m²)

⁵ For laboratory examinations, the reports/results issued by our hospital within 7 days prior to surgery can be accepted.

⁶ For imaging examination, the reports/results issued by our hospital within 1 month prior to surgery can be accepted.

⁷ In case of air leakage and bleeding at the staple line, the investigators may decide whether the subject will receive relevant examinations according to actual needs.

⁸ The duration from the time of thoracotomy to chest closure is recorded.

⁹ Right upper lobe, right middle lobe, right lower lobe, left upper lobe and left lower lobe.

- loss, whether the transection and anastomosis are successful, whether the cutting staple line is complete, whether there is leakage in the cutting staple line, whether there is bleeding in the cutting staple line, whether there is air leakage in the staple line, whether there is conversion to thoracotomy ¹⁰, number of stitches for staple line repair ¹⁰, whether there is staple line dehiscence ¹⁰, whether there are other complications ¹⁰, evaluation on the stapler operating performance;
- (2) Use of device ¹¹: models and specifications of electric-endoscopic stapler, reload models and specifications, device defect;
 - (3) Postoperative follow-up observation: whether there is anastomotic line bleeding ¹², whether there are other complications ¹², duration of chest tube drainage (days), whether there is anastomotic infection;
 - (4) Routine observation items: adverse events/serious adverse events, concomitant and combined medication, protocol deviations.

8.2.6 Device using specification

According to the results of random allocation, the investigators apply the study medical device to the subjects who enter the test group, and the control medical device to the subjects who enter the control group.

8.2.7 Visit schedule

Subjects should receive 4 visits before the end of the study: before surgery (Day -7 ~ 0), during surgery (Day 0), on Day 7 ± 2 after surgery, and on Day 30 ± 5 after surgery.

8.2.7.1 1st visit 1/before surgery (Day -7 ~ 0)

- (1) The ICF form shall be signed;
- (2) Screening subjects according to the inclusion and exclusion criteria;
- (3) Randomized grouping;
- (4) Recording demographic data, current medical history and past medical history are recorded;
- (5) Examining vital signs;
- (6) Performing preoperative imaging examination (CT/MRI examination);
- (7) Performing laboratory examinations (blood routine, coagulation routine, blood biochemistry, heart examination, pulmonary function test, arterial blood gas analysis, blood pregnancy test) ;

¹⁰ If yes, please specify the reason and whether it is related to the stapler. If the investigators use other devices for repair, it is necessary to record the product name and the quantity used, and finally replace the number of repair needles with the quantity used.

¹¹ Truthfully record the model and amount used in the clinical trial.

¹² If corresponding symptoms occur, the reasons for confirming the event should also be briefly described, so as to facilitate the review. The diagnosis is confirmed by thoracoscopy or second surgical exploration.

(8) Recording concomitant and concomitant medications;

(9) Recording protocol deviations;

8.2.7.2 2nd visit 2/during surgery (Day 0)

(1) Examining vital signs;

(2) Intraoperative observation records;

(3) Use of device;

(4) Recording concomitant medication/treatment;

(5) Monitoring of adverse events and serious adverse events;

(6) Recording device defects;

(7) Recording protocol deviations;

8.2.7.3 3rd visit 3/after surgery (Day 7 ± 2 days)

(1) Examining vital signs;

(2) Conducting laboratory examinations (blood routine, coagulation routine, blood biochemistry);

(3) Conducting postoperative imaging examinations (angiography, thoracoscopy) ⁵;

(4) Postoperative follow-up observation;

(5) Recording concomitant and combined medications;

(6) Recording adverse events and serious adverse events;

Recording protocol deviations;

8.2.7.4 4th visit/after surgery (Day 30 ± 5 days)

(1) Examine vital signs;

(2) Conducting postoperative imaging examinations (angiography, thoracoscopy) ⁵;

(3) Postoperative follow-up observation;

(4) Recording concomitant and combined medications;

(5) Monitoring adverse events and serious adverse events;

(6) Recording protocol deviations;

8.2.8 Concomitant medication/treatment management

If a subject takes medication or therapeutic measures that caused the AE for various reasons during the study, such situation shall be recorded in the case report form (CRF) by investigators for reference.

(1) Prohibited therapies: Products with similar efficacy to the study medical device and the control device are prohibited for use by all subjects during the clinical trial.

(2) Allowable treatment: The treatment with other drugs or therapeutic measures required due to underlying diseases, concomitant diseases or adverse events should be recorded in the original medical records in

detail, including the name, dosage, usage of the drug and the name, usage and frequency of the therapeutic measures.

- (3) The generic name, type, dose, frequency, route of administration and days of administration of other drugs used by the subjects as well as the name, method and frequency of treatment should be recorded in detail. The information is used to evaluate the possible bias on clinical trial results.

8.2.9 Case report form completion

Case Report Form (CRF) is a document designed according to the clinical protocol to record the data of each subject during the trial. The correct completion of case report form will help to make a correct judgment on the trial results. Each investigator should complete the case report form carefully and pay attention to the following points:

- (1) Required fields on the CRF must not be left blank;
- (2) When completing the CRF, be sure to make the unit of test results consistent with the unit filled in;
- (3) All data on the CRF should be derived from original documents and be consistent with the original documents;
- (4) The investigator should complete and sign the CRF within 2 weeks after each subject completes relevant examinations, which is for overall monitoring by the clinical research associate;
- (5) If there is any uncertainty in the filling process, please communicate with the sponsor in a timely manner.

8.3 Monitoring Plan

The sponsor shall, in accordance with the requirements of *Good Clinical Practice for Medical Devices*, select the qualified monitor to perform monitoring responsibilities for all the participating institutions. The monitoring frequency of each participating institution will be adjusted according to the speed of subject recruitment and subject number in each institution participating in the study. The monitor will complete the following responsibilities:

- (1) Confirming that the clinical trial institution has been equipped with appropriate conditions before the trial, including qualified staffing and training, complete laboratory equipment, good working condition, expected to have a sufficient number of subjects, and investigators familiar with the trial requirements;
- (2) Monitoring whether the clinical trial institutions and investigators abide by relevant regulations, this specification and clinical trial protocol before, during and after the trial;
- (3) Confirming that each subject has signed the ICF before participating in the clinical trial, and that the inclusion of subjects as well as the progress of the trial are well known; monitoring and recording the failure follow up, non-conducted trial, non-conducted examination by the investigators, and whether the errors and omissions are corrected; for the revised ICF, confirming that the clinical trial process is not

completed and the affected subjects have re-signed the ICF;

- (4) Confirming that all the CRFs are correctly filled in and consistent with the original data; that all the errors or omissions have been corrected or explained, which have been signed and dated by the investigators; and that the disease type, total number of cases and gender, age and therapeutic effect in each trial have been confirmed and recorded;
- (5) Confirming that the subject withdrawing from the clinical trial or non-compliant with the requirements set forth in the ICF is recorded and discussed with the investigators;
- (6) Confirming that all adverse events, complications and other device defects are recorded, and serious adverse events and device defects that may cause serious adverse events are reported and recorded within the specified time;
- (7) Monitoring the supply, use, maintenance, transportation, reception, storage, distribution, processing and recovery of investigational medical device samples;
- (8) Supervising the regular maintenance and calibration of relevant equipment during the clinical trial;
- (9) Ensuring that all the documents related to the clinical trial received by the investigators are of the latest version;
- (10) Submitting a written report to the sponsor after each monitoring, which shall include the name of the monitor, date of monitoring, time of monitoring, place of monitoring, content of monitoring, name of the investigator, completion status of the project, existing problems, conclusions, and corrections to errors and omissions, etc.

9. Statistical Considerations

9.1 Statistical Design, Methods and Analysis Procedure

This study is a prospective, multi-center, stratified block randomization, incomplete blinding, parallel positive control and non-inferiority test. All statistical analysis methods are detailed in the statistical analysis plan. The first draft of the statistical analysis plan will be formed after the trial protocol and case report forms are finalized and will be finalized before database lock.

9.2 Sample Size Calculation

9.2.1 Total sample size

164 subjects are planned to be enrolled in this clinical trial. The sample size determination process is as follows:

This clinical study is a parallel positive control clinical validation design. The "device cutting and anastomosis success rate" is adopted as the basis for sample size estimation. Based on literature review ^[2], a

total of 160 patients were included in the study. All patients underwent thoracoscopic surgery with Johnson & Johnson™ Echelon linear cutter stapler. During the surgery, 8 patients were converted into thoracotomy. Only 1 case of conversion to thoracotomy was resulted from linear cutter stapling failure. The other reasons were: intraoperative pulmonary artery bleeding in 3 cases, pulmonary vein bleeding in 1 case, lymph node metastasis or lymph node calcification and adhesion in 2 cases, and pulmonary artery anatomical malformation in 1 case. Therefore, when the similar device was used for VAT anatomical lobectomy (segmentectomy), the success rate of device cutting and anastomosis was up to 99.38% (1/160). The clinical investigators participating in this trial determined that the "device cutting and anastomosis success rate" of this clinical trial was conservatively estimated to be 97%, with a clinical non-inferiority margin of -10% based on the relevant literatures and considering the clinical practice. The parameters were set as follows: $\alpha = 0.025$ (one-sided), $\beta = 0.2$, $\delta = -10\%$, and the ratio of test group and the control group was 1:1. The formula of sample size estimation for non-inferiority clinical validation is adopted, and the sample size required in each group is at least 65 cases. The calculation formula and correction formula are as follows:

$$n_3 = \frac{(\pi_A(1 - \pi_A) + \pi_B(1 - \pi_B))(Z_{1-\beta} + Z_{1-\alpha})^2}{((\pi_A - \pi_B) - d)^2}$$

$$n_4 = \frac{n_3 \left\{ 1 + [1 + 4/(\delta n_3)]^{\frac{1}{2}} \right\}^2}{4}$$

Considering that the number of subjects who may drop out or be lost to follow-up during clinical validation is about 20%, the sample size of each group is expanded to 82 cases, and the total sample size required for the two groups is 164 cases.

9.2.2 Number of cases for clinical trial of each disease and reasons for determination

For the patients accepting anatomical lobectomy (segmentectomy) in this trial, the number of cases for each disease in the clinical trial is not specified.

9.2.3 Minimum and maximum number of subjects in each clinical trial institution and reason

A total of 164 subjects are planned to be enrolled in this clinical trial, and a total of 5 clinical trial institutions are planned to participate in this trial. In the trial, the sample balance in each site should be maintained as far as possible. In combination with the total sample size and the number of participating clinical trial institutions, considering different number of disease sources in each center and possible difficulties in individual centers, each clinical trial institution should recruit at least 16 subjects but no more than 60 subjects.

9.3 Significance Level and Power of Clinical Trial

Significance level α (one-sided) = 0.025 with 80% power.

9.4 Expected Dropout Rate

A 20% dropout rate is expected during sample size estimation.

This drop-out rate refers to the proportion of subjects who finally can not be included in the primary analysis. Drop-out subjects refer to those confirmed by the principal investigator with serious protocol violation (affecting the main evaluation/efficacy evaluation). The following conditions may be included, but not limited to: the subject violates the inclusion and exclusion criteria; the subject does not receive any study device; the subject receives the concomitant medication/device prohibited by the protocol; the data of main evaluation indicators are missing; the subject fails to obtain any data after randomization; the subject fails to complete the follow-up visit according to the regulations.

9.5 Qualification/Disqualification Criteria for Clinical Trial Results

The non-inferiority hypothesis is tenable if the lower limit of 95% confidence interval of the difference in the primary effectiveness evaluation indicator "device cutting and anastomosis success rate" between the test group and the control group is greater than the non-inferiority critical value (-10%).

The secondary effectiveness evaluation indicators and safety evaluation indicators are tested by difference test, and the significance level is set at 0.05. If the test result of the difference in any secondary evaluation indicator and safety evaluation indicator between the test group and the control group is $p > 0.05$, then there is no significant difference in this indicator between the test group and the control group.

9.6 Criteria and Reasons for Terminating Trial Based on Statistical Reasons

No interim analysis will be performed for this study, therefore, no standard for terminating the trial for statistical reasons is developed.

9.7 Statistical Methods of Data, Together with Processing Methods of Missing, Unused or Wrong Data and Unreasonable Data

9.7.1 Statistical method for data

Statistical analysis is performed using the SAS 9.1 software.

Statistical description on all data, including demographic data, baseline data, all efficacy indicators and all safety data shall be performed. Measurement data will be statistically described using mean, standard deviation, minimum, maximum, median, 25th percentile and 75th percentile; enumeration data will be statistically described using frequency and percentage.

The primary effectiveness evaluation indicator (device cutting and anastomosis success rate) is evaluated by the confidence interval method.

The secondary evaluation indicators are evaluated by the difference test. Measurement data are compared using two independent samples t-test. Enumeration data are compared using Pearson's χ^2 test or Fisher's

exact test.

The safety evaluation indicators are evaluated by the difference test. Measurement data are compared using two independent samples t-test. Enumeration data are compared using Pearson's χ^2 test or Fisher's exact test. The adverse events and device defects occur in this trial shall be shown in list.

The stapler operational performance is evaluated by the difference test, with the Pearson's χ^2 test or Fisher's exact test.

The central effect is considered in the analysis process. The general linear model is used for measurement data, and CMH is used for enumeration data.

9.7.2 Missing, unused, or incorrect data (including dropouts and withdrawals) and unreasonable data

All missing, unused or incorrect data (including dropouts and withdrawals) and unreasonable data will be discussed and finally determined by the investigators and biostatistician in the data review stage. The basic statistical principles for the processing of these data are as follows:

- (1) The details of each drop-out case shall be described, and the drop-out between the test group and the control group shall be compared using Pearson χ^2 test or Fisher exact test;
- (2) Missing data at baseline shall not be estimated;
- (3) Missing values of main evaluation indicator (success rate of device cutting and anastomosis) shall not be estimated;
- (4) Missing values of secondary evaluation indicators shall not be estimated;
- (5) The wrong and unreasonable data shall be regarded as missing values;
- (6) Missing values for all safety measures shall be estimated using the worst case;
- (7) Sensitivity analysis shall be performed on the safety data set for all the safety evaluation indicators after treatment.

9.8 Procedures for Reporting Deviations from Original Statistical Plan

In general, only the statistical analysis contents pre-specified in the statistical analysis plan can be presented in the clinical trial report of this trial. The increased need for statistical analysis due to various indeterminate reasons will be considered exploratory only.

9.9 Selection Criteria and Reasons for Subjects Included in Analysis

- (1) Full analysis set (FAS): According to the basic principles of intention to treat analysis, all the subjects who have received at least one treatment with the investigational medical device and had at least one post-baseline observation data are included in the FAS. The FAS will be used as the primary population for baseline data (and efficacy evaluation) in this study.
- (2) Per-protocol set (PPS): It refers to the subjects in the FAS who do not seriously deviate from the trial

protocol, have good compliance and have no missing primary efficacy indicators. The PPS will be used as the primary population for efficacy analysis in this study.

- (3) Safety set (SS): All subjects who have received at least one treatment with the investigational medical device and had at least one safety evaluation. The SS will be used as the primary population for safety analysis in this study.

9.10 Exclusion of Special Information During Hypothesis Validation and Reasons

Not applicable.

10. Data Management

10.1 Database Design

The CRF designer is responsible to annotate the final version of the blank CRF and describe in detail the necessary elements (e.g., variable name, length, type, range of the database variable, etc.).

The database designer determines the requirements for database establishment and database structure according to CRF and annotated CRF, and designs the database.

The database must be thoroughly tested to ensure that the contents and structure of the database are consistent with those of the CRF and annotated CRF, that the specified data are appropriately entered into the database, and that the entered data points are correctly exported into the corresponding database variables.

The database can be officially released for use only after it is approved.

10.2 Data Receipt and Data Entry

Data may be received by mail, courier with tracking for confidentiality, hand delivery by the monitor, web entry, or other electronic means. When the paper CRF data is received, it shall be recorded in the data receiving record form. After the completion of data receiving, the CRF shall be scanned and backed up; the data received via network or other electronic means shall be stored in the corresponding documents and records to confirm the data source.

Data entry can be double-entered by blind verification (two persons independently enter the data, and the inconsistency is resolved at the second entry), and double-entered by interactive verification (the second entry resolves the inconsistency between the two entries, and pays attention to the value of the first entry), manual review of single entry and other forms.

Prior to data entry, the data manager shall train the entry clerk or investigators on system operation or matters specific to this study.

After the completion of data entry, the database data shall be checked with original data to ensure that the database data are consistent with original data.

10.3 Management of Data Verification and Query

Data verification includes manual inspection of data, and computer verification designed for identification error, invalid data range, data integrity, protocol violation and consistency inspection. Data administrator shall complete relevant procedures of preparing and approving the Data Verification Plan before the first data receipt to determine the data verification content. Programmer will program and test the verification procedure according to the finalized Data Verification Plan.

Based on the results of manual inspection and computer verification, the data administrator will resolve the doubtful data in the form of data query form, data processing agreement and self-evident correction, and update the database according to the answers or descriptions in relevant documents.

10.4 Database Lock and Data Archiving

The clinical trial database must be closed and locked to ensure the consistency and integrity of the data results, data analysis and data submitted to the official authorities.

When all data have been received and all data cleaning has been completed, all queries have been answered, all unanswered questions have been explained, the data review report has been confirmed, the division of statistical analysis dataset has been clarified, and the data quality control has been completed (the database quality is up to standard), the data administrator shall apply for database locking to the project team members. After being jointly confirmed by the investigators, the sponsor and project team members, the data editing authority of the database will be taken back to lock the database. After the database is locked, the data administrator shall transfer the data to the statistical analysis department in a timely manner. If data correction is required, the review and approval of Database Unlock Application Form shall be completed, and the database shall be unlock, so that corresponding personnel can enter the database to correct the data.

The data administrator shall establish the Data Management Project Folder, and shall be responsible for archiving the data or record files when they meet the conditions for archiving. Upon completion of the study, the data administrator will revoke the database access of relevant personnel, create a copy of the database, and archive the study documents and data management documents.

11. Feasibility Analysis

11.1 Probability Analysis of Success

- (1) According to the relevant clinical trials, it shows that such product is safe and reliable; and such product has been used in clinical practice for many years, and its safety and reliability have been verified.
- (2) The company producing the study medical device has rich experience. The product quality is reliable and meets the requirements of clinical trials; the device has passed the product registration inspection designated by NMPA (original CFDA).
- (3) The clinical institutions undertaking this trial have complete instruments, equipment and technical

resources.

- (4) The investigators are excellent academic leaders in related disciplines in China, all of whom are physicians at or above the deputy director level, have rich clinical trial experience, have received GCP and relevant training, and can carry out treatment in strict accordance with the clinical trial protocol.

11.2 Probability Analysis of Failure

- (1) Too many subjects withdrawing from the trial due to SAE;
- (2) Too many subjects lost to follow-up;
- (3) The IFU and label information of test product can not effectively instruct the operator;
- (4) Unskilled operation of investigators.

11.3 Conclusion

In conclusion, this product can be used for stapling lung tissues for the benefit of lung patients. This study has great feasibility and high possibility of success.

12. Quality Control of Clinical Trial

12.1 Quality Assurance

12.1.1 Definition

Quality assurance is defined as planned and systematic activities. These activities are established to ensure that the study, generation, citation (recording) and reporting of data are conducted in compliance with *Good Clinical Practice* (GCP) and applicable regulatory requirements.

12.1.2 Audit

Audit refers to a systematic independent inspection of clinical trial-related activities and documents organized by the sponsor to determine whether the implementation of such activities and the recording, analysis and reporting of data comply with the clinical trial protocol, standard operating procedures, GCP and relevant applicable management requirements.

12.1.3 Inspection

Inspection refers to the activities of the food and drug regulatory authority in reviewing and supervising relevant documents, facilities, records and other aspects of a clinical trial.

12.1.4 Clinical trial monitoring

The sponsor shall be responsible for comprehensively tracking and monitoring the implementation of clinical

trial, ensuring that the trial is conducted in compliance with GCP and relevant regulations, and is in compliance with the trial protocol. NMPA (original CFDA) may also inspect the clinical trial during or after the trial is completed, and promptly notify the sponsor if investigators receive such notification.

The monitor will make periodic visits to the conduct and completion of the trial. Monitors will check the completeness of case records, the accuracy of case report form content, verify trial data, check the investigator's compliance with the trial protocol and the GCP, and ensure the accuracy of distribution, storage and counting of the study device.

12.2 Quality Control

12.2.1 Definition

Quality control is defined as operational techniques and activities performed under a quality assurance system, such as monitoring, which is intended for verifying that the study-related activities meets the requirements.

Quality control should be applied to each stage of data processing to ensure that all data are credible and correctly processed.

12.2.2 Study monitoring

The authorized and qualified monitor will visit the clinical trial centers regularly according to the monitoring plan to verify the source data, the investigator's compliance with the protocol and regulations, and to assist the investigators.

12.2.3 Laboratory quality control

Each clinical trial center laboratory should establish unified standards for experimental test indicators, standard operating procedures and quality control procedures. Each experimental test item must adopt the national legal unit of measurement, and the test report items must be complete (including date, test item, test result and normal range), and relevant personnel should sign the report. Special inspection items must be specially assigned to be tested.

12.2.4 Quality control for other inspections

Each clinical trial center should establish unified evaluation criteria, standard operating procedures and quality control procedures for the evaluation of vital signs and damage. Each clinical trial center should evaluate the subject's condition according to the unified standard and standard operating procedures.

12.2.5 Qualification of study personnel

The investigators participating in the clinical trial must have the professional expertise, qualification and ability of clinical trials, and pass the qualification review. The personnel requirements shall be relatively fixed.

12.2.6 Study staff training before study initiation

The sponsor is responsible for the training of study personnel before the start of the trial to help the clinical study personnel get a full understanding and understanding of the overall situation of the trial, protocol, CRF, etc.

13. Ethical Issues and Informed Consent of Clinical Trial

13.1 Ethical Considerations

This clinical trial must be carried out in accordance with the *Declaration of Helsinki*, relevant medical device clinical trial regulations and laws in China. The study protocol must be approved by the Ethics Committee (EC) of the clinical trial institute before the study is initiated.

Before each subject is enrolled in the study, the study physician is responsible for explaining the contents of the trial to the subject in a verbal or written language understandable to the subject and carefully answering the subject's questions. Each subject must be given a written subject informed consent prior to enrollment. It is the responsibility of the study physician to obtain informed consent prior to each subject entering the study and to retain it in the study file.

This study is conducted based on the *Declaration of Helsinki*, *Good Manufacturing Practice* (GMP) for medical devices, *Good Clinical Practice* for medical devices, applicable laws and regulations in the People's Republic of China.

Based on ethical considerations, in accordance with GCP, relevant laws and regulations applicable within the territory of the People's Republic of China, the design and implementation of this study put the safety of subjects first on the requirements that the study medical device must be tested firstly when used under clinical conditions. Therefore, in this pre-clinical study, the biological safety of the materials used is tested firstly, and the purpose is to minimize the potential risks in the use of the materials. In order to fundamentally elucidate safety issues, various criteria for this study are also identified.

13.2 Approval of Study Protocol

Before the clinical trial, the investigators should submit the study protocol, informed consent form and other relevant documents to the medical ethics committee of the hospital where the study site is responsible for. The clinical trial can only be initiated after obtaining the approval of the Ethics Committee. Any amendment

to the study protocol must be approved by the Ethics Committee before implementation. Serious adverse events during the clinical trial should be reported to the Ethics Committee in written form in a timely manner.

13.3 Informed Consent Process and Informed Consent Form Text

Before each subject is enrolled in this study, the investigators are responsible for introducing the contents and description of the informed consent form to him/her or his/her guardian completely and comprehensively in written form.

The informed consent form should be written in a language understandable to the subject or his/her guardian. The informed consent form shall contain no content that may cause the subject to give up his/her legal rights and interests, and that exempt the clinical trial institution, the investigators, the sponsor or his/her agent from taking the responsibilities.

Each subject must be given a written subject informed consent prior to enrollment. It is the responsibility of the investigators to obtain informed consent before each subject enters the study and to retain it in the study file.

For an incapacitated subject, if the Ethics Committee agrees in principle and the investigators believe that the subject's participation in the clinical trial is in his/her own interests, he/she can also enter the clinical trial, but his/her guardian should sign and date it before the trial.

When the subject or his/her guardian is unable to read, a witness shall be present during the process of informed consent. After the informed consent form is explained in detail, the witness reads the informed consent form, which shall be consistent with the oral informed consent. After the subject or his/her guardian gives oral consent, the witness signs and dates the informed consent form, and the signature of the witness shall be on the same day as that of the investigator.

See the Informed Consent Form for details.

13.4 Personal Data and Data Protection

- (1) All data obtained in clinical trials are subject to data protection. The investigators should not disclose the subject's name and other personal data (excluding date of birth/age and gender).
- (2) The CRFs and other documents transmitted to the sponsor will not contain names, but only the study code of the subject.
- (3) Similarly, the data used for statistical evaluation can only be performed under the study code of the subject. Only the investigator can identify the subject's name and other personal information through the study code.
- (4) During the study, if the subject name needs to be identified for medical reasons, all relevant personnel are obligated to keep it confidential.

14. Provisions for Reporting Adverse Events and Device Defects

14.1 Adverse Event

14.1.1 Definition of adverse events

Adverse event (AE): refers to any adverse medical event occurring in the clinical trial process, no matter whether it is related to the investigational medical device or not.

14.1.2 Adverse event grading

- (1) Mild: Symptoms inquired.
- (2) Medium: Symptoms that are actively described but can be tolerated.
- (3) Severe: Intolerable symptoms with objective manifestations.

14.1.3 Reporting and handling of adverse events

All the adverse events and serious adverse events observed in the clinical trial should be clearly filled in the adverse event column of case report form. The duration of adverse event (i.e., the date of occurrence and disappearance of adverse event), severity, relationship with the clinical trial, and corresponding treatment and other handling measures should be recorded. Subjects occurring such event shall be followed up to end of the trial or the adverse event outcome.

14.1.4 Correlation between adverse event and product

The correlation between all adverse events and study interventions must be judged by the investigators. The investigators then assessed based on temporal relationship and clinical experience. The followings are the evaluation criteria and judgment results of correlation:

- ① The adverse event occurs within a reasonable time after the application of the device;
- ② The adverse event can not be explained by other reasons;
- ③ The adverse event is relieved or disappeared after the device is not applied;
- ④ The adverse event reappears after reapplication of the device.

Table 3. Correlation evaluation table

Standard \ Results	①	②	③	④
Definitely related	+	+	+	+
Probably related	+	+	+	?
Possibly related	+	-	?	?

Unlikely related	-	-	?	?
Unrelated	-	-	-	-
"+" indicates conformity, "-" indicates non-conformity, "?" Indicates unable to judge				

14.1.5 Handling protocol of adverse events that may occur in this clinical trial

- (1) Pulmonary vascular injury: Accidental injury of pulmonary vessels during operation may cause massive hemorrhage, which is life-threatening. The reasons include ① anatomical variations; ② tight adhesions; ③ improper operation; ④ poor exposure. Once it occurs, the vessel rupture should be immediately compressed, but the force should be appropriate to avoid further injury. At this point, if the surgical field is not exposed enough, it should be enlarged. The proximal and distal ends of the ruptured vessels are then carefully dissected. After the distal vascular occlusion, the yarn ball can be removed, the hematocele can be sucked up. The rupture can be seen, and the continuous or bedding suture can be performed with atraumatic suture. Sometimes, the vessels are dissected and blocked from the pericardium. After pulmonary vascular injury, do not panic. Do not use common vascular clamp. You should ask experienced doctors on stage to help deal with. If the rupture is not large, 1 ~ 3 atraumatic Allis forces can be used to clamp quickly after sucking up the blood, and the rupture can also be directly sutured. It is not necessary to free the proximal and distal ends of the breach.
- (2) Pneumothorax: It occurs mostly in patients with bullae, and rupture of bullae causes pneumothorax. In patients undergoing extensive mediastinal lymphadenectomy, mediastinal pleural rupture can also lead to pneumothorax. After the occurrence of pneumothorax, pulmonary ventilation is more and more difficult. Finally respiratory insufficiency occurs, which is life-threatening. Once this complication is found during surgery, the gas in the pleural cavity should be immediately evacuated, or a chest drainage tube should be placed by expanding the mediastinal rupture. The incidence of this complication is low, with a rate of 0.8% reported in the literature.
- (3) Arrhythmia and myocardial ischemia: It mostly occur in patients with a history of heart disease, but no preoperative cardiac manifestations, intraoperative temporary physiological disorders can cause arrhythmia and myocardial ischemia. In order to prevent and reduce such complications, ① careful preoperative evaluation of cardiac function should be performed to prepare the population for medical or other treatment; ② some factors of cardiac dysfunction should be avoided during surgery, such as hypokalemia, hypovolemia, hypovolemia, tachycardia and acidosis, which should be corrected immediately once the above factors occur; ③ excessive extrusion and stimulation of the heart should not be performed during surgery. Once arrhythmia and myocardial ischemia are caused by extrusion or stimulation, the operation should be immediately stopped, and the operation should be continued after

- the cardiac function recovers; ④ Conduct ECG monitoring during the surgery. Once arrhythmia and myocardial ischemia are found, the corresponding drug treatment should be rapidly given, and the cardiologist should be invited for assistance when necessary.
- (4) Postoperative intrathoracic hemorrhage: Patients with intrathoracic hemorrhage after pneumonectomy who are forced to re-enter the chest for hemostasis account for about 1% of pneumonectomy. The reasons include: ① Bleeding or oozing at the site of pleural adhesion transection, generally mostly at the top of the thoracic cavity; ② Bleeding after thoracic wall vascular injury, such as bleeding from the intercostal artery or internal thoracic artery. It is not easy to stop spontaneously because the bleeding comes from the systemic circulation with high pressure; ③ Pulmonary macrovascular injury, mostly resulted from loose ligatures, fierce blood loss, which is often too late to rescue. Treatment: In any of the following cases, early thoracotomy and hemostasis should be pursued without hesitation, and adequate amounts should be prepared to supplement hypovolemia.
- (5) Cardiac hernia: Cardiac hernia may occur after pericardiotomy or partial resection without suturing or repairing. This complication is rare, but dangerous, with a mortality rate as high as 50%. It occurs mostly after pneumonectomy, but the occurrence of cardiac hernias after lobectomy has also been reported. Typical clinical manifestations are sudden hypotension, tachycardia, and cyanosis. Causes are intrathoracic vacuum aspiration, endotracheal suctioning, severe cough, postural changes, and positive pressure ventilation. Diagnosis is extremely difficult, mainly based on vigilance and experience. Emergency chest X-ray is of great help in judging right hernia, and the heart can be seen to move from the original to the right side, but it is difficult to judge left hernia. Right hernia not only causes torsion of the superior chamber and vein, but also causes distortion and obstruction of the outflow tract. Left hernias are true strangulated hernias that can seriously affect filling and myocardial supply, and electrocardiograms may present with manifestations similar to myocardial infarction. Once the possibility of cardiac hernia is considered in clinical practice, the patient should be immediately placed on the healthy side, and individual patients have the possibility of cardiac repositioning. If the condition does not improve, bedside thoracotomy should be performed decisively for cardiac reduction and pericardial defect repair. There are various methods for the repair of right pericardial defects, including epicardial and pericardial fixation, artificial material or autologous tissue repair. Left pericardial defects may not be repaired if they extend down to the diaphragm.
- (6) Cardiac tamponade (cardiac tamponade): The pericardium is opened during pneumonectomy, and the treatment of bleeding point may be missed, resulting in blood accumulation in the pericardium. Hypotension, elevated central venous pressure, odd pulse, and heart failure may occur when hemopericardium reaches a certain degree. Ultrasound and radiography can confirm the diagnosis.

Treatment should rapidly drain the hemopericardium, open it from the original incision in the chest, or make an incision under the xiphoid process.

- (7) Arrhythmia: Arrhythmia often occurs after pneumonectomy in patients over 60 years of old. The incidence ranges from 20% to 30% after pneumonectomy and from 15% to 20% after lobectomy. Among all arrhythmias, atrial fibrillation is the most common, followed by sinus tachycardia, and atrial flutter, extrasystoles, nodal rhythm, chronic arrhythmias, and bigeminy may also occur. Paroxysmal atrial tachycardia with block, multifocal atrial tachycardia, tachycardia, sick sinus syndrome and atypical tachycardia are relatively rare. More than half of the arrhythmias events occur in the first 24 hours after surgery, with a peak period on postoperative Day 2 to 3. For the treatment of arrhythmia, the cause, such as improving hypoxia, appropriate sedation and analgesia, correcting water and electrolyte imbalance and maintaining acid-base balance, should be firstly removed. Then, depending on the type of arrhythmia, different drugs or other measures should be applied. Digitalis preparations can be used for tachyarrhythmias, and the general dosage for adults needs to reach 0.8 to 1.2 mg to be effective. Verapamil (verapamil) is effective in terminating tachyarrhythmias, with a slow intravenous bolus of 5 to 10 mg for the first time and repeated after 10 to 15 minutes if necessary. After the control of supraventricular arrhythmias, it is changed to verapamil orally (40~80 mg/time, 3 times/day) for maintenance. For atrial arrhythmias that fail to respond to medical treatment and have hemodynamic disturbances, synchronous direct current cardioversion should be used. The first choice for drug treatment of tachycardia is lidocaine (50~100 mg, intravenous bolus), followed by continuous intravenous infusion maintenance at 1~2 mg/min. Bradycardia may be treated with atropine or intravenous isoproterenol. When third-degree atrioventricular block or sick sinus syndrome occurs, the installation of artificial cardiac pacemaker should be considered.
- (8) Myocardial ischemia and myocardial infarction: The incidence of silent myocardial ischemia after pneumonectomy is about 3.8%. Patients with coronary heart disease and myocardial infarction are prone to experience this situation, which often occurs on the 2nd to 4th day after surgery. Therefore, close cardiac monitoring should be performed after pneumonectomy. Once diagnosed, enteric-coated aspirin (160 to 325 mg/day) may be given. Some beta blockers are used appropriately to prevent myocardial infarction and death. The incidence of myocardial infarction after pneumonectomy is about 1.2%, and patients with preoperative diagnosis of coronary heart disease are prone to experience this situation, with a mortality rate as high as 50% to 75%. Once the diagnosis is confirmed, the cardiology consultation should be urgently requested to assist in treatment.
- (9) Orthostatic hypoxemia: After pneumonectomy or lobectomy (usually right upper lobe or right upper lobe and middle lobe), the patient has no or slight dyspnea in supine position, and the saturation is normal or

slightly lower than normal. However, when the patient sits up or stands up, the dyspnea or dyspnea is aggravated, and the saturation becomes abnormal or further decreased. This is called "orthostatic hypoxemia".

- (10) Pulmonary edema after pneumonectomy: Pulmonary edema after pneumonectomy, especially after right pneumonectomy, should be diagnosed if progressive dyspnea, cyanosis, tachycardia and dysphoria, pink foamy sputum and moist rales in the lungs occur. Although the incidence is not high (about 2% to 5%), the mortality rate is high (7% to 80%). The mechanism of this complication is not very clear, but clinical observations and experiments have demonstrated that excessive perioperative fluid infusion is an important cause. As fluid filtered out of the pulmonary capillaries exceeds the capacity for lymphatic reflux, fluid begins to accumulate in the peribronchiolar space. The lungs become stiff, and the work of breathing increases. When the peribronchial space is completely filled with water, the alveoli are also quickly affected, so hypoxemia or even death occurs. The treatment includes oxygen inhalation, fluid restriction, application of morphine and diuretics. In severe cases, endotracheal intubation and mechanical respiration should be performed again with PEEP ventilation mode.
- (11) Respiratory insufficiency: Incision pain, chest wall activity limitation, cough and weakness, airway secretion retention after pneumonectomy may cause atelectasis and inflammation, coupled with flatulence and other factors will cause respiratory insufficiency, manifested as shortness of breath, difficulty, rapid pulse, restlessness, elderly patients and even confusion. Blood gas tests show a decrease in PO_2 and PCO_2 in the early stage, followed by an increase in PCO_2 . Chest X-ray shows partial or total atelectasis of the remaining lungs. Once this disease is suspected, secretions in the airway should be immediately removed in addition to oxygen inhalation. If fiberoptic bronchoscopy suction does not work, tracheotomy should be considered and mechanical breathing should be performed if necessary. This complication often occurs in people with poor lung function. Therefore, pulmonary function assessment is very important. Antibiotics, bronchodilators, smoking cessation, chest physiotherapy, postural drainage and appropriate exercise can improve pulmonary function and reduce the occurrence of postoperative respiratory insufficiency.
- (12) Massive atelectasis: The incidence of massive and severe atelectasis after lobectomy has been reported to be about 7.8%, and the incidence after right upper or right upper middle lung resection is higher than that after right lower or left lobectomy (15.5%, 6%). Patients present with shortness of breath, varying degrees of cyanosis, shortness of breath, accelerated heart rate, and even high fever. Careful physical examination and chest X-ray can help in the diagnosis. Once the diagnosis is confirmed, tracheal and bronchial aspiration and lavage should be performed immediately.
- (13) Lobar torsion and gangrene: Lobar rotation of 180° or more along the bronchovascular pedicle, called

lobar torsion, is mostly caused by excessive intraoperative traction and inversion of the lobe. In addition, if the oblique and horizontal fissures are fully developed, the right middle lobe is prone to torsion after resection of the right upper lobe or right lower lobe. After lung lobe torsion, it causes vascular occlusion, which in turn causes pulmonary infarction and pulmonary gangrene. Once this disease is suspected, the chest should be opened immediately. If the twisted lobe remains viable, it can be reduced and fixed in its normal position. Otherwise, lobectomy should be performed.

- (14) Pulmonary infarction after pulmonary angioplasty: Patients who require pulmonary angioplasty at the same time for lobectomy or bronchial sleeve lobectomy may have thrombosis due to angulated bending or stenosis of pulmonary vessels after surgery, thus causing pulmonary infarction. Patients have low-grade fever, increased density of the lung parenchyma on chest X-ray, and no bronchial obstruction or stenosis can be seen on fiberoptic bronchoscopy. The bronchial mucosa is initially normal and later becomes congested, edematous, and cyanotic. Pulmonary angiography or pulmonary perfusion scan can help in the diagnosis. Once the diagnosis is established, the remaining pneumonectomy should be performed.
- (15) Postoperative pneumonia: The occurrence of pneumonia is closely related to airway secretion retention and atelectasis. The incidence of postoperative pneumonia is 7% in lobectomy and 6.6% in pneumonectomy, respectively. Aspiration pneumonia most commonly occurs in patients with bronchopleural fistula. The condition is prone to develop acute respiratory distress syndrome, and the mortality rate can be 40%. Patients accepting postoperative mechanical respiratory assistance also often suffer from pneumonia. Most pneumonias are caused by bacteria, while viruses, fungi, and mycoplasmas may also be pathogens. In treatment, it is mainly based on the antibiotics application, nutritional support and respiratory tract management.
- (16) Prolonged air leak: Air leaks from the rough surface of the lung after lobectomy or segmentectomy are common. With the expansion of the remaining lung and the elimination of the residual cavity, air leakage generally stops 2 to 3 days after surgery. If there is still air leakage more than 7 days after surgery, it is called "prolonged air leakage". The incidence after pneumonectomy is 15.2%. Although prolonged air leaks rarely cause other complications, the length of hospital stay will be prolonged by 5 to 6 days. The management of prolonged air leaks varies according to the condition, while the experience of thoracic surgeons also plays a significant role, with no single and fixed treatment mode. Increase or decrease the attraction force of closed thoracic drainage system, and install another thoracic drainage tube. Use medical glue to block the small bronchi of air leaking lung segment. Spraying fibrin glue on the rough surface of lung or spraying talc on thoracic cavity can be applied. The second thoracotomy is rarely require. However, if there is an infection in the thoracic cavity and there is a lot of air leakage,

bronchopleural fistula may occur. After the diagnosis of bronchopleural fistula is confirmed, reoperation is required.

- (17)Bronchopleural fistula: The incidence of bronchopleural fistula after pneumonectomy or lobectomy ranges from 1% to 4%. Residual tumor at the bronchial stump, preoperative radiotherapy, diabetes, positive sputum for pulmonary tuberculosis, and right pneumonectomy are high risk factors for the development of fistulas. Closed thoracic drainage should be performed promptly after the diagnosis of bronchopleural fistula, and systemic antibiotics should be given. Before closed thoracic drainage, the patient should be kept in the affected lateral decubitus position to prevent intrathoracic effusion from filling into the lung. Patients with small fistulas (about 3 mm) have the possibility of self-healing after closed drainage. Fibrin glue injection via fiberoptic bronchoscopy can promote the healing of small fistulas. In patients with large fistulas, reoperation may be considered within 7 days after surgery. The bronchial stump may be resutured, embedded with muscle flap, pericardium, or greater omentum, and lavage with postoperative thoracic antibiotic solution to make it sterile. Bronchopleural fistulas more than 7 days are mostly associated with severe thoracic infection and are no longer suitable for re-suturing of the broken end, and thoracic fenestration should be performed 2 weeks after closed thoracic drainage when the mediastinum has been fixed. Fistulas can heal in 30% of patients after fenestration. If the fistula is not healed, thoracoplasty or exposure of the main bronchus through a midline incision should be taken into consideration, which is cut and closed at the root of the main bronchus. Bronchopleural fistula after lobectomy is sometimes solved by total resection of the remaining lung.
- (18)Persistent residual cavities: After lobectomy, especially in the elderly, persistent residual cavities often appear in the pleural space, most of which are located at the top of the chest and a few above the diaphragm. Most postoperative residual cavities are asymptomatic, and disappear after a few months with expansion of the remaining lung, rise of the diaphragm, displacement of the mediastinum, collapse of the chest wall, and organization of residual pleural fluid. Therefore, no special treatment is required. A small number of patients with residual cavity infection, and even complicated by bronchopleural fistula, will present with symptoms. They should be treated as empyema and bronchopleural fistula. It is estimated that residual cavities will occur after surgery, especially after upper lobectomy. The parietal pleura can be dissected to form a tent to convert the intrapleural residual cavity into extrapleural.
- (19)Pleural effusion: After lobectomy or smaller lung resection, it is common to have a small amount of effusion at the bottom of the chest. Over time, the effusion gradually absorbs and organizes and may not be treated. However, when a large amount of effusion occurs, the cause should be identified. Thoracentesis should be performed or closed drainage should be installed again. Otherwise, there will be effusion infection and consequences affecting lung function.

- (20) Empyema: The incidence is 1% to 3%. Empyema may occur simultaneously with bronchopleural fistula or alone. Intraoperative thoracic contamination, secondary thoracotomy for hemostasis or repair of bronchial fistula, and postoperative mechanical respiration are risk factors for concurrent empyema. Clinically, the patient is febrile, with toxic symptoms, increased white blood cell count, and elevated C-reactive protein levels. The diagnosis can be established by evacuation of pus by thoracentesis. Once the diagnosis of empyema is established, a series of measures for the treatment of empyema, such as thoracic drainage and the use of antibiotics, should be taken immediately.
- (21) Chylothorax: Chylothorax after pneumonectomy is rare, with an incidence of about 0.05% to 0.74%. If the thoracic duct is injured during pneumonectomy, or the lymphatic vessels connect the mediastinal lymph nodes and thoracic duct, or the lymphatic vessels of the lung, this disease may occur after surgery. Chylothorax may complicate left or right, lobar or total, upper or lower lobectomy. Chylothorax after lobectomy is mainly characterized by pleural fluid becoming milk-like after eating and increased pleural fluid drainage. Chylothorax after pneumonectomy is characterized by increased pleural effusion, mediastinal displacement, dyspnea, and even decreased blood pressure. Milk-like fluid can be aspirated by thoracentesis. Lymphangiography can be done in conditional hospitals, either to show thoracic duct fistulas or not. Those who do not show a fistula are often cured conservatively. Once the diagnosis of chylothorax is established, fasting, intravenous nutritional support, and maintenance of patent thoracic drainage should be performed immediately. In patients with lobectomy, the fistula is 50% likely to be closed and cured by itself. Conservative cure of chylothorax after pneumonectomy is unlikely, and if the daily pleural fluid exceeds 300 mL, low thoracic duct ligation should be performed by thoracotomy, or closure with fibrin glue. The daily drainage fluid of chylothorax after lobectomy and then thoracotomy is more than 500ml.
- (22) Esophageal injury: When pneumonectomy is performed for pulmonary inflammatory lesions or lung cancer, if mediastinal anatomy and mobilization are difficult, it may damage the esophagus, with an incidence of 0.5% reported in the literature, and more than 90% events occur on the right side. Preoperative gastric tube placement may reduce the chance of esophageal injury. When esophageal injury is suspected intraoperatively, air or methylene blue can be injected into the esophagus. Once the diagnosis is confirmed, it should be repaired immediately. If missed diagnosis occurs during surgery, mediastinal infection and empyema will occur after surgery.
- (23) Wound infection and wound dehiscence: Wound infection and dehiscence after pneumonectomy are relatively rare. In case of incision infection, there are corresponding symptoms and signs, and antibiotics should be drained and used in time. Early incision dehiscence can be sutured again; late dehiscence, mostly as a result of infection, should be debrided, drained, and sutured after infection control.

- (24) Subcutaneous emphysema: The residual gas in the pleural cavity after pneumonectomy can enter the soft tissue of the chest wall through an intercostal incision when the patient changes position or coughs, forming subcutaneous emphysema. If there is no bronchopleural fistula, the subcutaneous emphysema is mild and confined to the vicinity of the incision. No special measures are generally taken.
- (25) Phrenic and recurrent laryngeal nerve injury: When lung cancer or metastatic lymph nodes invade the phrenic nerve, the phrenic nerve has to be cut or removed to remove the tumor, which is the most common cause of pulmonary resection complicated by phrenic nerve injury. This was followed by severe adhesions on the separated mediastinal surface as well as phrenic nerve injury when clamping or coagulation of bleeding points near the phrenic nerve. Phrenic nerve injury affects the ventilatory function of patients, and special care should be taken for those with marginal pulmonary function.
- (26) Spinal cord injury: During pneumonectomy, if there is bleeding from the spinal costal horns, blood may flow into the spinal canal, causing spinal cord compression and paraplegia. If bleeding is controlled with hemostatic materials, these materials have the potential to slip into the spinal canal to compress the spinal cord. If bleeding continues at the spinal costal angle, neurosurgeons should be asked to help enlarge the hole and stop bleeding and oozing with bipolar electrocautery and other special methods.
- (27) Dura mater tear - arachnoid pleural fistula: This complication often occurs during lung cancer surgery. When the tumor invades the posterior chest wall and spinal costal angle, it may tear the dura mater during separation, resulting in arachnoid pleural fistula, or cerebrospinal fluid leakage. During pneumonectomy. If clear fluid is found to flow out from the posterior end of the incision and near the spino-costal angle, this disease should be highly suspected, and neurosurgeons should be consulted to enlarge the foramen, or another laminectomy should be performed to directly suture the dural cleft. If this disease is not found during surgery, it is not easy to confirm the diagnosis in the early postoperative period because cerebrospinal fluid mixed in the blood is not easy to identify. The possibility of this complication should be vigilant when there are gradually less blood components in the pleural fluid and the pleural fluid becomes clear, while the total amount is still many. Cranial copy and CT can show pneumocephalus, and radionuclide myelography can further help in the diagnosis. Neurosurgical consultation should be requested as early as possible after the diagnosis is confirmed. Although some cases can heal spontaneously, preparation should be made for posterior laminectomy. Patients with cerebrospinal fluid leakage may be asymptomatic, but also may be with symptoms, such as headache, seizures, mental status changes, meningeal irritation. There are even death cases reported. Bacterial meningitis is rare.
- (28) Peritumoral vascular embolism: Tumor emboli can be detached from the pulmonary veins and flow into branch vessels of the aortic arch and femoral arteries. Some other embolized vessels reported in the

literature are the cerebral, upper extremity, and mesenteric arteries. 84% of embolizations are caused by surgery, 16% are spontaneous, and the mortality rate is as high as 50%. When possible, embolectomy should be performed. The method of intraoperative prophylaxis is mobilization, ligation, and transection of the pulmonary veins at the proximal end of the intrapericardial tumor thrombus, and if necessary, resection of a small portion of the left atrium.

- (29)Pulmonary embolism: The incidence of pulmonary embolism after pneumonectomy ranges from 1% to 5%. Most emboli originate from veins in the lower extremities or pelvis, and a small proportion from the pulmonary artery stump. The pulmonary artery stump is prone to thrombus formation due to its ruffling and blind tubular shape, where blood slowly swirls. Right embolization after pneumonectomy is rare because of the stump of the right pulmonary artery. The mortality rate of pulmonary embolism is as high as 50%, often with sudden onset and death. Lower limb infusion should be avoided. Patients are encouraged for early physical activities. Pulmonary artery treatment with stapler can reduce the occurrence of this complication.
- (30)Deep Vein Thrombosis: Ziomek et al. (1993) reported this complication. Of the 77 patients with lung cancer, 4 had preoperative deep vein thrombosis and 11 had postoperative deep vein thrombosis, for an overall incidence of 19%. Of the 17 patients who received preoperative oral aspirin or ibuprofen, none developed thrombosis after surgery, and it was considered that patients with adenocarcinoma, large tumors, late TNM stage and relatively large surgery were prone to deep venous thrombosis and embolism.
- (31)Renal failure: Renal failure after pneumonectomy is rare. However, 197 cases of pneumonectomy have been reported. Renal failure occurs in 15% of patients after operation. The elderly over 70 years old are an important cause of postoperative death.
- (32)Cerebrovascular accident: For patients with lung cancer, a distinction should be made between a true "stroke" or a "brain metastasis" that is not diagnosed preoperatively. CT and MRI can help differentiate.
- (33)Gastrointestinal bleeding: Gastrointestinal bleeding after pneumonectomy is rare, often accompanied by other fatal complications (e.g., severe respiratory failure, sepsis, etc.), which mostly manifests before the end of life.
- (34)Pyogenic empyema: After pneumonectomy, the chest residual cavity can be infected and form empyema. Bacteria come via blood from other parts of the body, or from insidious bronchopleural fistulas. The clinical presentation of this complication is not obvious at the beginning, and typical symptoms and signs of empyema appear later.
- (35)Hemorrhagic empyema: It has been reported in the literature that several years after pneumonectomy, a gradually growing hematoma appears in the thoracic stump with blood from granulation tissue in the

stump. Hematomas can cause mediastinal displacement and compression of the unaffected lung.

Re-operation to remove the hematoma can cure this disease.

- (36) Fungal empyema: Late fungal empyema after pneumonectomy is mostly caused by *Aspergillus fumigatus*. In treatment, the best way is thoracoplasty. Open drainage should not be advocated because reoperation, such as muscle flap packing or thoracoplasty, is also required after open drainage.
- (37) Device-related adverse event: The test device may cut to nerves and great vessels during use. The inability to deploy or retract the reload may occur during the deployment of the reload. There may also be battery or motor failures that render the device inoperable. In such cases, the necessary treatment should be performed according to the actual situation, such as conversion to thoracotomy or device replacement.

14.1.6 Handling of adverse events that may occur in this clinical trial

The above are foreseeable adverse events of this study. However, when the subjects suffer from adverse events, the investigators may treat adverse events according to clinical diagnosis and treatment regulations as needed and record them faithfully.

14.2 Serious Adverse Event

14.2.1 Definition of serious adverse event

Serious adverse event refers to events occur in the process of clinical trial that causes death or serious deterioration of health, including fatal disease or injury, permanent damage to the body structure or body function, requiring hospitalization or prolonging the hospitalization, requiring medical or surgical intervention to avoid the permanent damage to the body structure or body function.

14.2.2 Handling of serious adverse events

- (1) In case of any serious adverse event in the clinical trial, the investigators should immediately take appropriate therapeutic measures for the subjects, report to the management department of medical device clinical trial of the clinical trial institution in written form, and notify the sponsor in written form.
- (2) The management department of medical device clinical trial shall submit a written report to corresponding ethics committee, food and drug regulatory authority, and health commission authority of the local province, autonomous region or municipality directly under the central government where the clinical trial institution is located within 24 hours.
- (3) For death-typed events, the clinical trial institutions and investigators should provide all required data to the ethics committee and the sponsor.

- (4) For serious adverse events and device defects that may cause serious adverse events, the sponsor shall report to the filed food and drug regulatory authority and the health commission authority at the same level within 5 working days after being informed, notify other clinical trial institutions and investigators participating in the trial, and timely notify the ethics committee of the clinical trial institution through its medical device clinical trial management department.

14.3 Reporting and Handling of Device Defects

Device defects refer to unreasonable risks that may endanger human health and life safety under normal use of a medical device during a clinical trial, such as labeling errors, quality problems and failures, etc.

Possible device defects: identification error, product quality problem; design defect; damage of disinfection package, etc.

All device defects are documented on the Device Defect Record form, and device defects that could have led to serious adverse events shall be reported and documented within 5 working days.

15. Device Management

15.1 Study Device Management

Clinical trial institutions shall establish a strict registration system for the distribution of study devices. The sponsor shall assign a special person to directly deliver the study medical devices and control medical devices to each clinical trial institution. Clinical trial institutions and the sponsor shall establish perfect device receiving formalities, and the records shall include the date, quantity, batch number and expiry date. Each clinical trial institution uses a special *Record Form for Use of Medical Devices in Clinical Trials* to register the name of the subject, date of use, and signature of the device administrator, etc.

15.2 Device Packaging

The external package and number of study medical devices shall be unified, and the number, name, specification and model, production date, shelf life and the words "investigational" shall be indicated on the package.

15.3 Recovery of Remaining Device

Upon termination of the trial, the sponsor will collect all unused study medical devices and control medical devices, and the sponsor is responsible for destroying the unused medical devices. All supplies must be accounted for at the end of the study.

16. Deviations from Clinical Trial Protocol and Provision for Amendments to Clinical Trial Protocol

16.1 Deviations from Clinical Trial Protocol

16.1.1 Definition of protocol deviation

Deviation refers to a situation where the requirements of the clinical trial protocol are not followed intentionally or unintentionally.

16.1.2 Protocol deviation grading

Important protocol deviations, which have one or more of the following attributes: (1) affecting the safety and rights and interests of subjects; (2) affecting the willingness of subjects to continue to participate in the trial; (3) affecting the quality and integrity of data;

Minor protocol deviations: protocol deviations without the 3 attributes listed above are minor protocol deviations.

16.1.3 Significant protocol deviations

- (1) The subject does not meet any of the inclusion criteria or meets any of the exclusion criteria but is included in the trial.
- (2) Subjects, who meet the criteria for discontinuation during the trial, do not withdraw from the trial. Conditions that should lead to premature discontinuation of the trial include: laboratory test becomes significantly abnormal; the patient withdraws consent; the patient is pregnant; the patient experiences an endpoint event or the patient does not meet criteria to enter the next phase of the study.
- (3) The examination on safety indicators, primary efficacy indicator or key secondary efficacy indicator are not conducted as required by the protocol, thus affecting the scientificity of the study. Including: ① The above examinations are performed outside the time window specified in the protocol, which constitute an important protocol deviation when the data of the primary or key secondary efficacy indicators become inappropriate for statistical analysis; ② The specific requirements of the protocol for examination operation are not met; ③ The examination of safety indicators is not performed, which makes the patient faces safety risks.
- (4) The subject receives incorrect investigational device to the extent that it affects the subject safety or statistical analysis. For example, the device model does not match the use.
- (5) Any serious violation of the principles of GCP in the implementation of the protocol. The main categories include: ① inappropriate process of obtaining informed consent, such as obtaining informed consent by unauthorized researchers; failure to obtain informed consent from patients again or serious delay in obtaining informed consent again after the update of informed consent version; irregular completion/signing of informed consent form (e.g., incomplete signature/date, selection of questions not completed by subjects themselves, signing untrue dates, etc.); existing evidence indicating that the

obtained consent does not meet the principles of complete notification, full understanding, and independent selection; ② improper management of investigational medical devices, such as transportation, receipt, distribution, use/counting of investigational medical devices are not recorded; expired investigational medical devices are still used; investigational medical devices are used for unqualified subjects; unqualified personnel use investigational medical devices for subjects; unsafe storage of investigational medical devices or conditions do not meet the protocol regulations; unqualified labeling of investigational medical devices, etc.

- (6) The principal investigator (PI) fails to perform his/her responsibilities. For example, the PI enrolls subjects before the trial is authorized by the sponsor; the PI fails to timely implement the new protocol after the protocol is amended; the PI fails to supervise the investigators, such as the use of unqualified investigators or the inadequate training. The PI fails to follow the provisions for the reporting of serious adverse events, for example, the adverse event reaches a serious grade but is not reported as a serious adverse event; it is not reported to the relevant parties within the specified time; the investigator often does not make a judgment on the causal relationship between the serious adverse event and the investigational medical device.

16.1.4 Minor protocol deviations

- (1) The visit/observation/examination is out of the time window, but it does not affect the effectiveness in the evaluation of primary efficacy indicators and key secondary efficacy indicators.
- (2) The data points or laboratory parameters, which are required to be observed based on the protocol, are missing, while which do not affect the primary efficacy or key secondary efficacy or safety indicator results.
- (3) The observation/evaluation is insufficient, while which does not affect the primary/secondary key efficacy or safety outcomes.

16.1.5 Conditions that do not constitute a protocol deviation

- (1) If the trial is prematurely discontinued (due to patient withdrawing consent, or patient deciding to discontinue participation in the trial for other reasons), undone examinations resulted from it does not constitute a protocol deviation;
- (2) Discontinuing use of the study device due to an adverse event does not constitute a protocol deviation.

16.1.6 Deviation control measures

16.1.6.1 Sponsor

- (1) The protocol shall be clearly described, and possible protocol deviations shall be considered. The protocol shall be designed in a way to minimize the occurrence of protocol deviations, and possible protocol deviations shall be simply and clearly identified and confirmed.
- (2) In the investigator meeting and monitor training of the study project, the implementation details of the protocol, and the reporting and monitoring of protocol deviations should be given intensive training.
- (3) There is a special protocol deviation control group that regularly reviews the protocol deviation records to identify potential trends and system errors; the correction of protocol deviation should be implemented. If it is necessary to immediately amend the protocol or study plan, the study project should be suspended in time, and the study can be continued after the protocol is amended and the retraining is completed.

16.1.6.2 Investigators

- (1) At the project evaluation stage before the start of the study, the study protocol shall be read carefully, and the enforceability of the protocol in the institution shall be fully discussed with the sponsor.
- (2) The investigators are responsible to explain the importance of compliance with the protocol to the subject at the time of informed consent.
- (3) The investigators shall comply with the protocol approved by the Ethics Committee. In addition, the occurrence of protocol deviation should be recorded and explained immediately when he/she recognizes it.

16.1.7 Retraining

In the following cases, the investigators must receive re-training, and may receive multiple trainings if necessary.

- (1) There are deviations in the understanding of the protocol;
- (2) There are new investigators participating in the trial;
- (3) The study plan changes;
- (4) The investigators are reminded of the protocol details.

16.1.8 Protocol deviation report

In case of any deviation from the clinical trial protocol that affects the rights, safety and health of the subjects or the scientificity of the clinical trial, including the requested deviation and the reported deviation; the investigators should timely report to the management department for medical device clinical trial of the clinical trial institution, and timely notify the sponsor and the Ethics Committee to determine whether the trial can continue to be conducted.

In order to protect the rights, safety and health of subjects, where the deviation occurs in an emergency can not be reported in a timely manner, it shall be reported in written form as soon as possible afterwards in

accordance with relevant regulations.

During the clinical trial, in case of revising the clinical trial protocol, informed consent form and other documents, requesting for deviation and resuming the suspended clinical trial, the implementation can only be continued after obtaining the written approval from the Ethics Committee.

The investigators should strictly follow the clinical trial protocol. Without the approval of the sponsor and the Ethics Committee, or failing to obtain the approval from NMPA (original CFDA) as required, the investigators should not deviate from the protocol or materially change the protocol. However, in the event of an emergency that requires immediate removal, such as when the subject is at an immediate risk, a written report may be made afterwards.

16.2 Amendments to Clinical Trial Protocol

Medical device study protocol should be based on the primary principle of protecting the rights, safety and health of the subjects to the maximum extent. It should be jointly designed and formulated by the medical institution responsible for the clinical trial and the sponsor, and should not be implemented until recognized by the Ethics Committee; if the investigators have any revision to the study protocol, informed consent form and case report form during the study, the clinical trial should not be continued until approved by the Ethics Committee.

17. Direct Accessed Source Data/Source Documents

The authorized monitor has the right to access and verify the source data/documents of the subjects, so as to judge whether the investigators conduct the clinical trial according to the requirements of the protocol, and whether the source data/documents are timely recorded in the subject medical records. The monitor should confirm that the source data/documents are traceable and verifiable. Meanwhile, the data filled in the CRF should correspond to the source data/documents.

When any personnel participating in the clinical trial question the clinical trial data, the investigators shall timely provide the source data/source documents for verification by relevant personnel and carefully reply to the relevant queries.

Source data/documents should be properly kept in accordance with the requirements of relevant laws and regulations. Prior to destruction, it shall be signed and confirmed by the sponsor, and comply with relevant laws and regulations.

18. Finance and Insurance

The sponsor of this clinical trial is responsible for providing the investigational medical device and control medical device required by this clinical trial, and providing the inspection item cost required by this clinical trial. The sponsor shall bear the treatment cost and corresponding economic compensation for the subjects who suffer from injury or death related to the clinical trial, except for the damage caused by the faults of

medical institutions and their medical staff in the diagnosis and treatment.

The specific requirements for finance and insurance are detailed in the relevant agreements.

Medical accidents refer to accidents in which medical institutions and their medical staff violate the laws, trial protocols and guidelines for medical and health management in medical activities, and negligence results in personal damage to patients.

19. Contents Covered in Clinical Trial Report

The clinical trial report shall be prepared in accordance with relevant requirements of GCP and with reference to *Template for Clinical Trial Report of Medical Devices*, including but not limited to:

- (1) General information;
- (2) Abstract;
- (3) Introduction;
- (4) Clinical trial purpose;
- (5) Clinical trial method;
- (6) Clinical trial content;
- (7) General clinical data;
- (8) Investigational medical device and control medical device or control diagnosis and treatment method;
- (9) Statistical analysis method and evaluation method adopted;
- (10) Clinical evaluation criteria;
- (11) Organizational structure of clinical trial;
- (12) Ethical statement;
- (13) Clinical trial results;
- (14) Adverse events identified in clinical trials and handling situations;
- (15) Analysis and discussion of clinical trial results, especially indications, applicable scope, contraindications and precautions;
- (16) Clinical trial conclusion;
- (17) Existing problems and suggestions for improvement;
- (18) List of trial personnel;
- (19) Other situations requiring explanation.

20. Confidentiality

The investigators must strictly keep confidential the study results, protocol and other data, and can not release them by themselves unless authorized by the sponsor in written form. If referenced, a written authorization must also be obtained from the sponsor in advance.

All personal information of subjects is confidential. Besides relevant study personnel, only the sponsor, ethics

committee members and relevant national/local food and drug administration personnel related to this study project can be allowed to access such information. The investigators and sponsor should keep the personal information of the subjects confidential.

21. Agreement on Publication of Study Results

The sponsor and investigators shall agree with the final study report.

The results of the study may be made publicly available as scientific literature and the results may be submitted to authorities. This clause is intended to protect trade secret materials.

All information about the study medical device (e.g., patent application, previously undisclosed manufacturing process provided to the investigators by the sponsor, basic scientific data, etc.) is considered confidential and its ownership belongs to the sponsor. The investigators may not use it for other purposes without the written permission of the sponsor.

Prior to publication or presentation of the results of this study, the investigators allow the sponsor to review the manuscript and comment within 30 days to confirm that the confidential information is not disclosed, and to supplement the relevant information. According to the generally accepted principles of scientific cooperation, the investigators should discuss the manuscript with the relevant personnel of the sponsor and reach a consensus before the manuscript is published.

22. Responsibilities Assumed by Each Party

22.1 Responsibilities of Sponsor

- (1) The sponsor is responsible for initiating, applying for, organizing and monitoring the clinical trial, and is responsible for the authenticity and reliability of the clinical trial.
- (2) The sponsor is responsible for organizing the development and revision of the investigator's brochure, clinical trial protocol, informed consent form, case report form, relevant standard operating procedures and other relevant documents, and organizing the training necessary for the clinical trial.
- (3) Before signing the clinical trial agreement with the clinical trial institution, the sponsor shall provide the clinical trial institution and investigators with the latest investigator's brochure and other relevant documents for them to decide whether they can undertake the clinical trial.
- (4) The sponsor shall not exaggerate the publicity of the mechanism and efficacy of the study medical device in organizing the development of clinical trial protocol.
- (5) During the clinical trial, when the sponsor obtains important information affecting the clinical trial, the investigator's brochure and relevant documents shall be revised in a timely manner, and submitted to the ethics committee for review and approval through the medical device clinical trial management department of the clinical trial institution.
- (6) The sponsor shall reach a written agreement with the clinical trial institution and the investigators on the

detailed rules related to the clinical trial.

- (7) The sponsor is responsible for the safety of the study medical device in the clinical trial.
- (8) Where the sponsor decides to suspend or terminate the clinical trial, it shall notify the management department of medical device clinical trial of all the clinical trial institutions within 5 days, and explain the reasons in written form.
- (9) The sponsor shall ensure that all investigators conducting clinical trials strictly follow the clinical trial protocol, and timely point out and correct the clinical trial institutions and investigators who fail to comply with relevant laws and regulations, this Practice and clinical trial protocol; if the situation is serious or does not change continuously, the trial shall be terminated and reported to the food and drug regulatory authority of the provinces, autonomous regions and municipalities directly under the central government where the clinical trial institutions are located, as well as NMPA (original CFDA).
- (10) The sponsor shall bear the treatment cost and corresponding economic compensation for the subjects who suffer from injury or death related to the clinical trial, except for the damage caused by the faults of medical institution and its medical staff in the diagnosis and treatment.
- (11) The sponsor should bear the responsibility of monitoring and verifying the clinical trial.
- (12) For serious adverse events and device defects that may cause serious adverse events, the sponsor shall report to the food and drug regulatory authority which files the trial, as well as the health commission at the same level within 5 working days after being informed. The sponsor shall notify other clinical trial institutions and investigators participating in the trial, and timely notify the Ethics Committee of the clinical trial institution through its medical device clinical trial management department.
- (13) The sponsor shall ensure that the clinical data in the electronic clinical database or remote electronic clinical data system are controlled and authentic, and form complete validation documents.
- (14) The sponsor shall ensure that the CRFs are rigorously and rationally designed, so that the coordinating investigators can obtain all data from the clinical trial institutions of each sub-site.

22.2 Responsibilities of Clinical Trial Institution and Investigators

- (1) Before accepting a clinical trial, the clinical trial institution should evaluate relevant resources according to the characteristics of the study medical device to determine whether to accept the clinical trial.
- (2) The clinical trial institution should properly keep the clinical trial records and basic documents according to the agreement with the sponsor.
- (3) The investigators responsible for the clinical trial should have the corresponding qualification.
- (4) Before a clinical trial, the management department for clinical trial of medical device of the clinical trial institution should cooperate with the sponsor in submitting an application to the Ethics Committee and submitting the relevant documents as required.

- (5) The investigators should ensure that relevant staff participating in the trial are familiar with the principles, applicable scope, product performance, operation methods, installation requirements and technical indicators of the study medical device; understand the pre-clinical study data and safety data of the study medical device; and master the prevention and emergency treatment methods of possible risks in the clinical trial.
- (6) The investigators should ensure that the study medical device is only used for the subjects of the clinical trial, and should not charge any fee.
- (7) The investigators should strictly follow the clinical trial protocol. Without the approval of the sponsor and the ethics committee, or failing to obtain the approval from NMPA (original CDFA) as required, the investigators should not deviate from the protocol or materially change the protocol. However, in the event of an emergency that requires immediate removal, such as when the subject is at an immediate risk, a written report may be made afterwards.
- (8) The investigator is responsible for recruiting the subject and communicating with the subject or his/her guardian.
- (9) The investigators or other personnel involved in the trial should not coerce or induce the subject to participate in the trial in other improper ways.
- (10) When the investigators find unexpected adverse events of the study medical device during the clinical trial, they should modify the relevant contents of the informed consent form with the sponsor. After reporting to the Ethics Committee for review and approval according to the relevant working procedures, the affected subjects or their guardians should re-sign and confirm the revised informed consent form.
- (11) The investigators shall be responsible for making medical decisions related to the clinical trial. In case of any adverse event related to the clinical trial, the clinical trial institution and the investigators shall ensure to provide adequate and timely treatment, as well as processing for the subjects.
- (12) In case of any serious adverse event in the clinical trial, the investigators should immediately take appropriate therapeutic measures for the subjects, report to the management department of medical device clinical trial of the clinical trial institution in written form, and notify the sponsor in written form.
- (13) The investigators should record all adverse events occur and device defects found during the clinical trial, work with the sponsor to analyze the causes of the events, form a written analysis report, propose the comments on continuing, suspending or terminating the trial, and report to the Ethics Committee by the management department for clinical trial of medical device of the clinical trial institution for review.
- (14) The investigators should ensure that the clinical trial data are accurately, completely, clearly and timely recorded in the CRF.
- (15) The clinical trial institution and investigators shall ensure that the data, documents and records generated

from the clinical trial are authentic, accurate, clear and safe.

- (16) The clinical trial institution and investigators should accept the monitoring and verification by the sponsor, as well as the supervision by the Ethics Committee; and provide all required records related to the trial.
- (17) When the clinical trial institution and investigators need to suspend or terminate the clinical trial when they find that the risks outweigh the possible benefits, or the results are enough to judge the safety and effectiveness of the study medical device, they should inform the subjects and ensure that the subjects receive appropriate treatment and follow-up, report according to the provisions, and provide a detailed written explanation.
- (18) The clinical trial institution and investigators shall report to the food and drug regulatory authority of the province, autonomous region or municipality directly under the central government where the sponsor is located if the sponsor violates relevant regulations or requirements to change the trial data and conclusions.
- (19) At the end of the clinical trial, the investigators should ensure the completion of various records and reports, and timely deliver the required clinical data to the clinical trial institution as required.
- (20) The investigators should provide relevant training for the authorized personnel, and form the corresponding documents.

Note: See the "Clinical Trial Contract" and other relevant documents for the responsibilities assumed by specific parties.

23. References

- [1] Yang Zhiguang, Lin Xingyu, Zhang Peng, et al. 112 cases of lung cancer treated by VAT lobectomy/segmentectomy [J]. Chinese Journal of Gerontology, 2014 (11): 3024-3026.
- [2] Jiang Wei, Xi Junjie, Wang Qun, et al. Discussion on lymph node dissection of VAT radical resection of lung cancer [J]. Chinese Journal of Minimally Invasive Surgery, 2012 (11): 969-972.

24. Supplementary Provisions

Meanings of the academic terms involved in this study protocol:

Medical device clinical trial institution: it refers to the medical institution undertaking medical device clinical trial recognized by NMPA (original CFDA) and NHFPC. Unless otherwise specified, "clinical trial institution" in this protocol refers to "medical device clinical trial institution".

Investigational/study medical device: it refers to the medical device to be applied for registration whose safety and effectiveness are confirmed or verified in the clinical trial.

Sponsor: it refers to the institution or organization initiating, managing and providing financial support for the clinical trial.

Investigator: it refers to the person responsible for conducting the clinical trial in the clinical trial institution.

If a trial is conducted by a team of individuals at a clinical trial institution, the investigator is the responsible leader of that team, also known as the principal investigator.

Ethics Committee: it refers to an independent institution set by the clinical trial institution to review the scientificity and ethicality of the medical device clinical trial project.

Management department of medical device clinical trial: it refers to the division or department set within a clinical trial institution to be responsible for the organization management and quality control of the medical device clinical trial.

Multi-center clinical trial: it refers to the clinical trial conducted in more than three (including three) clinical trial institutions in accordance with the same clinical trial protocol.

Subject: it refers to an individual recruited for the clinical trial of a medical device.

Informed consent: it refers to the process that the subject voluntarily participates in the clinical trial after being informed of all aspects of the clinical trial, as evidenced by the signed and dated informed consent form.

Informed consent form: it refers to the supporting document for the subject voluntary participation in the clinical trial.

Monitoring: it refers to the activities that the sponsor designates special personnel to evaluate and investigate the clinical trial institution and investigator; verify, record and report the data during the clinical trial so as to ensure that the clinical trial carried out can comply with the clinical trial protocol, standard operating procedures, *Good Clinical Practice* for medical devices and relevant applicable regulatory requirements.

Monitor: it refers to the special personnel assigned by the sponsor to monitor the medical device clinical trial project.

Audit: it refers to a systematic and independent inspection of clinical trial-related activities and documents organized by the sponsor to determine whether the implementation of such activities, recording, analysis and reporting of data comply with the clinical trial protocol, standard operating procedures, *Good Clinical Practice* for medical devices and relevant applicable regulatory requirements.

Auditor: it refers to the personnel entrusted by the sponsor to audit the medical device clinical trial project.

Inspection: it refers to supervision and management activities conducted by regulatory authorities on relevant documents, facilities, records and other aspects of clinical trials

Inspector: it refers to the personnel assigned by regulatory authorities to inspect the medical device clinical trial project.

Deviation: it refers to the situation where the requirements of the clinical trial protocol are not followed

intentionally or unintentionally.

Case report form: it refers to the document designed in accordance with the clinical trial protocol to record all the information and data of each subject obtained during the trial.

Endpoint: it refers to an indicator used to assess a clinical trial hypothesis.

Source data: it refers to all information in original records and their approved copies of clinical findings, observations and other activities in a clinical trial that can be used for reconstruction and evaluation of the clinical trial.

Source document: it refers to the printed, visible or electronic documents containing source data.

Adverse event: it refers to any adverse medical event occurring in the clinical trial process, no matter whether it is related to the investigational medical device or not.

Serious adverse event: it refers to any adverse event that results in death or serious deterioration in the health of the subject during the clinical trial, including a life-threatening illness or injury, a permanent impairment of a body structure or a body function, in-patient hospitalization or prolongation of existing hospitalization, medical or surgical intervention to prevent permanent impairment of a body structure or a body function; results in fetal distress, fetal death, or a congenital abnormality or birth defect, etc.

Device defect: it refers to unreasonable risks that may endanger human health and life safety under normal use of a medical device during a clinical trial, such as labeling errors, quality problems and failures, etc.

Standard operating procedure: it refers to the standard and detailed written procedure proposed to effectively implement and complete each work in the clinical trial.

Clinical data: it refers to the safety and performance information obtained from the relevant literature or clinical application of the medical device.

Investigator Statement

I agree that:

1. This clinical trial should be conducted in strict accordance with the *Declaration of Helsinki*, current laws and regulations of China, and the requirements of the study protocol.
2. All required data shall be recorded accurately on the case report form (CRF), and the clinical trial report shall be completed on time.
3. The investigational medical device is only used for this clinical trial. The receipt and use of the investigational medical device will be completely and accurately recorded during the clinical trial, and the records will be kept.
4. The monitor, inspector and regulatory authority authorized or dispatched by the sponsor will be allowed to monitor, verify and inspect the clinical trial.
5. The terms of clinical trial contract/agreement signed by all parties shall be strictly performed.

I have read the clinical trial protocol, including the above statement, and I agree all contents above.

Comments of sponsor

Agree.

Fengh Medical Co., Ltd.



3202811906439

Signature (seal): Zhang Xinghua

Date: 2019.03.25

Comments of investigator

Agree.

Signature: Yu Shusheng

Date: 2019.03.27

Comments of medical device clinical trial institution

Agree.



Signature (seal):

Date: 2019.04.03