

CLINICAL INVESTIGATION PLAN

CLINICAL INVESTIGATION TITLE:	A Post-Market Clinical Evaluation of the ReUnion Reverse Shoulder Arthroplasty (RSA) System
DEVICE NAME:	ReUnion RSA System
CLINICAL INVESTIGATION DESIGN:	Post-Market, Multicenter, Prospective, Non-Randomized
INDICATIONS:	This clinical investigation will adhere to the indications and contraindications for the ReUnion RSA System as are detailed in the device's Instructions for Use.
REGULATORY STATUS:	510(k) Clearance received on 27Dec2013
CLINICAL INVESTIGATION PLAN PHASE:	Post-Approval Clinical Investigation
SPONSOR:	Stryker Orthopaedics 325 Corporate Drive Mahwah, NJ 07430
AUTHOR:	Lindsay Mattfolk
INVESTIGATORS:	Investigators' information is on file at the Sponsor
MEDICAL EXPERT:	Dr. Helmut Mueckter, M.D.
COMPLIANCE STATEMENT:	This clinical investigation will be conducted in compliance with the Clinical Investigation Plan International Conference of Harmonisation Good Clinical Practice (ICH-GCP), and all other applicable regulatory requirements, including the retention of essential documents. Investigators will be trained on the clinical investigation devices and surgical techniques prior to implanting clinical investigation subjects.
CONFIDENTIALITY STATEMENT:	This Clinical Investigation Plan contains confidential information and its use is limited to investigational staff intending to conduct the clinical investigation, Institutional Review Boards (IRBs) and any others charged with reviewing the clinical investigation.
VERSION:	1
DATE:	12Oct2018

Approval Page









APPROVERS			
Role	Name	Signature	Date
<i>Medical Expert</i>	Dr. Helmut Mueckter, MD	  Name des Unterzeichners: Helmut Mueckter Signiergrund: Ich habe dieses Dokument geprüft Signaturzeit: 2018-10-15 15:53:27Z (UTC)	15-Okt-2018 8:53 AM F
<i>Clinical Research Head (CRH)</i>	Georgia Mitchell	  Signer Name: Georgia Mitchell Signing Reason: I approve this document Signing Time: 2018-10-15 14:55:59Z (UTC)	15-Oct-2018 7:56 AM F
<i>Regulatory Affairs (RA)</i>	Jemin Dedania	  Signer Name: Jemin Dedania Signing Reason: I approve this document Signing Time: 2018-10-16 13:56:12Z (UTC)	16-Oct-2018 6:56 AM F
<i>Statistician</i>	Claudia Beimel	  Signer Name: Claudia Beimel Signing Reason: I approve this document Signing Time: 2018-10-16 07:42:53Z (UTC)	16-Oct-2018 12:42 AM

Table of Contents

1. LIST OF ABBREVIATIONS	4
2. SYNOPSIS	4
3. GENERAL INFORMATION AND ADMINISTRATIVE STRUCTURE	5
4. PRODUCT INFORMATION	6
5. RISKS AND BENEFITS	6
6. INTRODUCTION	6
7. CLINICAL INVESTIGATION DESIGN	7
8. OBJECTIVE	7
9. SELECTION OF CLINICAL INVESTIGATION POPULATION	7
10. CLINICAL INVESTIGATION EVALUATIONS, PROCEDURES AND ASSESSMENTS	9
11. STATISTICAL METHODS	10
12. CLINICAL INVESTIGATION PLAN DEVIATIONS	17
13. ADVERSE EVENTS	17
14. REOPERATIONS	19
15. ETHICS	19
16. DATA COLLECTION PROCESS	20
17. CLINICAL INVESTIGATION MONITORING	20
18. DATA HANDLING AND RECORD KEEPING	20
19. REPORTS	20
20. COMPLETION OF THE CLINICAL INVESTIGATION	20
21. ESSENTIAL DOCUMENTS	21
22. PUBLICATION POLICY	21
23. REFERENCES	21
24. CLINICAL INVESTIGATION PLAN SIGNATURE PAGE	23
25. DOCUMENT VERSION HISTORY	24

1. List of Abbreviations

<u>Acronym</u>	<u>Definition</u>
ADE	Adverse Device Event
AE	Adverse Event
ASES	American Shoulder and Elbow Surgeons
CI	Confidence Interval
CIP	Clinical Investigation Plan
CRF	Case Report Form
eCRF	Electronic Case Report Form
EDC	Electronic Data Capture
ICF	Informed Consent Form
IFU	Instructions for Use
IRB	Institutional Review Board
ITT	Intent-to-Treat
LTFU	Lost to Follow-Up
PP	Per Protocol
RSA	Reverse Shoulder Arthroplasty
SADE	Serious Adverse Device Effect
SAE	Serious Adverse Event
SAP	Statistical Analysis Plan
TSA	Total Shoulder Arthroplasty
UADE	Unanticipated Adverse Device Effect

2. Synopsis

Title	A Post-Market Clinical Evaluation of the ReUnion Reverse Shoulder Arthroplasty (RSA) System
Treatment	ReUnion RSA System
Design	<ul style="list-style-type: none"> • Post-Market • Multicenter • Prospective • Non-Randomized
Objectives	<p>The objective of this clinical investigation is to demonstrate the safety and efficacy/performance of the ReUnion RSA System.</p> <p>Efficacy/performance of the procedure will be measured the American Shoulder and Elbow Surgeons (ASES) Shoulder Score.</p> <p>Safety of the ReUnion RSA System will be demonstrated through reporting of device-related intra-operative and post-operative Adverse Events (AEs).</p>
Target Population	80 subjects are to be enrolled in this clinical investigation. Enrolled subjects will be assessed at Pre-Operative, Operative/Discharge, and at 6 Weeks, 6 Months, 12 Months, 24 Months and annually thereafter up to 10 years following the index procedure.

Endpoints	<p><u>Primary Endpoint:</u></p> <p>The objective of the clinical investigation is to demonstrate the efficacy/performance of the device at 24 months, as measured by the ASES Shoulder Score.</p> <p><u>Secondary Endpoints:</u></p> <ul style="list-style-type: none"> • Safety: Incidence of device-related intra-operative and post-operative AEs • Efficacy/Performance: Implant survivorship will be monitored
Inclusion Criteria:	<ol style="list-style-type: none"> Subject is willing to sign the informed consent. Subject is willing and able to comply with post-operative scheduled clinical assessments. Subject is male or non-pregnant female and 18 years or older at the time of surgery. Subject's joint has gross rotator cuff deficiency, a functional deltoid muscle and is anatomically and structurally suited to receive the implant and subject has one or more of the following: <ul style="list-style-type: none"> • Painful, disabling joint disease of the shoulder resulting from degenerative arthritis or rheumatoid arthritis • Failed previous shoulder joint replacement
Exclusion Criteria:	<ol style="list-style-type: none"> Subject has an active or suspected latent infection in or about the shoulder joint. Subject has mental or neuromuscular disorder which would create an unacceptable risk of prosthesis instability, prosthesis fixation failure or complications in postoperative care. Subject has bone stock compromised by disease, infection or prior implantation which cannot provide adequate support and/or fixation to the prosthesis. Subject has anticipated activities which would impose high stresses on the prosthesis and its fixation. Subject is obese such that he/she produces a load on the prosthesis which can lead to failure of fixation of the device or to failure of the device itself. Subject has concomitant disease(s) which may significantly affect the clinical outcome. Subject has traumatic or pathologic fracture of the proximal humerus.

3. General Information and Administrative Structure

3.1. SPONSOR

Stryker Orthopaedics
325 Corporate Drive
Mahwah, NJ 07430

3.2. KEY SPONSOR PERSONNEL

Lindsay Mattfolk
Clinical Study Manager
Role: Clinical Investigation Manager
Email: lindsay.mattfolk@stryker.com
Phone: +1 617 417 2956

Hanna Schlyter
Senior Director Clinical Affairs Trauma & Extremities
Role: Clinical Affairs Head

Dr. Helmut Mueckter, M.D.
Role: Medical Expert

3.3. EDC SYSTEM

Name: IBM Clinical Development
Email: icdhelp@us.ibm.com
Phone: +1 888 500 4247

4. Product Information

All components of the ReUnion RSA System have been cleared and approved for sale and use in the United States prior to starting the clinical investigation. 510(k) Clearance was received on 27Dec2013. This system is to be used only for indications for which it has been approved. Please see the approved Instructions for Use (IFU) and Operative Technique manuals for a detailed description of the medical device(s) and instrumentation as well as the intended use information.

Medical device product traceability will be achieved by capturing the implant lot number.

5. Risks and Benefits

This prospective, multicenter clinical investigation is designed to examine the safety and efficacy/performance of the ReUnion RSA System in accordance with the approved IFU, labelling and instrumentation. The potential risks to subjects are described in the approved IFU and Operative Technique manuals.

Potential benefits resulting for the ReUnion RSA System over other devices and procedures as demonstrated by superior scoring on the outcome survey and positive results on other clinical evaluation measurements would suggest affirmative clinical efficacy.

6. Introduction

The ReUnion RSA System is designed as a reverse total shoulder endoprosthesis to address unrepairable gross rotator cuff tear, rotator cuff arthropathy with pseudoparalysis of the shoulder joint, advanced arthritic and rheumatic disorders affecting the shoulder joint, and failed previous shoulder arthroplasty. The ReUnion RSA System is comprised of a Humeral Cup, Humeral Insert, Glenosphere, Glenoid Baseplate and Screws. It is used with the humeral stems from the ReUnion Total Shoulder Arthroplasty (TSA) System. The intended purposes of the ReUnion RSA System are to achieve pain relief, improvement of range of motion and restoration or improvement of the shoulder function while ensuring long-term replacement of the shoulder joint with sufficient stability of all endoprosthesis components.

7. Clinical Investigation Design

This investigation is a prospective, multicenter clinical investigation. It is anticipated that a total of eighty (80) subjects will be enrolled at approximately 4-7 sites. Neither subjects nor investigators are blinded to treatment and the clinical investigation does not include a contemporaneous control. The clinical investigation has been designed to follow the surgeon's standard of care for joint arthroplasty subjects, which entails clinical evaluation on a regular ongoing basis, or as needed should the subject become symptomatic in the treated joint. The enrollment period is expected to occur over 12 months.

7.1. CLINICAL INVESTIGATION RATIONALE

This clinical investigation will evaluate the safety and efficacy/performance of the ReUnion RSA System.

8. Objective

8.1. PRIMARY ENDPOINT

The primary endpoint of the clinical investigation is to demonstrate non-inferiority of the device to the selected literature controls, as measured by the ASES Shoulder Score at 24 Months post-operative.

8.2. SECONDARY ENDPOINT

In addition to the principal endpoint, information on the following outcomes will be assessed up to 10 years after the index procedure:

Safety: Incidence of device-related intra-operative and post-operative AEs

Efficacy/Performance: Implant survivorship will be monitored

9. Selection of Clinical Investigation Population

Subjects participating in this clinical investigation will be recruited from the investigator's standard subject population, where all patients presenting for primary or revision total shoulder replacement will be evaluated for clinical investigation participation based on the eligibility criteria listed below.

9.1. INCLUSION CRITERIA

- a. Subject is willing to sign the informed consent.
- b. Subject is willing and able to comply with postoperative scheduled clinical and radiographic evaluations.
- c. Subject is male or non-pregnant female and 18 years or older at the time of surgery.
- d. Subject's joint has gross rotator cuff deficiency, a functional deltoid muscle and is anatomically and structurally suited to receive the implant and subject has one or more of the following:
 - Painful, disabling joint disease of the shoulder resulting from degenerative arthritis or rheumatoid arthritis
 - Failed previous shoulder joint replacement

9.2. EXCLUSION CRITERIA

- a. Subject has an active or suspected latent infection in or about the shoulder joint.

- b. Subject has mental or neuromuscular disorder which would create an unacceptable risk of prosthesis instability, prosthesis fixation failure or complications in postoperative care.
- c. Subject has bone stock compromised by disease, infection or prior implantation which cannot provide adequate support and/or fixation to the prosthesis.
- d. Subject has anticipated activities which would impose high stresses on the prosthesis and its fixation.
- e. Subject is obese such that he/she produces a load on the prosthesis which can lead to failure of fixation of the device or to failure of the device itself.
- f. Subject has concomitant disease(s) which may significantly affect the clinical outcome.
- g. Subject has traumatic or pathologic fracture of the proximal humerus.

9.3. WITHDRAWAL CRITERIA

If during the clinical investigation a subject must be withdrawn prematurely, then the procedures outlined in this section must be followed. These procedures should not interfere with the initiation of any new treatments that are necessary to treat a subject's condition. Information on all withdrawn subjects will be documented.

Subjects may be withdrawn from the clinical investigation for any of the following reasons:

- a. **Subject Withdrawal:** A subject may voluntarily withdraw from the clinical investigation at any time and for any reason. The subject should be asked when possible, and without any form of coercion, the reason for his/her decision. If the participant withdraws from the clinical investigation completely, then data collected up until the point of withdrawal will be included in the final analysis.
- b. **Lost to Follow-Up (LTFU):** A subject will be considered LTFU after all reasonable efforts have been made to contact the subject and request his/her continued participation in the clinical investigation. All attempts to contact the subject must be documented and should include at least two attempts to contact the subject by phone and one attempt via a certified letter. Data collected up until the point where the subject is LTFU will be included in the final analysis.
- c. **Removal of Device or AE/Incident:** The discontinuation of a subject's participation in the clinical investigation due to the removal of the ReUnion RSA System or AE/incident that prohibits his/her continued participation must be fully explained. All available information concerning the removal of the device or AE/incident should be provided. Data collected up until the point of removal or AE/incident will be included in the final analysis.
- d. **Death:** The discontinuation of a subject's participation in the clinical investigation due to death must be fully explained. All available information concerning the death or AE should be provided. Removal of a subject from continued follow-up in the clinical investigation due to death will not be considered a device failure unless the death is directly caused by, or attributable to, the ReUnion RSA System. Data collected up until the point of death will be included in the final analysis.
- e. **Other:** A subject may be withdrawn by the investigator if he/she believes it is in the best interest of the subject, or if it is determined by the IRB that a subject's continued participation in the clinical investigation represents an unacceptable risk to the subject. The Sponsor must be notified immediately if this occurs. All data collected up until the point of withdrawal or IRB determination will be included in the final analysis.

A subject may also be withdrawn if the subject is non-compliant with the clinical investigation procedures or visits, or if a selection criteria violation is noted after the subject received the clinical investigation treatment and it is determined that the subject should be discontinued. All data collected up until the point of withdrawal will be included in the final analysis.

10. Clinical Investigation Evaluations, Procedures and Assessments

10.1. METHOD OF ASSIGNING SUBJECTS

No specific methods (e.g. randomization, blinding, or stratification) for assigning subjects are used in this clinical investigation plan (CIP). Consecutive subjects at each site meeting all eligibility criteria will be enrolled in this clinical investigation.

10.2. PROCEDURES

Subjects in the clinical investigation will undergo placement of the ReUnion RSA System. Please see the approved IFU and Operative Technique Manuals for a detailed description of the medical device(s) and instrumentation, intended use information and associated risk. Any additional clinically indicated procedures are permitted as deemed necessary by the investigator.

10.3. FOLLOW-UP EVALUATIONS

Subjects in this clinical investigation will be evaluated at Pre-Operative, Operative/Discharge, and at 6 Weeks (4 weeks – 8 weeks), 6 Months (23 weeks – 29 weeks), 12 Months (48 weeks- 56 weeks), 24 Months (100 weeks – 108 weeks) and annually thereafter. The follow-up evaluations will include assessment of device-related AEs/incidents, radiographs and ASES Shoulder Score.

Investigative site personnel will contact subjects prior to their scheduled follow-up evaluations to encourage compliance with clinical investigation visits and participation.

If a subject misses a visit and is outside of the visit window, then every effort should be made to collect data instead of noting the visit as missed.

10.4. SCHEDULE OF EVENTS

Assessment	Pre-Operative	Operative/ Discharge	6 Weeks ^{a, b} (+/- 2 weeks)	6 Months ^{a, b} (+/- 3 weeks)	12 Months ^{a, b} (+/- 4 weeks)	24 Months ^{a, b} (+/- 4 weeks)	Annually ^b
Informed Consent	X						
Demographics & Medical History	X						
Inclusion/Exclusion	X						
Physical Exam	X		X ^c	X ^c	X ^c	X ^c	X ^c
Surgical Procedure		X					
ASES Shoulder Score	X		X ^d	X ^d	X ^d	X ^d	X ^d
Image Evaluation ^{e, f}	X	X	X	X	X	X	X
Subject Disposition ^g			X	X	X	X	X
Device-Related AEs/Incidents & Reoperations will be collected throughout the course of the clinical investigation.							

	<ul style="list-style-type: none"> a. Follow-up visit schedule to reflect Institutions' Standard of Care practices b. If the subject misses a visit and is outside of the visit window, then every effort should be made to collect data instead of noting visit as missed. c. Evaluation may be collected when subject presents in-clinic for study visit. d. Evaluation can be collected via phone. e. Radiograph collection should follow Institution's Standard of Care practices and no additional x-rays should be made for study purposes. f. CT scans may be collected if part of Institution's Standard of Care practices. g. Subject Disposition assessment will occur at any time point for subject withdrawal prior to the completion of the clinical investigation.
--	---

Table 1: Schedule of Events

11. Statistical Methods

The 24 months post-operative results for subjects implanted with the ReUnion RSA System will be compared to a historical group and results reported by respective clinical outcome data in the scientific literature. The benchmark sources and values will serve as the control group for the ReUnion RSA System subjects.

Hypotheses are developed to allow for a comparison of 24 months post-operative ASES score results and 24 months effectiveness/performance between these two populations.

The Statistical Analysis Plan (SAP) lists all variables/questions within this clinical investigation. Therefore, no additional "evaluation" chapter is required nor needed for this CIP.

Data will be captured via IBM Clinical Development electronic data capture (EDC) system and statistical analysis will be performed using IBM SPSS. All statistical hypotheses tests will be with confidence levels (1- α) of 95% and power (1- β) of 80%. The significance level (α) is 0.05 and the beta-value (β) is set to 0.20. Therefore, p-values ≤ 0.05 will indicate statistical significance.

Results will be presented using summary tables and optionally supported by graphs. For detailed information per variable, see SAP.

The primary endpoint of the clinical investigation is to demonstrate non-inferiority of the device to the selected literature controls, as measured by the ASES Shoulder Score at 24 Months post-operative. Specifically, the proposed study will test the null hypothesis that the population mean is 71.27 points. The criterion for significance (alpha) is set at 0.05. With the sample size of 10 cases, the clinical investigation will have power of 80.3% to yield a statistically significant result. This computation assumes that the population from which the sample will be drawn has a mean of 71.27 points with a standard deviation of 19.40 points. The observed value will be tested against a theoretical value (lower limit -0) of 51.87 points. To account for an estimated appropriate overall drop-out rate, Stryker intends to enroll the calculated number of subjects, multiplied by the number of cleared indications (four) to reflect the underlying subject population adequately.

11.1. DETERMINATION OF SAMPLE SIZE


The determination of sample size is based on benchmark sources and values. To account for an estimated appropriate overall drop-out rate, Stryker intends to enroll the calculated number of subjects, multiplied by the number of cleared indications to reflect the underlying subject population adequately. For details, please see Table 2.

Benchmark and Objectives for Clinical Investigation

Endpoint		The non-inferiority (equal or better) of the ASES Score in relationship to the officially cleared indications in comparison to respective clinical outcome data in the scientific literature.			
Estimated drop-out rate		56% (confirmed by Medical Expert, see Cuff et al. [1])			
Benchmark Sources & Values ASES [points]					
Source		n	Mean	Std. Dev.	Comments No.
No.	Title				
1	Nolan et al. [2] 2011	71	76.1	NA	NA
2	Naveed et al. [3] 2011	49	65	NA	NA
3	Kelly 2 nd et al. [4] 2012	30	71.8	12.7	Revision Arthroplasty in failed arthroplasty
4	Patel et al. [5] 2012	28	66.2	21.3	Revision Arthroplasty in failed arthroplasty
5	Cuff et al. [6] 2012	76	75	NA	NA
6a	Wiater et al. [7] 2014	37	84.9	18.9	Cemented humeral stem
6b	Wiater et al. [7] 2014	64	77.1	18.3	Cement-less humeral stem
7a	Black et al. [8] 2014	33	74	23.8	NA
7b	Black et al. [8] 2014	32	69.7	18	NA
8a	Poon et al. [9] 2014	27	70	NA	Concentric
8b	Poon et al. [9] 2014	23	68	NA	Eccentric
9	Kiet et al. [10] 2015	53	77	19	NA
10	Steen et al. [11] 2015	24	79.9	20.2	NA
11	Wagner et al. [12] 2015	154	65.0	NA	Revision Arthroplasty in failed arthroplasty
12a	Crosby et al. [13] 2017	73	55	25	Stem exchanged
12b	Crosby et al. [13] 2017	29	68	20	Stem retention
13	Vourazeris et al. [14] 2017	86	77.7	NA	NA
14a	Otto et al. [15] 2017	32	58.6	NA	NA
14b	Otto et al. [15] 2017	35	40.8	NA	Revision Arthroplasty in failed arthroplasty
15a	Mollon et al. [16] 2017	428	84.1	17.1	Patients without scapular notching
15b	Mollon et al. [16] 2017	48	78.1	21.8	Patients with scapular notching
16	Holschen et al. [17] 2017	44	65.3	20.9	Revision Arthroplasty in failed arthroplasty
17a	Friedmann et al. [18] 2017	340	86.7	15.8	Patients with subscapularis repair
17b	Friedmann et al. [18] 2017	251	82.1	18.2	Patients without subscapularis repair
18	Williams et al. [19] 2017	17	65.6	NA	Revision Arthroplasty in failed arthroplasty
Identified Cleared Indications					
No.	Indication				
1	Gross rotator cuff deficiency				
2	Degenerative arthritis				
3	Rheumatoid arthritis				
4	Revision of previously failed shoulder joint replacement				
Explorative Meta-Analysis - ASES (single group)					
Acceptance Criteria					
Confidence Interval (CI)				0.95 (95%) two-sided	
Software Used				IBM SPSS V20	

Case Processing Summary						
	Cases					
	Valid		Missing		Total	
	N	Percent	N	Percent	N	Percent
Mean ASES ReUnion_RSA_CER	25	100,0%	0	0,0%	25	100,0%
Std.Dev. ASES ReUnion_RSA_CER	15	60,0%	10	40,0%	25	100,0%

Descriptives				Statistic
Mean ASES ReUnion_RSA_CER	Mean			71,2680
	95% Confidence Interval for Mean	Lower Bound		67,0288
		Upper Bound		75,5072
	Median			71,8000
	Std. Deviation			10,26985
	Minimum			40,80
	Maximum			86,70
	Interquartile Range			12,45
Std.Dev. ASES ReUnion_RSA_CER	Mean			19,4000
	95% Confidence Interval for Mean	Lower Bound		17,7028
		Upper Bound		21,0972
	Median			19,0000
	Variance			9,393
	Std. Deviation			3,06478
	Minimum			12,70
	Maximum			25,00
	Interquartile Range			3,30



Red line: Pooled mean ASES
Blue lines: 95% CI of pooled mean ASES
Black line: Pooled median pooled ASES
Box: Interquartile Range
Green Line: Pooled mean ASES minus pooled std. dev. ASES (19.40)
 $71.27 - 19.40 = 51.87$ points


Benchmark n	25
Benchmark Mean	71.27 (pooled mean ASES [points])
Benchmark Std. Dev.	19.40 (pooled std. dev. ASES [points])
Benchmark Value Non-Inferiority Margin (-θ)	Pooled mean ASES minus pooled std. dev. ASES (19.40) 71.27 – 19.40 = 51.87 points
Software Used	IBM SPSS Sample Power V3.0
IBM SPSS Sample Power Output <p>One goal of the proposed clinical investigation is to test the null hypothesis that the population mean is 71.27 points. The criterion for significance (alpha) has been set at 0.05. The test is 2-tailed, which means that effects in both directions will be interpreted. With the proposed sample size of 10 cases, the clinical investigation will have power of 80.3% to yield a statistically significant result. This computation assumes that the population from which the sample will be drawn has a mean of 71.27 points with a standard deviation of 19.40 points. The observed value will be tested against a theoretical value (constant, non-inferiority margin) of 51.87 points.</p>	
IBM SPSS Sample Power Output – Screenshot 	
<p>Estimated overall drop-out rate is 56% which leads to the requirement of enrolling additional 6 subjects into the clinical investigation.</p>	
Sample Size	Overall number of subjects to be enrolled: 16 subjects (rounded up to 20 subjects)
Overall Sample Size (multiplied by number of indications)	80 subjects

Table 2: Sample Size Justification

In conclusion, the calculated number of subjects to be enrolled (ten) plus the estimated overall drop-out rate of 56%, leads Stryker to enroll 20 subjects per cleared indication (four) to reflect the underlying subject population adequately. As a result, an enrollment target of 80 subjects in total will be aspired.

11.2. ANALYSIS POPULATIONS

It is expected during this clinical investigation only one population for ReUnion RSA System will exist and all subjects will be analyzed “Per Protocol” (PP). However, it cannot be fully avoided that in theory subjects might need to be excluded from the PP population. In this occasion, there will be two groups being fully analyzed to ensure transparency and avoid bias.

The groups are defined as follows:

- **Intent-to Treat Population**

The Intent-to-Treat (ITT) Population is defined to be all enrolled subjects. An enrolled subject is a subject that has signed informed consent, all screening procedures have been successfully completed, is eligible and can receive treatment. The ITT population will not be analyzed for the annual reports and will only be included in the final report.

- **Per Protocol Population**

The PP Population is defined to be all subjects in the ITT Population with no major protocol violations. The protocol deviations that will exclude a subject are as follows:

- The subject does not receive the ReUnion RSA System
- The subject does not meet all eligibility criteria
- The subject has a protocol deviation that is considered likely to affect subject outcomes.

After the clinical investigation has been completed, a review of the data will be conducted to determine which subjects are to be excluded from the PP population.

11.3. ANALYSIS AND EVALUATION

11.3.1. Statistical Analysis

Evaluation elements are defined as the questions on the CRF/eCRF. The SAP lists all evaluation elements and secondary elements which will be based on calculations between two or more evaluation elements.

All quantitative variables, including those based on calculations (secondary elements), will be analyzed with a case summary evaluation before the detailed characteristics and parameters can be evaluated. A case summary contains a listing of the number of valid cases/values, missing cases/values (if any) and total cases/values in the specific analysis. In general, as central position parameter for quantitative variables the mean, median and mode will be analyzed. As variation parameter the standard deviation, 95% confidence interval of the mean, interquartile range and range (based on maximum and minimum) will be calculated. All quantitative variables will be assessed for normality using the Shapiro-Wilk test. For optional visualization of quantitative variables, box-and-whisker plots will be used. Additional analyses like skewness and kurtosis measures or standard errors also are optional.

All qualitative variables, including those based on summaries (secondary elements), will be analyzed listing the proportions, frequencies, column and row totals, and missing proportion, if any.

The SAP reflects this approach and specifies the variables characteristics (quantitative or qualitative) in detail together with the related analysis strategy. This also includes calculation and summaries based on primary elements and the required analysis.

11.3.2. Primary Analysis / Endpoint

The objective of the clinical investigation is to demonstrate the non-inferiority (equal or better) of the ASES Shoulder Score in relationship to the officially cleared indications in comparison to respective clinical outcome data in the scientific literature.

Data collection of ASES Shoulder Score will start pre-operatively and will be collected according to schedule in Table 1. This will be repeated annually in all subjects who have the prosthesis with full or partial implant survival (including all subjects without removal of all endo-prosthesis components).


The 24 months post-operative results for subjects implanted with ReUnion RSA System will be compared to a historical group and results reported by respective clinical outcome data in the scientific literature. The benchmark sources and values will serve as the control group for the ReUnion RSA System subjects.

Higher ASES Shoulder Score results are linked to better subject results and vice versa.

The clinical investigation endpoint is non-inferiority to the control, meaning the clinical investigation result should be equal or better than the control. In this clinical investigation, an equal or better ASES score result means equal or more (\geq). As only results from samples will be captured, results are mostly estimates of the true population parameter. These estimates vary by a certain area, where it is expected that the true population parameter falls within. Based on this, it is required to specify a lower limit for the acceptable difference or zone of indifference, denoted as $-\theta$.

Hypotheses are developed to allow for a comparison of the 24 months post-operative ASES score effectiveness / performance between the two underlying populations. The 24 months post-operative ASES Shoulder Score is the primary endpoint of this clinical investigation. Hypothesis tests will be one-sided with a significance level α of 5%.

Hypothesis	Equations	Interpretation
Null (H_0)	$A - B < -\theta$	Central tendency of A is inferior to the central tendency of B.
	ReUnion RSA System – Control (Benchmark) $< -\theta$	
Alternative (H_1)	$A - B \geq -\theta$	Central tendency of A is non-inferior to the central tendency of B.
	ReUnion RSA System – Control (Benchmark) $\geq -\theta$	



Possible Evidence (p)	Possible Decisions	Possible Conclusions – ASES score
$p\text{-value} > \alpha (0.05)$	Fail to reject null hypothesis (H_0)	ReUnion RSA System $<$ Control (Benchmark) Insufficient evidence to reject the null hypothesis ($H_0: A - B < -\theta$) at the pre-determined significance level of 5%.
$p\text{-value} \leq \alpha (0.05)$	Reject null hypothesis (H_1)	ReUnion RSA System \geq Control (Benchmark) Sufficient evidence to reject the null hypothesis ($H_0: A - B < -\theta$) at the pre-determined significance level of 5%.

To test non-inferiority, the 24 months mean ASES Shoulder Score result of the ReUnion RSA System group will be compared to the mean estimate of the control group (71.27 points).

To be able to identify an acceptable difference or zone of indifference ($-\theta$), the pooled standard deviation of the ASES Shoulder Score result at 24 months post-operative in the control group (Benchmark) was used as lower limit (pooled standard deviation of control is 19.40 points).

The lower maximum acceptable difference ($-\theta$) is 51.87 points (mean of control - θ or $71.27 - 19.40 = 51.87$ points).

Based on the underlying distribution of the data and the result of the normality assessment, either the parametric one-sample t-test or the non-parametric one-sample sign test will be used to compare the 24 months post-operative ASES Shoulder Score results of the ReUnion RSA System against the value of 51.87 points.

11.3.3. Secondary Endpoints

The incidence of device-related AEs and implant survivorship will be assessed up to 10 years after the index procedure and monitored through collection and analyses. These analyses will be part of the annual and final reports.

Furthermore, time to (earliest) device-related AEs will be analyzed as well. For analysis of the time to the (earliest) device-related AEs as well as the time to secondary procedure (revision, removal, reoperation), the Kaplan-Meier method will be used. The time between surgery and the last available assessment will be used together with the time between date of surgery and the date of secondary procedure. Considered variables, the level of measurement and the planned analysis steps are listed in detail in the SAP.

11.3.4. Additional Analyses

Additional analyses are outlined in the subsequent sections. Analysis details (variables, level of measurement, planned steps) are listed in-depth in the SAP.

- **Mortality**
For analysis of the time to death or mortality, the Kaplan-Meier method will be used. The times between surgery and the last available assessment will be used together with the times between date of surgery and the date of death. This analysis will be part of the annual and final reports.
- **Total ASES Shoulder Score – Within subject changes by visit**
The within subject score changes of the ASES Shoulder Score from visit to visit will be analyzed to help identify the changes on the subject level. This analysis will be part of the annual and final reports.

11.4. MISSING DATA / SAP DEVIATIONS

The intent is to collect as complete a dataset as possible. Nevertheless, in some situations missing data cannot be avoided. The reports and tables therefore will show the number and percentage of missing cases for each analyzed variable in relation to the enrolled cases for each post-operative assessment.

Any deviations from the SAP will be listed in the annual or final reports.

11.5. REPORTS

11.5.1. Interim Analysis and Reports

Interim analyses will be performed on a yearly basis. The progress of the clinical investigation will be reported together with the interim results on the variable level according to the analysis plan.

The analysis of the primary endpoint will be part of the related interim / annual report when all subjects have completed the 24 months post-operative including the ASES Shoulder Score.

11.5.2. Final Analysis and Reports

The full final report with complete analysis including progress and conduct reporting will be created at the end of this clinical investigation.

12. Clinical Investigation Plan Deviations

A CIP deviation is a departure from the approved CIP that is not implemented or intended as a systemic change. All CIP deviations are recorded and reported to each site's IRB in accordance with the respective site's IRB policies.

13. Adverse Events

As this CIP is being carried out to satisfy the post-market requirements to support safety and efficacy/performance according to the European Medical Device Regulation (EU MDR), categorization and definition of device-related adverse events will follow the guidelines outlined in the EU MDR as "incident" reporting.

13.1. DEFINITIONS

- **An Adverse Device Effect (ADE)** is defined as any untoward or unintended response to the clinical investigation treatment; and/or a medical response which may have a causal relationship to the treatment.
- **An Incident** is defined as any malfunction or deterioration in the characteristics or performance of a device made available on the market, including use-error due to ergonomic features, as well as any inadequacy in the information supplied by the manufacturer and any undesirable side-effect.
- **A Serious Adverse Device Effect (SADE)** is defined as any ADE that results in consequences characteristic of a SAE or might lead to the consequences if suitable action or intervention is not taken; causes considerable interference with the subject's usual activities and may be incapacitating or life-threatening, including those events resulting in a subject's disability or permanent damage, or required intervention to prevent disability or permanent damage; results in a life-threatening illness or injury; and/or results in death (fatal).
- **A Serious Incident** is defined as any incident that directly or indirectly led, might have led or might lead to any of the following:
 - the death of a patient, user or other person;
 - the temporary or permanent serious deterioration of a patient's, user's or other person's state of health;
 - a serious public health threat
- **An Unanticipated Adverse Device Effect (UADE)** is defined as an AE not described in the informed consent, CIP or device labeling which has resulted in any of the consequences of a SAE or which might have led to any of the consequences of a SAE if suitable action had not been taken, intervention had not occurred, or if circumstances had been less opportune.

13.2. ADVERSE EVENT SEVERITY

The severity of all AEs is assessed by the Investigator utilizing the following categories:

- **Mild:** The AE is transient and easily tolerated by the subject.
- **Moderate:** The AE causes the subject discomfort and interrupts the subject's usual activities.
- **Severe:** The AE causes considerable interference with the subject's usual activities and may be incapacitating or life-threatening, including those events resulting in a subject's disability or permanent damage, and/or required intervention to prevent permanent disability or damage.
- **Life-Threatening:** The AE results in a life-threatening illness or injury.
- **Fatal:** The AE results in death.

13.3. RELATIONSHIP TO THE DEVICE

Only events considered possibly, probably or definitely related to the device will be captured for this clinical investigation.

13.4. ADVERSE EVENT/INCIDENT REPORTING

In the event a SADE, UADE or serious incident occurs, the Investigator is required to notify the Sponsor within 48 hours of being made aware of the event. The Investigator also is required to notify their IRB in accordance with the policies of their local laws and regulations.

13.5. FORESEEABLE ADEs, SADEs and INCIDENTS

ADEs, SADEs and incidents which may be expected as part of the surgical intervention include:

- *Perioperative complications*
 - Malpositioning of the humeral or glenoidal component
 - Oversizing of implant components
 - Undersizing of implant components
 - Intraoperative fracture of the humerus or glenoid
 - Humeral cement leakage, if applicable
 - Nerve injury, mild (minor motor or sensory loss, or spontaneous recovery)
 - Nerve injury, severe (significant motor or sensory loss or requiring surgical revision)
 - Vessel injury
 - Tendon injury
 - Wound complications (e.g. hematoma, wound healing disturbances)
 - Superficial infection
 - Deep infection
 - Deep vein thrombosis
 - Pulmonary embolism
- *Complications in the follow-up period*
 - Implant dislocation
 - Other subluxation or instability, symptomatic
 - Implant component dissociation, humeral cup from humeral stem
 - Implant component dissociation, humeral insert from humeral cup
 - Implant component dissociation, glenosphere from baseplate
 - Implant fretting or crevice corrosion
 - Implant breakage/wear, humeral or glenoidal
 - Implant loosening, humeral
 - Implant loosening, glenoidal with or without screw breakage
 - Radiographic lucency, humeral or glenoidal
 - Pain related to the implant, severe
 - Late infection (e.g. hematogenous or protracted)

- Periprosthetic fracture of the humerus or glenoid
- Stress fracture of the acromion or the scapular neck
- Stress fracture of the coracoid
- Scapular notching, asymptomatic or symptomatic
- Stiffness
- Heterotopic ossification, asymptomatic or symptomatic
- Healing disturbances of subscapularis tenotomy, if applicable
- Non-union of lesser tuberosity osteotomy, if applicable
- Malunion of lesser tuberosity osteotomy, if applicable

14. Reoperations

Reoperations and reason(s) for reoperations will be collected throughout the course of the clinical investigation. A reoperation may include but not limited to irrigation and debridement, revision surgery and/or implant removal.

15. Ethics

This clinical investigation is to be conducted according to International Conference of Harmonisation of Good Clinical Practice (ICH-GCP), applicable regulations, institutional research policies and procedures, Declaration of Helsinki and in compliance with the CIP. Investigators will be trained on the clinical investigation devices and surgical techniques prior to implanting clinical investigation subjects.

This CIP and any amendments will be submitted to a properly constituted independent ethics board, in agreement with local legal prescriptions, for formal approval of the clinical investigation conduct. The decision of the ethics board concerning the conduct of the clinical investigation will be made in writing to the Investigator before commencement of this clinical investigation. Clinical investigations shall not begin until the governing regulatory authority has provided full, unconditional approval. Off-label use of the ReUnion RSA System is not permitted.

15.1. INSTITUTION REVIEW BOARD (IRB)

IRB approval will be obtained at each of the investigational sites prior to enrolling clinical investigation subjects at that site. In addition, any SAEs and UADEs that meet the reporting criteria of the IRB, will be reported to the IRB. During the clinical investigation, the Investigator should promptly provide written reports to the IRB of any changes that affect the conduct of the clinical investigation and/or increase risk to the subjects, unless otherwise submitted by the Sponsor.

15.2. INFORMED CONSENT

The Investigator, or qualified clinical investigation personnel designated to perform this task, will explain the nature of the clinical investigation to the subject, and answer all questions regarding participation in this clinical investigation. Prior to any clinical investigation procedures being performed, the informed consent form (ICF) will be reviewed, signed and dated by the subject, and by the person administering the informed consent. A copy of the ICF will be given to the subject, and the original will be placed in the subject's clinical investigation records. Subjects will need to sign updated versions of the ICF if required by the Investigator's IRB during the clinical investigation.

16. Data Collection Process

The Sponsor will collect clinical data for this clinical investigation utilizing eCRFs through an EDC system. All data entered in the eCRFs are supported by source documentation. All clinical data is entered into the EDC system by designated personnel at each of the Investigator sites.

16.1. RADIOGRAPHS

All radiographs shall be de-identified of personal health information. Radiographs will be uploaded as DICOM images into the EDC system. The radiologic analysis shall be based at minimum on an axillary lateral radiograph and an anteroposterior radiograph. Additional radiographs (e.g., 30° to 40° posterior oblique radiographs in internal and external rotation) may be used where available. For measurements, all digital radiographs shall be sized to 100%, based on the diameter of the glenosphere or other suitable reference sizes.

17. Clinical Investigation Monitoring

It is the responsibility of the Investigator to oversee the safety of the clinical investigation at his/her site, to include the careful assessment and appropriate reporting of AEs/incidents as noted above as well as the implementation of site data safety. The Sponsor, or designee, will monitor the sites to ensure informed consent has been documented appropriately, to ensure the information documented on the completed CRFs match the medical records and to resolve any differences. The Sponsor will take all steps necessary to ensure data integrity. The Sponsor also will review significant new information, including UADEs and ensure that such information is provided to all Investigators, their IRBs, and applicable regulatory authorities. Additionally, a quality assurance check will be performed to ensure the investigator is complying with the CIP and applicable regulations in the collection of all clinical investigation data.

18. Data Handling and Record Keeping

Information about clinical investigation subjects will be kept confidential. In the event a subject revokes authorization to collect or use protected health information, the Investigator, by regulation, retains the ability to use all information collected prior to the revocation of subject authorization. The Health Insurance Portability and accountability Act (HIPAA) will apply to ensure data protection and document anonymization. Records are to be stored in a secure location. Retention of records shall be maintained through the clinical investigation duration as well as specified years following the clinical investigation completion as required by local regulatory authority.

19. Reports

Analysis will be performed, and interim reports will be prepared on a yearly basis. Upon the completion of all subjects' final post-operative assessment, data freeze will occur, and the final report will be prepared.

20. Completion of the Clinical Investigation

The Investigator will conduct this clinical investigation in compliance with the CIP and will complete the clinical investigation within the timeframe specified in the contract. Continuation of the clinical investigation beyond this time must be mutually agreed upon in writing by both the Investigator and Stryker. The Investigator will provide a summary of the clinical investigation results in accordance with the IRB guidelines.

Stryker may terminate this clinical investigation prematurely, either in its entirety or at this site, for reasonable cause provided that written notice is submitted a reasonable time in advance of the intended termination. The Investigator also may terminate the clinical investigation at their site for reasonable cause, after providing written notice to Stryker a reasonable time in advance of the intended termination. If Stryker terminates the clinical investigation for safety reasons, then it will immediately notify the Investigator by telephone and subsequently provide written instructions for clinical investigation termination.

21. Essential Documents

All essential documentation will be stored as specified under the Sponsor's Standard Operating Procedures.

22. Publication Policy

Refer to the clinical investigation agreement for the publication policy.

23. References

1. Cuff DJ, Pupello DR, Santoni BG, Clark RE, Frankle MA. Reverse Shoulder Arthroplasty for the Treatment of Rotator Cuff Deficiency: A Concise Follow-up, at a Minimum of 10 Years, of Previous Reports. *J Bone Joint Surg Am.* 2017 Nov 15; 99(22):1895-1899.
2. Nolan BM1, Ankersen E, Wiater JM. Reverse total shoulder arthroplasty improves function in cuff tear arthropathy. *Clin Orthop Relat Res.* 2011 Sep;469(9):2476-82.
3. Naveed MA1, Kitson J, Bunker TD. The Delta III reverse shoulder replacement for cuff tear arthropathy: a single-centre study of 50 consecutive procedures. *J Bone Joint Surg Br.* 2011 Jan;93(1):57-61.
4. Kelly JD 2nd, Zhao JX, Hobgood ER, Norris TR. Clinical results of revision shoulder arthroplasty using the reverse prosthesis. *J Shoulder Elbow Surg.* 2012 Nov;21(11):1516-25.
5. Patel DN, Young B, Onyekwelu I, Zuckerman JD, Kwon YW. Reverse total shoulder arthroplasty for failed shoulder arthroplasty. *J Shoulder Elbow Surg.* 2012 Nov;21(11):1478-83.
6. Cuff D, Clark R, Pupello D, Frankle M. Reverse shoulder arthroplasty for the treatment of rotator cuff deficiency: a concise follow-up, at a minimum of five years, of a previous report. *J Bone Joint Surg Am.* 2012 Nov 7;94(21):1996-2000
7. Wiater JM, Moravek JE Jr, Budge MD, Koueiter DM, Marcantonio D, Wiater BP. Clinical and radiographic results of cementless reverse total shoulder arthroplasty: a comparative study with 2 to 5 years of follow-up. *J Shoulder Elbow Surg.* 2014 Aug;23(8):1208-14.
8. Black EM, Roberts SM, Siegel E, Yannopoulos P, Higgins LD, Warner JJ. Reverse shoulder arthroplasty as salvage for failed prior arthroplasty in patients 65 years of age or younger. *Shoulder Elbow Surg.* 2014 Jul;23(7):1036-42.
9. Poon PC, Chou J, Young SW, Astley T. A comparison of concentric and eccentric glenospheres in reverse shoulder arthroplasty: a randomized controlled trial. *J Bone Joint Surg Am.* 2014 Aug 20;96(16):e138.
10. Zafra M, Uceda P, Flores M, Carpintero P. Reverse total shoulder replacement for nonunion of a fracture of the proximal humerus. *Bone Joint J.* 2014 Sep;96-B(9):1239-43.
11. Steen BM, Cabezas AF, Santoni BG, Hussey MM, Cusick MC, Kumar AG, Frankle MA. Outcome and value of reverse shoulder arthroplasty for treatment of Glenohumeral osteoarthritis: a matched cohort. *J Shoulder Elbow Surg.* 2015 Mar 11. pii: S1058-2746(15)00043-9.

12. Wagner ER, Houdek MT, Elhassan BT, Sanchez-Sotelo J, Cofield RH, Sperling JW. What Are Risk Factors for Intraoperative Humerus Fractures During Revision Reverse Shoulder Arthroplasty and Do They Influence Outcomes? *Clin Orthop Relat Res*. 2015 Jul 11.
13. Crosby LA, Wright TW, Yu S, Zuckerman JD. Conversion to Reverse Total Shoulder Arthroplasty with and without Humeral Stem Retention: The Role of a Convertible-Platform Stem. *J Bone Joint Surg Am*. 2017 May 3;99(9):736-742.
14. Vourazeris JD, Wright TW, Struk AM, King JJ, Farmer KW. Primary reverse total shoulder arthroplasty outcomes in patients with subscapularis repair versus tenotomy. *J Shoulder Elbow Surg*. 2017 Mar;26(3):450-457.
15. Otto RJ, Clark RE, Frankle MA. Reverse shoulder arthroplasty in patients younger than 55 years: 2- to 12-year follow-up. *J Shoulder Elbow Surg*. 2017 May;26(5):792-797.
16. Mollon B, Mahure SA, Roche CP, Zuckerman JD3. Impact of scapular notching on clinical outcomes after reverse total shoulder arthroplasty: an analysis of 476 shoulders. *J Shoulder Elbow Surg*. 2017 Jul;26(7):1253-1261.
17. Holschen M1, Franetzki B2, Witt KA2, Liem D3, Steinbeck J2. Is reverse total shoulder arthroplasty a feasible treatment option for failed shoulder arthroplasty? A retrospective study of 44 cases with special regards to stemless and stemmed primary implants. *Musculoskelet Surg*. 2017 Aug;101(2):173-180.
18. Friedman RJ1, Flurin PH2, Wright TW3, Zuckerman JD4, Roche CP5. Comparison of reverse total shoulder arthroplasty outcomes with and without subscapularis repair. *J Shoulder Elbow Surg*. 2017 Apr;26(4):662-668.
19. Williams PN1, Trehan SK2, Tsouris N2, Dines JS2, Dines DM2, Craig EV3, Gulotta LV2, Warren RF2. Functional Outcomes of Modular Conversion of Hemiarthroplasty or Total to Reverse Total Shoulder Arthroplasty. *HSS J*. 2017 Jul;13(2):102-107.

24. Clinical Investigation Plan Signature Page

ReUnion RSA System

I have read this Clinical Investigation Plan and agree that this clinical investigation is ethical. I agree to conduct this clinical investigation in accordance with this Clinical Investigation Plan, as well as all applicable regulations and guidelines. I agree to maintain the confidentiality of all information received or developed in connection with this Clinical Investigation Plan.

Signature of Investigator

Date of Signature

Name of Investigator (Printed)

25. Document Version History

Version	Effective Date	Description	Revised/Created by
1	12Oct2018	Initial version	Lindsay Mattfolk

Certificate Of Completion

Envelope Id: D3DE7B77CB2948FC801BC84171E5CDCE	Status: Completed
Subject: Please DocuSign: 12OCT2018_ReUnion RSA_Protocol Development.pdf	
Source Envelope:	
Document Pages: 24	Signatures: 4
Certificate Pages: 5	Initials: 0
AutoNav: Enabled	
Enveloped Stamping: Disabled	
Time Zone: (UTC-08:00) Pacific Time (US & Canada)	Envelope Originator: Lindsay Mattfolk 325 Corporate Drive New Jersey, NC 07430 lindsay.mattfolk@stryker.com IP Address: 24.94.18.137

Record Tracking


Status: Original 10/15/2018 7:10:45 AM	Holder: Lindsay Mattfolk lindsay.mattfolk@stryker.com	Location: DocuSign
---	--	--------------------

Signer Events

Signer Events	Signature	Timestamp
Claudia Beime claudia.beime@stryker.com Stryker Security Level: Email, Account Authentication (Required)	<i>Claudia Beime</i> Signature Adoption: Pre-selected Style Signature ID: CF9D5AE3-52C5-45BA-B2C9-FF408F88CDAF Using IP Address: 165.225.72.71 With Signing Authentication via DocuSign password With Signing Reasons (on each tab): I approve this document	Sent: 10/15/2018 7:16:30 AM Viewed: 10/16/2018 12:37:54 AM Signed: 10/16/2018 12:42:59 AM


Electronic Record and Signature Disclosure:

Not Offered via DocuSign


Georgia Mitchell georgia.mitchell@stryker.com Security Level: Email, Account Authentication (Required)	 Signature Adoption: Pre-selected Style Signature ID: B2566DB4-9823-4B21-8287-8F1BD6240617 Using IP Address: 69.141.177.192 With Signing Authentication via DocuSign password With Signing Reasons (on each tab): I approve this document	Sent: 10/15/2018 7:16:30 AM Viewed: 10/15/2018 7:55:06 AM Signed: 10/15/2018 7:56:04 AM
--	---	---

Electronic Record and Signature Disclosure:

Accepted: 10/15/2018 7:55:06 AM
ID: 247762e6-8dac-4154-a9b1-a68c98d82aeb

Signer Events	Signature	Timestamp
Helmut Mueckter helmut.mueckter@stryker.com Security Level: Email, Account Authentication (Required)	 Signature Adoption: Pre-selected Style Signature ID: C24B3BDD-5E94-46ED-8398-FBC9B3EE30A7 Using IP Address: 165.225.80.49 With Signing Authentication via DocuSign password With Signing Reasons (on each tab): Ich habe dieses Dokument geprüft	Sent: 10/15/2018 7:16:31 AM Viewed: 10/15/2018 8:52:22 AM Signed: 10/15/2018 8:53:42 AM

Electronic Record and Signature Disclosure:
 Accepted: 4/27/2018 7:50:14 AM
 ID: 11dda966-5c40-496b-a7d0-1b88be367a95

Jemin Dedania jemin.dedania@stryker.com Security Level: Email, Account Authentication (Required)	 Signature Adoption: Pre-selected Style Signature ID: 5B8D8CF9-5DB1-4E5E-A800-D11CDCAB128D Using IP Address: 69.74.233.66 With Signing Authentication via DocuSign password With Signing Reasons (on each tab): I approve this document	Sent: 10/15/2018 7:16:32 AM Viewed: 10/16/2018 6:54:40 AM Signed: 10/16/2018 6:56:15 AM
--	---	---

Electronic Record and Signature Disclosure:
 Accepted: 10/16/2018 6:54:40 AM
 ID: 9a092ab3-fa45-4d6f-a605-01bbf03d9206

In Person Signer Events	Signature	Timestamp
Editor Delivery Events	Status	Timestamp
Agent Delivery Events	Status	Timestamp
Intermediary Delivery Events	Status	Timestamp
Certified Delivery Events	Status	Timestamp
Carbon Copy Events	Status	Timestamp
Notary Events	Signature	Timestamp
Envelope Summary Events	Status	Timestamps
Envelope Sent	Hashed/Encrypted	10/15/2018 7:16:33 AM
Certified Delivered	Security Checked	10/16/2018 6:54:40 AM
Signing Complete	Security Checked	10/16/2018 6:56:15 AM
Completed	Security Checked	10/16/2018 6:56:15 AM
Payment Events	Status	Timestamps
Electronic Record and Signature Disclosure		

ELECTRONIC RECORD AND SIGNATURE DISCLOSURE

From time to time, Stryker Corporation - Trauma & Extremities - Part 11 (we, us or Company) may be required by law to provide to you certain written notices or disclosures. Described below are the terms and conditions for providing to you such notices and disclosures electronically through your DocuSign, Inc. (DocuSign) Express user account. Please read the information below carefully and thoroughly, and if you can access this information electronically to your satisfaction and agree to these terms and conditions, please confirm your agreement by clicking the 'I agree' button at the bottom of this document.

Getting paper copies

At any time, you may request from us a paper copy of any record provided or made available electronically to you by us. For such copies, as long as you are an authorized user of the DocuSign system you will have the ability to download and print any documents we send to you through your DocuSign user account for a limited period of time (usually 30 days) after such documents are first sent to you. After such time, if you wish for us to send you paper copies of any such documents from our office to you, you will be charged a \$0.00 per-page fee. You may request delivery of such paper copies from us by following the procedure described below.

Withdrawing your consent

If you decide to receive notices and disclosures from us electronically, you may at any time change your mind and tell us that thereafter you want to receive required notices and disclosures only in paper format. How you must inform us of your decision to receive future notices and disclosure in paper format and withdraw your consent to receive notices and disclosures electronically is described below.

Consequences of changing your mind

If you elect to receive required notices and disclosures only in paper format, it will slow the speed at which we can complete certain steps in transactions with you and delivering services to you because we will need first to send the required notices or disclosures to you in paper format, and then wait until we receive back from you your acknowledgment of your receipt of such paper notices or disclosures. To indicate to us that you are changing your mind, you must withdraw your consent using the DocuSign 'Withdraw Consent' form on the signing page of your DocuSign account. This will indicate to us that you have withdrawn your consent to receive required notices and disclosures electronically from us and you will no longer be able to use your DocuSign Express user account to receive required notices and consents electronically from us or to sign electronically documents from us.

All notices and disclosures will be sent to you electronically

Unless you tell us otherwise in accordance with the procedures described herein, we will provide electronically to you through your DocuSign user account all required notices, disclosures, authorizations, acknowledgements, and other documents that are required to be provided or made available to you during the course of our relationship with you. To reduce the chance of you inadvertently not receiving any notice or disclosure, we prefer to provide all of the required notices and disclosures to you by the same method and to the same address that you have given us. Thus, you can receive all the disclosures and notices electronically or in paper format through the paper mail delivery system. If you do not agree with this process, please let us know as described below. Please also see the paragraph immediately above that describes the consequences of your electing not to receive delivery of the notices and disclosures electronically from us.

How to contact Stryker Corporation - Trauma & Extremities - Part 11:

You may contact us to let us know of your changes as to how we may contact you electronically, to request paper copies of certain information from us, and to withdraw your prior consent to receive notices and disclosures electronically as follows:

To contact us by email send messages to: rebecca.gibson@stryker.com

To advise Stryker Corporation - Trauma & Extremities - Part 11 of your new e-mail address

To let us know of a change in your e-mail address where we should send notices and disclosures electronically to you, you must send an email message to us at rebecca.gibson@stryker.com and in the body of such request you must state: your previous e-mail address, your new e-mail address. We do not require any other information from you to change your email address..

In addition, you must notify DocuSign, Inc to arrange for your new email address to be reflected in your DocuSign account by following the process for changing e-mail in DocuSign.

To request paper copies from Stryker Corporation - Trauma & Extremities - Part 11

To request delivery from us of paper copies of the notices and disclosures previously provided by us to you electronically, you must send us an e-mail to rebecca.gibson@stryker.com and in the body of such request you must state your e-mail address, full name, US Postal address, and telephone number. We will bill you for any fees at that time, if any.

To withdraw your consent with Stryker Corporation - Trauma & Extremities - Part 11

To inform us that you no longer want to receive future notices and disclosures in electronic format you may:

- i. decline to sign a document from within your DocuSign account, and on the subsequent page, select the check-box indicating you wish to withdraw your consent, or you may;
- ii. send us an e-mail to rebecca.gibson@stryker.com and in the body of such request you must state your e-mail, full name, IS Postal Address, telephone number, and account number. We do not need any other information from you to withdraw consent.. The consequences of your withdrawing consent for online documents will be that transactions may take a longer time to process..

Required hardware and software

Operating Systems:	Windows2000? or WindowsXP?
Browsers (for SENDERS):	Internet Explorer 6.0? or above
Browsers (for SIGNERS):	Internet Explorer 6.0?, Mozilla FireFox 1.0, NetScape 7.2 (or above)
Email:	Access to a valid email account
Screen Resolution:	800 x 600 minimum
Enabled Security Settings:	<ul style="list-style-type: none">•Allow per session cookies•Users accessing the internet behind a Proxy Server must enable HTTP 1.1 settings via proxy connection

****** These minimum requirements are subject to change. If these requirements change, we will provide you with an email message at the email address we have on file for you at that time providing you with the revised hardware and software requirements, at which time you will

have the right to withdraw your consent.

Acknowledging your access and consent to receive materials electronically

To confirm to us that you can access this information electronically, which will be similar to other electronic notices and disclosures that we will provide to you, please verify that you were able to read this electronic disclosure and that you also were able to print on paper or electronically save this page for your future reference and access or that you were able to e-mail this disclosure and consent to an address where you will be able to print on paper or save it for your future reference and access. Further, if you consent to receiving notices and disclosures exclusively in electronic format on the terms and conditions described above, please let us know by clicking the 'I agree' button below.

By checking the 'I Agree' box, I confirm that:

- I can access and read this Electronic CONSENT TO ELECTRONIC RECEIPT OF ELECTRONIC RECORD AND SIGNATURE DISCLOSURES document; and
- I can print on paper the disclosure or save or send the disclosure to a place where I can print it, for future reference and access; and
- Until or unless I notify Stryker Corporation - Trauma & Extremities - Part 11 as described above, I consent to receive from exclusively through electronic means all notices, disclosures, authorizations, acknowledgements, and other documents that are required to be provided or made available to me by Stryker Corporation - Trauma & Extremities - Part 11 during the course of my relationship with you.